LER2583 Ages 8+

## Agar **Preparation & Use**

Agar is a gelatinous material extracted from certain marine algae. It is used for gelling foods and as a bacteria culture medium. Nutrient agar will resemble a solution something like Jell-O® with beef soup. It is not edible and should not be ingested.

## Contents:

2 plastic petri dishes 2 packets of Nutrient Agar (Use one packet per petri dish)

## Instructions:

- 1. Boil 2 tablespoons of distilled water.
- Add entire contents of one package of Nutrient Agar to the boiling water.
- 3. solution until the agar completely dissolved. The solution will start to resemble a jelly mixture.
- 4. Pour the agar solution into one petri dish.
- 5. Refrigerate to gel (about 30 minutes). After the solution turns into a waxy 6.
- gel, turn the dish over and draw pie-cut sections on the bottom and up the sides. Number each section.
- Re-draw the sections on paper to record which sample is in each section.
- Next, use a clean toothpick to collect 8. a sample of what you want to test by rubbing it on the sample's surface, then rub it on one section of the agar.
- Cover the dish and leave it in a warm, 9. draft-free area for 24 to 48 hours.
- 10. Leaving the dish open and exposed to the air is another way to collect samples floating around on breeze. You can even sneeze on gel and see what grows.
- 11. By turning the dish upside down it allows for less moisture to collect.
- best to use a magnifier to See your dealer see the results. magnifiers available various from Learning Resources.



## WARNING! NOT A FOOD PRODUCT. DO NOT INGEST.



Agar made in USA

Free Manuals Download Website

http://myh66.com

http://usermanuals.us

http://www.somanuals.com

http://www.4manuals.cc

http://www.manual-lib.com

http://www.404manual.com

http://www.luxmanual.com

http://aubethermostatmanual.com

Golf course search by state

http://golfingnear.com

Email search by domain

http://emailbydomain.com

Auto manuals search

http://auto.somanuals.com

TV manuals search

http://tv.somanuals.com