# Leica Abbe Mark II Refractometer

Models 13104800/13104810/  
13104940/13104950

Instruction Manual  
Gebrauchsanleitung  
Manuel d'utilisation  
Manual de instrucciones

**Leica**  
MICROSYSTEMS


# Leica ABBE MARK II Refractometer

## Instruction Manual

### Contents

1.0	INTRODUCTION .....	3
2.0	UNPACKING .....	3
3.0	POWER REQUIREMENTS .....	3
4.0	SET UP .....	3
4.1	Equipment Installation .....	3
5.0	CONTROLS .....	3
5.1	Variable Conditions in Taking Readings .....	4
5.2	Calibrating With Distilled Water .....	4
6.0	OPERATING INSTRUCTIONS .....	4
7.0	LEICA ABBE MARK II PLUS REFRACTOMETER .....	5
7.1	Computer Interface .....	5
7.2	Setting the Baud Rate .....	5
7.3	Printer Interface .....	5
7.4	Custom Programming .....	5
8.0	REFRACTOMETER USE, TECHNIQUES AND REFRACTION PRINCIPLES .....	6
8.1	High Index/High Dispersion Samples .....	6
9.0	CALIBRATION WITH TEST GLASS .....	6
10.0	ILLUMINATOR BULB REPLACEMENT .....	7
10.1	Fuses .....	8
11.0	CARE AND CLEANING .....	8
12.0	SERVICING .....	8
13.0	SPECIFICATIONS .....	8
13.1	The following specifications refer only to the LEICA ABBE MARK II PLUS Refractometer .....	9

### Explanation of Warning Symbols

 5X20 250mA T/250V (Type 10480/10494)

 **CAUTION:** Risk of electric shock

 **CAUTION:** (Refer to accompanying documents)

## 1.0 INTRODUCTION

**NOTE:** Please read these instructions carefully before attempting to use the instrument.

This Instruction Manual is designed specifically for the LEICA ABBE MARK II Refractometer and describes proper set up, use and routine maintenance of the instrument. **Sections 7.0 - 7.4 are applicable only to the LEICA ABBE MARK II Plus Refractometer.**

The LEICA ABBE MARK II Refractometer is a durable bench-top instrument. It provides efficient, reliable and repeatable readings in three modes:

- % Solids (°Brix)
- % Solids (°Brix) temperature compensated or (BX-TC) and
- Refractive index ( $n_D$ )

The LEICA ABBE MARK II Plus Refractometer can accommodate a custom program on the % Solids (°Brix) scale. This can be accomplished with pre-programming or retrofitting at the factory.

## 2.0 UNPACKING

The LEICA ABBE MARK II Refractometer has been carefully packaged for shipping. Before setting aside the packing materials, check the box(es) to assure that all items ordered have been received. Among the accessories in the refractometer box will be:

- a refractive index test glass (Catalog No. 10451)
- a 1/16" ( No. 4) Allen wrench (Catalog No. X30600)
- a bottle of 1- Bromonaphthaline for measuring solids (Catalog No. 10452)
- a dust cover (Catalog No. 10450-140)

Check the instrument and model number against the original order before discarding any packing materials. The original packing materials may be saved for future storage or transportation of the refractometer.

## 3.0 POWER REQUIREMENTS

Models 13104800/13104940	115~ +/-10%, 50/60 Hz
Models 13104810/13104950	230~ +/-10%, 50/60 Hz

- Input power 25P max
- Indoor use only
- Altitude up to 2000m
- Temperature range 5°C to 40°C
- Pollution degree 2
- Installation category 2
- Max relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.

## 4.0 SET UP

The LEICA ABBE MARK II Refractometer is easy to use. Very little training is needed to set up and operate the instrument to obtain temperature and refractive index measurements.

### 4.1 Equipment Installation

The LEICA ABBE MARK II Refractometer must be used on a hard, stable and level surface such as a laboratory table or work bench.

 **CAUTION:** As a safety precaution, the power cord has been grounded to the instrument base. Never use a 2-prong adapter between the power cord and the power source; it will render the grounding feature ineffective.

**CAUTION:** Use of this equipment other than in the manner specified in this manual may result in reduction of safety protection.

## 5.0 CONTROLS

THE LEICA ABBE MARK II Refractometer controls are identified in Figure 1.

**NOTE:** The mode selector switch no longer functions as the on/off switch. The ON/OFF switch is located on the lower left front cover of the instrument above the Leica logo.

### Power Switch:

The Power Switch is located in the lower right corner of the front panel (Figure 1). | = on 0 = off

The mode selector switch is used to select measurement modes:

The modes are:

- % Solids-TC - Temperature Compensated % Solids as Sucrose ( $^{\circ}$ Brix) per International Commission for Uniform Methods of Sugar Analysis Refractive Index Tables.
- $n_D$  - Refractive Index
- % Solids - Non-Temperature Compensated % Solids as Sucrose

The reticle adjustment access hole permits vertical adjustment of the crosshair during calibration procedures.

Depress the read display button to obtain measurement values.

Depressing the temperature display button will provide a reading of the prism temperature.

Readings are shown digitally in the display window.

### 5.1 Variable Conditions in Taking Readings

Distilled water readings above 20°C can be obtained in the refractive index % Solids-TC mode, assuming the instrument is properly calibrated. Correct readings cannot be obtained in the non-temperature compensated % Solids mode when the water sample temperature is above or below 20°C. Refractive index ( $n_D$ ) readings for distilled water can be obtained from 0°C to 80°C (see Chart 1 for values between 18-30°C).

Refractive index is inversely proportional to temperature. As temperature increases, the solution density decreases and the refractive index or % Solids reading will decrease. As temperature decreases the solution density will increase and the refractive index (% Solids) will increase. The % Solids-TC mode automatically corrects the reading for temperature based on the temperature coefficient of sucrose in solution. This coefficient is typical of most aqueous solutions.

### 5.2 Calibrating With Distilled Water

To calibrate:

1. Place a distilled water sample on the measuring prism surface and let it stand undisturbed for 3-5 minutes to obtain temperature stability. Turn the mode selector switch to the refractive index ( $n_D$ ) position.
2. Focus the eyepiece until the crosshair clearly appears.
3. Rotate the adjustment control clockwise until the shadow line nearly meets the crosshair. Use the dispersion correction wheel to obtain the sharpest possible color-free shadow line. Adjust the

shadow line to meet the crosshair. Press the Read button and record the  $n_D$  value, then press the temperature button and record the temperature.

4. If the instrument is properly calibrated, the reading will be 1.3330 less 0.0001 for each degree above 20°C or plus 0.0001 for each degree below 20°C.

If the instrument is not calibrated, determine the proper reading (Chart 1) and rotate the adjustment control and depress the READ display button until the display reads correctly. (The shadow line will no longer intersect the crosshair on the reticle.)

5. Insert a 1/16" Allen wrench into the reticle adjustment access hole (Figure 2) located in the dispersion correction wheel. Turn the wrench to bring the crosshair into proper position on the shadow line.

**NOTE:** The final movement made with the adjustment correction control wheel should always be in a clockwise direction. The shadow line will move upward.

6. Turn the mode selector to the BX-TC position. Press the READ button. An LED reading of 000.0 or 0.1 is normal. If a negative sign and a decimal point appear, no reading can be given.

Go on to the next step.

If a "No Reading" condition persists after repeated READ attempts, rotate the adjustment control until the display reads 000.0 or 0.1. Only a very slight adjustment will be necessary. Readjust the reticle crosshair to the shadow line.

The display should now read 000.0.

Repeat the steps until the display reads 000.0.

## 6.0 OPERATING INSTRUCTIONS

To begin operation:

1. Open the prism assembly cover (Figure 4) and remove the protective lens tissue. Clean the measuring prism surface with alcohol and then distilled water. Wipe with a clean soft cloth or lens tissue, but do not wipe the measuring prism surface while it is dry. Any residue left on the prism can affect the accuracy of readings.
2. Turn the mode selector to the desired mode position.
3. Apply a sample solution to the measuring prism surface using a glass dropper. Close the prism cover and position the illumination arm so that the exposed face of the upper prism is fully illuminated.
4. Turn the dispersion correction wheel so that the crosshair adjustment access hole is at the six o'clock position.


Figure 1


Figure 2

5. Rotate the adjustment control counterclockwise to position the shadow line at the bottom of the field of view.
6. Rotate the eyepiece to bring the crosshair into focus.
7. Move the shadow line to the crosshair reticle with the coarse adjustment control.
8. Rotate the dispersion correction wheel to eliminate any red or green color at the edge of the shadow line.
9. Turn the adjustment control to center the shadow line to the crosshair. The shadow line must be perfectly centered to obtain an accurate reading (Figure 3a).
10. Depress the READ button. The value of the test sample will be digitally indicated in the display window. Depressing the TEMP button will activate a temperature sensing device located in the measuring prism. The display will digitally indicate the actual temperature of the measuring prism and sample.

**NOTE:** If the LEICA ABBE MARK II Refractometer is to be used in an environment where foreign materials may contaminate the display buttons, use the splash guard accessory (Catalog No. 13104590). It will prevent instrument damage and possible instrument failure.

## 7.0 LEICA ABBE MARK II PLUS REFRACTOMETER

**NOTE:** Sections 7.0 through 7.4 are applicable only for the LEICA ABBE MARK II Plus Refractometer - Model 13104940/13104950.

The LEICA ABBE MARK II Refractometer provides computer interface, custom programming and printer interface capabilities.

The instrument becomes the LEICA ABBE MARK II Plus Refractometer (Model 13104940/13104950) with the addition of the RS232 port. The LEICA ABBE MARK II Plus can be ordered as:

- an original purchase
- a factory upgrade or
- a retrofit to an existing LEICA ABBE MARK II Refractometer.

### 7.1 Computer Interface

Turn on the LEICA ABBE MARK II Plus Refractometer with the switch located on the instrument's front panel. When the refractometer is first turned on, the software version (e.g., 1.6) will appear in the read display for a few seconds. The refractometer is ready to use when this message disappears.

By interfacing the LEICA ABBE MARK II Plus to a computer, sample reading results can be transferred to IBM-PC compatible data-handling software programs with communications

capabilities (e.g., Microsoft Works®, Framework II™ and ProComm®)¹. Software can be used to store, graph, analyze and print data. Leica also offers LEICA ABBE MARK II Plus Refractometer interface software (MARKPLUS) specifically designed for use with the instrument.

The MARKPLUS software allows the instrument user to reprogram the % Solids scales from the default "Brix" (% Solids) scale to a user defined scale.

### 7.2 Setting the baud rate

The baud rate of the LEICA ABBE MARK II Plus is set to 9600. To change the baud rate remove the cap on the left side of the instrument (Figure 3b). Turn the instrument on and using the eraser end of a pencil press the post switch which is centered in the hole. The display will show the current baud rate. To change the rate press the read button on the front of the instrument panel until the appropriate baud rate is displayed. Choose the baud rate from the following options: 300, 1200, 2400, 4800, 9600. Once the appropriate baud rate is displayed press the switch again. Turn off the instrument and replace the cap. The new baud rate will be permanently retained in the instrument memory.

### 7.3 Printer Interface

The LEICA ABBE MARK II Plus Refractometer can be interfaced to a serial port printer for printed data results using the 9-pin RS232 output. A parallel port printer may be used with a serial to parallel converter.

### 7.4 Custom Programming

LEICA ABBE MARK II Plus Refractometer interface software (MARKPLUS) may be configured for customized calibration curves and temperature compensation. (The standard LEICA ABBE MARK II Refractometer corrects to 20°C only.) However, MARKPLUS must be configured for each specific application or custom programming. Leica's MARKPLUS is designed to accept a Table of Refractive Index, Temperature and Concentration values and then generate custom scale data for instrument programming. Data is entered from a text file. Contact your Leica representative for more information.


Figure 3a


Figure 3b

## 8.0 REFRACTOMETER USE TECHNIQUES AND REFRACTION PRINCIPLES

Refractive Index measurements taken on the LEICA ABBE MARK II and LEICA ABBE MARK II Plus Refractometers are affected by temperature. Different values can be obtained when the same sample is measured at different temperatures. For example: as the temperature of a sugar solution sample rises, its density and refractive index decrease. This causes an apparent decrease in concentration when using the non-temperature compensated % Solids mode. The % Solids-TC (Temperature Compensated) mode automatically corrects the reading to the correct value at 20°C.

For every 1°C change in temperature, an aqueous solution containing sugar will change readings by approximately 0.08%. This is the average value for samples of 0 to 85% at 15°C to 45°C. As temperature rises the % Solids reading decreases. As temperature decreases, the % Solids reading increases.

The LEICA ABBE MARK II Refractometer temperature will rise approximately 1°C every 10 minutes for 1 hour until it stabilizes. An instrument that has been turned on for more than 1 hour in a 20°C room could read 26°C. The rise in temperature is principally due to heat generated by the power transformer inside the instrument housing.

The sample placed on the measuring prism surface of a LEICA ABBE MARK II Refractometer that has been turned on for more than an hour in a 20°C environment will assume the temperature of the prism (*i.e.*, approximately 26°C). Although the sample is now warmer and would indicate a lower concentration when measured in the uncompensated Brix mode, the BX-TC mode automatically compensates for the temperature difference and displays the concentration of the solution as if its temperature were 20°C. Without compensation for the 6° rise, the readout of a 70% solution would be 0.5% less (6°C x 0.08%) or 69.5%.

The measuring prism surface temperature can be lowered by running cooling fluid through the water jackets in the prism housing. However, the heat generated inside the instrument housing still affects the prism temperature. For this reason it is necessary to run fluid slightly below the desired temperature. A thermometer placed in the fluid line may measure several degrees below the actual displayed temperature of the measuring prism surface displayed when the temperature button is pressed. The opposite is true if the prism temperature has to be elevated above the ambient temperature.

The fluid circulation system allows the prism to be adjusted for temperatures up to 80°C. This may be necessary when using samples that solidify at lower temperatures.

Measuring viscous samples, such as high Brix sugar syrups, may be more difficult than measuring products with lower total dissolved solids.


Figure 4

Stir the sample on the measuring prism surface with a rubber spatula or an angled piece of clear plastic tubing. Completely cover the surface making sure all the air pockets between the sample and the measuring prism surface are removed.

For removing carbonation from carbonated soft drink samples, open and close the prism cover two or three times after loading. Doing this provides the best possible reading line and the best way to judge where the shadow line intersects the crosshair reticle.

**NOTE:** Remember the accuracy of the instrument is only as good as the clarity of the shadow line.

### 8.1 High Index/High Dispersion Samples

High refractive index samples with high dispersion properties may require additional use of Catalog #13K30650 589nm bandpass filter accessory. This accessory will insure that the shadowline will be sharp for these samples.

## 9.0 CALIBRATION WITH TEST GLASS

Refer to Chart 1, Refractive Index of Triple Distilled Water at Various Temperatures to compare the refractive index reading with the calibrated temperature reading.

To calibrate the LEICA ABBE MARK II Refractometer use:

- a refractive index test glass (Catalog No. 13104510)
- a bottle of 1- Bromonaphthaline (Catalog No. 13104520)
- a 1/16" (No. 4) Allen wrench (Catalog No. 13X30600)
- lens tissue

The refractive index test glass is imprinted with the refractive index scale. The bottom and front surfaces of the glass are polished. It is important that these surfaces are very clean when beginning calibration procedures.

**To calibrate:**

1. Turn the "Mode" selector switch to the refractive index position.
2. Apply a single drop of 1-Bromonaphthaline to the surface of the measuring prism surface.
3. Place the refractive index test glass on the measuring prism surface with the imprinted side facing up. Position the polished end of the test glass toward the illuminator - about 1/8" (2mm) from the edge (Figure 4).
4. Align the illuminator so that the light is directed onto the polished face of the test glass. Cover the illuminator lens with lens tissue to diffuse the light.
5. Adjust the illuminator arm, dispersion correction wheel, eyepiece and adjustment control so that the shadowline is centered, in focus, free of color. This is the same procedure used with a

Chart 1: Refractive Index of Triple Distilled Water at Various Temperatures

Temp °C	Refractive Index
18	1.3332
19	1.3331
20	1.3330
21	1.3329
22	1.3328
23	1.3327
24	1.3326
25	1.3325
26	1.3324
27	1.3323
28	1.3322
29	1.3321
30	1.3319

Source: Molar Refraction of Organic Compounds

liquid sample.

6. Depress the READ button and record the digital reading which appears in the display window.

If the reading differs from the value printed on the test glass by more than  $\pm 0.0001$ , repeat the measurement several times. Clean the measuring prism and test glass surfaces with alcohol, followed by water after each test.

If an error persists, make the following adjustment:

7. Turn the shadow line adjustment control slightly and depress the READ button. Repeat this procedure until the display agrees with the test glass value.

The shadowline will now be above or below the center of the crosshair reticle.

8. Insert the Allen wrench through the reticle access hole in the dispersion correction wheel. Adjust the reticle until the shadow line passes through the center of the crosshair (Figure 2 & Figure 3a).

The instrument is now calibrated.

## 10.0 ILLUMINATOR BULB REPLACEMENT

**DANGER:** The glass envelope of the bulb may be extremely hot. Do not attempt to change the bulb before it is completely cooled or without wearing protection.

**WARNING:** Unplug the refractometer power cord before attempting to replace the illuminator bulb.

The illuminator bulb can be replaced by removing the condenser sleeve (Figure 6). Always use the suggested replacement bulb (Catalog No. 13104620, pkg. of 3) for maximum efficiency and safety. After inserting a new bulb, the socket can be adjusted to center the light on the measuring prism with the condenser sleeve in place. If necessary, push down on the spring-loaded condenser sleeve lock pin and slide back and forth until the pin snaps into place.


Figure 5

## 10.1 Fuses

 **Warning:** Unplug the Refractometer Power Cord before attempting to check or replace fuses.

 **CAUTION:** Replacement of fuses is not user serviceable. Fuses must be replaced by qualified service personnel. The fuse shall be replaced only with the same size, type and rating.

The LEICA ABBE MARK II Refractometer electrical circuitry is protected by two (2) fuses. To replace the fuses turn the Adjustment Control Knob (Figure 1) two turns clockwise. This raises the elevator assembly to provide clear access to the fuses. Remove the adjustment control knob with a 1/8" Allen wrench. Remove the two Phillips head screws which hold the right panel cover to the bottom of the instrument housing. Slide the panel cover down as far as possible to facilitate its removal. With the right panel cover removed the fuses will be visible on the bottom plate of the instrument (Figure 5b). Replace the fuses and reassemble.

## 11.0 CARE AND CLEANING

The LEICA ABBE MARK II Refractometer is a precision optical instrument designed for years of trouble-free operation. The interior design has been carefully constructed for protection from dust or spills. The exterior is painted with a durable enamel. The stainless steel refractive prism assembly is sealed with corrosion-resistant epoxy.

The use of strong solvents such as N, N-Dimethylformamide or N, N-Dimethylacetamide, cresols, phenols and other tar acids or strong acetic solutions should be avoided. Materials such as Tetrahydrofuran or some lacquer thinners and their components can deteriorate the prism seal over a long period of time or when applied in elevated temperatures.

An accumulation of dust, solvents or sample, particularly on the prism assembly, may eventually penetrate the refractometer and affect accuracy or accelerate wear.

When the instrument is not in use, the measuring prism assembly should be kept closed with a folded sheet of lens tissue between the prisms. The prism surfaces should be cleaned with distilled water, methanol or an appropriate solvent immediately after each sample is used.

Warm distilled water is recommended for removing sugar residues and methanol for oil-based samples. The prism may be wiped with lens cleaning tissue, but it should not be wiped with a hard, dry cloth. After wiping the sample from the prism, rinse the prism again with distilled water.

Figure 5b

**NOTE:** Do not wipe the measuring prism surface dry and do not wipe it with anything if it is dry. Do not contact the prism or the surrounding seal with metal or with sharp or pointed objects.

The glass parts of the prism assembly are relatively hard and have excellent resistance to staining. Strong acids or bases should not be measured with this instrument, nor used as cleaning agents.

To clean the eyepiece, snap out the rubber insert at the end of the adjustable eyepiece tube. Clean the lens with a cotton swab or lens tissue lightly dampened with acetone or methanol. Replace the rubber insert.

## 12.0 SERVICING

The LEICA ABBE MARK II and LEICA ABBE MARK II Plus Refractometers should be serviced or repaired by factory-trained service and repair technicians. If operating difficulty is experienced contact Leica's Refractometer Customer Service Department at 716-686-3000 or fax 716-686-3085. Outside the United States contact your local Leica representative.

## 13.0 SPECIFICATIONS

**LEICA ABBE MARK II Refractometer Model 13104800/13104810**

**LEICA ABBE MARK II Plus Refractometer**

**Model 13104940/13104950**

<b>Display</b>	Direct reading LED
<b>Modes of Operation</b>	Refractive index ( $n_D$ ), %Solids (°Brix) and temperature compensated Brix - BX-TC
<b>Modes</b>	<b>Range</b>
$n_D$	1.3200-1.7000
<b>Brix</b>	0-95
<b>BX-TC at 20°</b>	0-95
<b>Temperature °C</b>	15°-50°C
<b>Power Requirements</b>	
Model 13104800	115V/60Hz
Model 13104810	230V/50Hz
Model 13104940	115V/50-60Hz
Model 13104950	230V/50-60Hz


Figure 6

<b>Electrical Approvals</b>	UL listed, CSA certified, VDE listed, CE Accuracy $\pm 0.0001n_D$ $\pm 0.1\%$ Solids ( $^{\circ}$ Brix)- 0.0 - 95% $\pm 0.1$ on custom scale
<b>Temperature Compensation</b>	Default mode corrects %Solids scale readings taken at 15 $^{\circ}$ C to 50 $^{\circ}$ C to a datum temperature of 20 $^{\circ}$ C; Mark II Plus Custom Mode corrects to user specified data
<b>Temperature Measurement</b>	Thermistor mounted in prism assembly
<b>Temperature Compensation Accuracy</b>	$\pm 0.1^{\circ}$ C
<b>Sample Types</b>	Transparent or translucent liquids and solids
<b>Dimensions</b>	6"W x 10.5"L x 14"H (152mm x 262mm x 356mm)
<b>Weight (net)</b>	16lbs (7.3kg)
<b>Sample Temperature Control</b>	Fittings on prism assembly for connecting temperature controlled water circulator
<b>Prism Assembly</b>	Stainless steel, glass sealed with solvent-resistant epoxy
<b>Optical Compensation</b>	Dual amici prism compensator for sharp shadow-lines and 589nm (Sodium D line) illumination

All models have passed RF Interference Testing, Class A

### 13.1 The following specifications refer only to the LEICA ABBE MARK II Plus Refractometer

<b>Display</b>	Direct digital readout of refractive index or concentration (% solids or custom scale) and temperature
----------------	--

#### Default Modes

1. Refractive Index 1.3200 to 1.7000 $n_D$
2. % Solids ( $^{\circ}$ Brix) 0.0 to 95.0%
3. % Solids ( $^{\circ}$ Brix) temperature compensated 0.0 to 95.0%  
Default % Solids ( $^{\circ}$ Brix) scales (0.0-85.0%) based 1974 I.C.U.M.S.A. Tables (% Sugar vs.  $n_D$  at 20 $^{\circ}$ C)

#### Custom Modes

- 1. 0.0 to 999.9 readout of concentration. User inputs concentration vs. refractive index data
- 2. 0.0 to 999.9 Temperature Compensated readout of concentration vs. refractive index vs. temperature data Custom Mode data are input using Leica software program (IBM PC compatible)

<b>Interface</b>	RS232 Serial
<b>Communication</b>	User selectable Baud rate, 300, 1200, 2400, 4800, 9600 bps - 8 data bits, one stop bit, no parity, Baud rate factory set at 9600. ASCII output with tab separators
<b>Connection</b>	DB9 Female
<b>Pin Assignments</b>	2 output, 3 input, 5 ground

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Framework II is a trademark of Microsoft Corporation.
ProComm is a registered trademark of Datastorm Technologies Inc.
- Escalas de % de sólidos por defecto (°Brix) (0.0-85.0%) basadas en Tablas I.C.U.M.S.A 1974 (% azúcar vs. nD a 20°C)
Framework II ist ein Warenzeichen der Microsoft Corporation
ProComm ist ein eingetragenes Warenzeichen von Datastorm Technologies Inc.
Framework II est une marque commerciale de Microsoft Corporation.