The system described herein is covered by one or more of the following U.S. Patents/ *Le système décrit aux présentes est protégé par un ou plusieurs des brevets américains suivants:* D367,109, D371,198, D428,150, 5,418,142, 5,515,170, 5,526,120, 5,563,031, 5,605,837, 5,780,304, 5,789,255, 5,922,530, 5,968,836, 6,045,567, 6,156,051, 6,197,040, 6,335,203 and/et 6,491,870. Use of the monitoring device described herein is protected under the U.S. patents 4,935,346/*L'utilisation du dispositif de surveillance décrit aux présentes est protégée par le brevet américain 4,935,346.*

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© 2003 LifeScan, Inc.
Milpitas, CA 95035

CA/OB/OTSS Rev. 08/2003

The OneTouch® SureStep® System is intended for use outside the body (in vitro diagnostic use). It should be used only for testing purposes and only with fresh whole blood samples. It should not be used for the diagnosis of diabetes.

Welcome to LifeScan

Congratulations! You've chosen a product made by LifeScan, a leading company in providing products to help you manage your diabetes. OneTouch® SureStep® blood glucose monitoring products are widely available without a prescription.

Blood glucose monitoring plays an important role in controlling your diabetes. A long-term study showed that keeping blood glucose levels close to normal can reduce the risk of complications by up to 60%. The results you get with the OneTouch® SureStep® System can help you and your Healthcare Professional monitor and adjust your treatment plan to gain better control of your diabetes.

Everything you need to know about using the OneTouch® SureStep® System is included in this owner's booklet and in the test strip package insert that comes in each box of test strips. It is important that you read both this booklet and the test strip package insert for critical information on how to use this product. If you have any questions, please feel free to call the toll-free number to speak with our Customer Care Line.

Remember to complete and mail the Warranty Registration Card.

Important phone numbers:

OneTouch® Customer Care Line:

CANADA 1 800 663-5521

www.OneTouch.ca

Meter Serial Number _____

Healthcare Professional _____

Pharmacist _____

Diabetes Educator _____

Other _____

CAUTION: Before using any product to test your blood glucose (blood sugar), read all instructions and practice the test. Do all quality control checks as directed and consult with a diabetes Healthcare Professional. These recommendations apply to all blood glucose monitoring systems and are supported by the Diabetes Educator Section of the Canadian Diabetes Association.

IMPORTANT: If you have followed the tips to ensure accuracy and your blood glucose results are still out of range, call your Healthcare Professional.

The OneTouch® SureStep® System may contain many small parts: test strip holder, battery door, batteries, lancets, control solution bottle, test strips, test strip bottle cap, etc. Keep the system out of the reach of small children because these parts may be dangerous if swallowed.

If using a 10% bleach solution to disinfect the meter, be sure to completely remove any bleach from the meter and test strip holder with a cloth dampened with water. Bleach will react with the test strips and cause inaccurate results. Failure to remove the excess bleach may result in an error message (E-5, E-6) or an inaccurately high result.

Units of Measure for Meter Readouts

A digital display showing the number 5.8 with a decimal point. Below the number, the unit mmol/L is displayed, and at the bottom, the date 9-17 is shown.

mmol/L

9-17

mmol/L
millimoles/liter
decimal point

A digital display showing the number 107. Below the number, the unit mg/dL is displayed, and at the bottom, the date 9-17 is shown.

mg/dL

9-17

mg/dL
milligrams/deciliter
no decimal point

Your OneTouch® SureStep® Meter can display results in two different units of measure, millimoles per liter (mmol/L) or milligrams per deciliter (mg/dL). The mmol/L unit is standard in Canada. The mg/dL unit is standard in the United States. Please contact your Healthcare Professional to determine the unit of measure you should use. If the incorrect unit of measure is used, you may get inaccurate results, which can lead to improper treatment.

Table of Contents

Learning		2
	Purpose	2
	What to Expect from Your System	3
	Warnings and Precautions	4
	Importance of Following Instructions	5
	OneTouch® SureStep® Test Strips	8
	Optional Colour Chart	9
	Display Symbols	10
Coding		13
	Coding the Meter	13
Checking		15
	Checking the System	15
	Understanding Control Solution Results	20
Testing		22
	Testing Your Blood	22
	Testing Tips for Accurate Results	23
	Test Procedures	24

OneTouch® Customer Care Line:

CANADA 1 800 663-5521

www.OneTouch.ca

Results		35
	Understanding Your Results	35
	How to Check if Your Meter is Working Correctly	37
	Comparing Meter and Laboratory Results	38
	Health Conditions Which Can Cause Out-of-Range Results	40
Maintenance		41
	Caring for Your Meter	41
	General Care	41
	Cleaning the Meter	41
	Changing the Batteries	44
	Storage Instructions	47
Troubleshooting		48
	Error Messages	48
Symbols		56
	Symbols	56
Settings		58
	Changing Meter Settings	58
	Entering Set-Up Mode	59
	Changing the Date and Time	59
	Setting the Time	60

OneTouch® Customer Care Line:

CANADA 1 800 663-5521

www.OneTouch.ca

Setting the Year	60
Setting the Month	61
Setting the Day	61
Display of Last Results	61
Unit of Measurement	61
Beep Signals	62
14- and 30- Day Averages	63
Time Format	64
Date Format	64
Exiting the Set-Up Mode	64
<hr/>	
Memory	65
Background Information on Meter Memory	65
Getting All Results from Meter Memory	65
First Time Turning on Meter	67
Uploading Memory Data	67
<hr/>	
Specifications	68
Specifications	68
Guarantee and Warranty	69
<hr/>	
Quick Reference	70

OneTouch® Customer Care Line:

CANADA 1 800 663-5521

www.OneTouch.ca

ONETOUCH[®]
SureStep[®]

Purpose

People with diabetes use the OneTouch® SureStep® Blood Glucose Monitoring System to measure the level of glucose (sugar) in their blood. The system is intended for testing whole blood samples at the time when the blood sample is obtained.

Why Test Your Blood Glucose

Diabetes is a disease where the body has trouble controlling glucose levels properly. Glucose levels may change greatly and quickly. Glucose levels that are too high or too low can lead to serious health issues. This test enables a person to know what his or her glucose is at the time it is measured.

In 1993 the National Institutes of Health finished the Diabetes Control and Complications Trial (DCCT.) This was a long-term study of people with type 1 diabetes. This

study found that you should keep your blood glucose close to the levels of people without diabetes. You can reduce the risk of problems with the eyes, kidneys, and nervous system by about 60%. You should be aware that keeping your blood glucose level close to the levels of people without diabetes increases your risk of hypoglycemic events¹ (low blood sugar).

How often you need to test your blood glucose depends on many factors. These include your age, if you have type 1 or type 2 diabetes, your medications, and your health. Your Healthcare Professional will guide you. It is important that you make testing part of your daily routine.

When Not to Use

Do not use OneTouch® SureStep® blood glucose results to determine whether or

not you have diabetes. Only your doctor can diagnose this disease.

What to Expect from your System

- ◆ Blood glucose meters are clinically accurate, but not a substitute for a laboratory result. Originating from laboratory methods, glucose meter systems provide a convenient, portable way to test inexpensively and quickly. However, meter results can be affected by factors and conditions that do not affect laboratory results in the same way. Keep in mind that testing conditions for laboratories are tightly controlled and highly automated. Your meter glucose result should agree with a laboratory result to within $\pm 20\%$ under normal conditions. A result within that range is considered clinically accurate when testing with the OneTouch® SureStep® Meter. However, results can differ by more than $\pm 20\%$ in some situations. See **Comparing Meter and Laboratory Results** for more information.
- ◆ Blood glucose meters increase your ability to frequently check your blood sugar at reasonable cost so you can see the effect different foods, exercise, and medications can make on your glucose level. You can keep your blood glucose in better control with these devices.
- ◆ Make sure that you are familiar with the limitations of the system before testing. See important information in your OneTouch® SureStep® Test Strip Package Insert. Know the factors that can affect the accuracy of your results.
- ◆ Make sure that you understand your results before acting on them, especially if your results do not reflect how you

feel. See **Understanding Your Results** on page 35.

Warnings and Precautions

- ◆ Read all the instructions, including the test strip package insert, and practice the test with control solution (see **Checking the System**) before testing your blood glucose. Do all quality control checks as directed, including checking the meter display symbols (see **Description of System**), and consult your Healthcare Professional before using this product.
- ◆ The system may contain small parts such as lancets, bottle caps, etc. Keep your system out of the reach of small children because these parts may be dangerous if swallowed. The test strip packaging contains a drying agent that may be harmful if inhaled or swallowed. It may also cause skin or eye irritation.
- ◆ Do not use the test strips after the expiration date printed on the bottle label or more than 4 months after first opening the bottle, whichever comes first. You may get inaccurate results if you use the test strips after these dates.
- ◆ Do not use strips that are bent, cut, or altered in any way.
- ◆ Do not use the control solution after the expiration date printed on the bottle label or more than 3 months after first opening the bottle, whichever comes first. Control solution results may not be reliable if you use the solution after these dates.
- ◆ Testing at a very high or very low temperature or relative humidity may affect results particularly if both the temperature and humidity are either high or low. If you test under these conditions, repeat

the test in a climate-controlled location as soon as possible.

- ◆ Excessive water loss (severe dehydration) may cause false low results.^{2,3} Severe dehydration, if untreated, can lead to life-threatening medical complications such as diabetic ketoacidosis. Some factors causing dehydration are vomiting, diarrhea, inadequate fluid intake, some drugs (for example, diuretics), sustained uncontrolled diabetes, and shock.⁴ If you have these symptoms, consult your Healthcare Professional immediately, even if your glucose reading is not unusually high.
- ◆ Because some diseases can be transmitted by blood, use care when disposing of used lancets and test strips.

NOTE: *Blood glucose monitoring requires the help of a Healthcare Professional. Work with your Healthcare Professional to set your own target blood glucose values, arrange your testing times, and discuss the meaning of your blood glucose results.*

Importance of Following Instructions

Your OneTouch® SureStep® Blood Glucose Monitoring System will provide accurate glucose results that can help you control your diabetes, if you follow the directions in this Owner's Booklet and the test strip package insert. Be sure to check the test strip package insert each time you open a new box of test strips for important updates. If you do not follow the instructions, you may get inaccurate results, which can lead to improper treatment.

The Complete OneTouch® SureStep® Blood Glucose Monitoring System

- OneTouch® SureStep® Meter (Batteries installed)
- Owner's Booklet
- Quick Reference Guide
- OneTouch® SureStep® Control Solution
- OneTouch® Penlet® Plus Adjustable Blood Sampler
- OneTouch® FinePoint™ Lancets
- Carrying Case
- Logbook
- OneTouch® SureStep® Test Strips (May be sold separately)

OneTouch® SureStep® Blood Glucose Meter

DATA PORT

Allows you to transfer the information stored in the meter's memory to a computer to view, analyze, and print.

DISPLAY

Shows your test results and symbols that guide you through the test.


LENS AREA

Contains the meter optics that reads the glucose level. Test strip holder covers the lens area.

SERIAL NUMBER (Example)

BATTERY COMPARTMENT

Holds two AAA alkaline batteries.

C BUTTON (Green)

Codes the meter, changes the meter settings, marks control solution test results, and scrolls through the memory.


POWER BUTTON (Blue)

Turns the meter on and off, accesses the memory, and navigates through the set-up mode.


CONTACT POINTS

Sense the position and orientation of the OneTouch® SureStep® Test Strip.

TEST STRIP HOLDER

The cover protects the lens area and the base holds the OneTouch® SureStep® Test Strip.

OneTouch® SureStep® Test Strips


When blood is applied to the pink test square, it is absorbed and a chemical reaction takes place. A blue colour forms in the confirmation dot on the back of the test

White Tip is inserted into the meter. The side where you apply blood must be face up when you insert the test strip.

Pink Test Square is where you apply a drop of blood.

White Pad absorbs the excess blood that may extend beyond the pink test square.

Do not apply blood directly onto the white pad.

strip and is read by the meter to determine the blood glucose level; the darker the colour, the higher the blood glucose level.

Black Tip helps the meter detect that the strip has been inserted properly.

Confirmation Dot on the back of the test strip is where you check to confirm proper blood application. When it turns completely blue, enough blood has been applied for an accurate test.

Optional Colour Chart

Directly after reading a meter result, you can compare the colour of the confirmation dot on the back of the test strip with the Colour Chart on the test strip bottle label. The chart shows two shades of blue, which correspond to a very low blood glucose (2.8 mmol/L) level and a very high blood glucose (19.4 mmol/L) level. It also shows an off-white shade for new, unused test strips.


Important: *The sample colour dots are only approximate. They can give you an indication of whether your glucose is very low (2.8 mmol/L) or very high (19.4 mmol/L). Results obtained using the sample colour dots are not a replacement for a OneTouch® SureStep® meter test result. Do not use Colour Chart to verify control solution tests.*

Display Symbols


The OneTouch® SureStep® Blood Glucose Meter uses symbols and error messages to guide you through a test and to help you determine if something is wrong.

Carefully check that all symbols shown on the right appear each time you turn on your meter. To check that all display symbols are working, press and hold down the blue power button. This display check also appears briefly each time you turn on the meter. Call the OneTouch® Customer Care Line at 1 800 663-5521 if any segments are missing.

Display Check


When you turn on the meter, it automatically recalls and displays your most recent glucose result and the date it was measured.

(Example)


The current code number is displayed.

The number on the display must match the code number on the OneTouch® SureStep® Test Strip vial.


(Example)


Apply blood to the centre of the pink test square on the test strip.

Check that the pink side of the test strip faces up and firmly push the test strip into the test strip holder until it stops. Insert the test strip tapered-end first.

The meter is measuring your glucose level.

The meter displays your test result in mmol/L or mg/dL together with the date and time of the test.

(Example)


(Example)

When the meter determines that your blood glucose level exceeds 27.8 mmol/L (500 mg/dL), it displays *HI* in virtually all cases.


(Example)


(Example)

Clean the meter and test strip holder.

Replace the two AAA alkaline batteries.

Coding the Meter


A code number is printed on each test strip vial label. The code numbers (ranging from 1 to 21) are used to calibrate the meter.

Each time you test, check that the code number on the meter display matches the code number on the test strip vial label.


CAUTION: *If the code number on the meter display does NOT match the code number on the test strip vial label, you may get inaccurate results.*

1. Press the blue power button to turn on the meter.

Make sure there is not a test strip in the meter when you turn it on.

Display Check


(Example)


The display check appears briefly. Carefully check that all symbols appear. The result of your last blood glucose or control solution test appears.

If you are turning on the meter for the first time, --- appears on the display.

- ◆ Next, a code number appears for a few seconds.


2. Compare the code number on the test strip vial with the code number on the meter display.

(Example)

If the code numbers match, you may begin testing.

3. If the code numbers do not match, press the green **C** button. Continue pressing until the code number on the display matches the code number on test strip vial. The meter will retain this code until you change it.

NOTE: If the meter advances past the code number before you are able to change it, continue pressing the green button until the code number appears.

CAUTION: Your meter has been preset to meet the needs of most users in the Canada. Refer to **Changing Meter Settings** if you need to verify or change the settings, such as date, time, or unit of measurement. See page 58 for a complete listing of features that can be changed.

Checking the System

Use OneTouch® SureStep® Glucose Control Solution to check if the meter and test strips are working together properly. A control solution test is performed exactly like a blood test except that control solution is used as a sample instead of your blood. This control solution is a liquid with a known amount of glucose in it. Practice the test with the control solution before testing your blood for the first time.

When to do a control solution test.

- ◆ Before testing your blood for the first time.
- ◆ Every time you open a new vial of test strips.
- ◆ If you drop the meter.

- ◆ To check that the meter and test strips are working together properly.
- ◆ If your results do not reflect how you feel.
- ◆ After you clean the meter.



OneTouch®
SureStep® Glucose
Control Solution

WARNING:

- ◆ *Do not use the control solution to clean the meter.*
- ◆ *Do not swallow or inject the control solution or put it in your eyes.*
- ◆ *The control solution is for use outside the body (in vitro use).*
- ◆ *The control solution contains a dye that stains clothing.*


1. Press the blue power button to turn on the meter. The display check appears briefly. Carefully check that all symbols appear. Make sure a test strip is not in the meter when you turn it on.
2. Compare the code number on the meter display with the code number on the test strip vial. If they do not match, press the green **C** button until the codes match.

3. Remove a test strip from the vial and replace the vial cap. Close it tightly.
 - ◆ Check the expiration date and discard date. Do not use the test strips after these dates.
 - ◆ Compare the colour of the confirmation dot on the test strip with the “Unused” colour dot shown on the Colour Chart on the test strip vial label. If the colour of the confirmation dot on the test strip is darker than that shown on the Colour Chart, do not use the test strip. The result may be inaccurate.

NOTE: *If opening the test strip vial for the first time, write the discard date (4 months from today) on the vial label and perform a control solution test following the instructions below.*

4. Shake the control solution vial. Apply 1 drop of control solution to the centre of the pink test square.


- ◆ The first time you open the control solution vial, write the discard date on the label in the space provided. The discard date for control solution is 3

months from the date you first open the vial.


- ◆ Check expiration date and discard date. Do not use the control solution after these dates.
- ◆ If the entire white pad surrounding the pink test square is soaked, you have applied too much Glucose Control Solution.

CAUTION: *Applying too much control solution may cause out-of-range results and may make the meter dirty.*

5. Turn the test strip over and check to see that the confirmation dot is completely blue. If it is not, discard the test strip. Repeat the test with a new test strip, applying a larger drop of control solution to the centre of the pink test square.


6. Firmly insert the test strip, confirmation dot side down, into the test strip holder until it stops.

CAUTION: *If you do not firmly insert the test strip into the meter until it stops, you may get an inaccurately low result.*

IMPORTANT: *Immediately insert the test strip after applying the control solution. If you do not insert the test strip within 2 minutes, discard the test strip and start over with a new test strip.*

A flashing clock symbol appears while the meter is reading your result. The control solution result appears in about 30 seconds. The wait may be shorter or longer depending upon the current temperature.

- Press the green **C** button to mark the test result as a control solution result (C E L) in meter memory. Marking the control result prevents it from being included in the 14- and 30-day blood glucose averages.

(Example)

NOTE: Pressing the green **C** button again while the glucose level is displayed removes the C E L marker.


- Compare the result with the control range printed on the test strip vial label.

OneTouch®
SureStep®
Control
Range


(Example)

CAUTION: The range printed on the test strip vial label is the expected range for control solution results. It is **NOT** a recommended range for your blood glucose level.

(Example)

- Record the control solution test result in your logbook. Place a C next to the result to mark it as a control solution test.

Understanding Control Solution Results

Control solution test results should fall within the printed range on the test strip vial label. If a control solution test result falls outside the range, repeat the control test. Control results that fall outside the range may indicate that:

- ◆ you are not following all test steps or performing them incorrectly.
- ◆ the code number on the meter does not match the code number on the test strip vial.
- ◆ the meter or test strip holder is dirty.
- ◆ the control solution is expired, or beyond its discard date, contaminated, or not well mixed.
- ◆ there is a problem with the test strip (expired, damaged, etc.).
- ◆ the meter has malfunctioned.
- ◆ the control solution test is performed outside of the temperature range 15–35°C (59–95°F).
- ◆ you have applied too much control solution.

If your control solution result falls outside the printed range, repeat the control solution test with a new test strip. If your second control solution result is also outside the printed range, call the OneTouch® Customer Care Line at 1 800 663-5521.


WARNING: *Do NOT use your system to test your blood when your control solution results fall outside the printed range. You may get inaccurate blood glucose results. Call the OneTouch® Customer Care Line at 1 800 663–5521 if this occurs.*

Testing your Blood

Choose a clean, dry work surface. Make sure you have all the items needed for a blood test as shown.


Meter

Sterile Lancet

OneTouch®
SureStep®
Test Strip

Lancing Device

Depth Indicator

Testing Tips for Accurate Results

There are several important steps you should follow before and during the test that will help you get accurate results.

Before Testing

- ◆ Wash and dry your hands.
- ◆ Use a new, sterile lancet.
- ◆ Check the test strip expiration date and discard date (4 months after first opening the test strip vial). Do not use the test strips after these dates.
- ◆ Compare the colour of the confirmation dot on the **Unused** test strip with the colour shown for an “Unused” colour dot on the Colour Chart on the test strip vial label. If the colour of the confirmation dot on the test strip is darker than the “Unused” colour, do not use the test strip.

- ◆ Make sure the code on the meter display matches the code on the test strip vial label.
- ◆ Do a test with OneTouch® SureStep® Glucose Control Solution to make sure the meter and test strip are working properly together.
- ◆ Make sure the test strip holder, lens area, and contact points on the meter are clean. See **Cleaning the Meter** for cleaning instructions.

During Testing

- ◆ Make sure you applied enough blood to the test strip by checking that the confirmation dot turned completely blue.
- ◆ Make sure you did not apply too much blood to the test strip. You have applied too much blood if the white pad is soaked and there is no white visible.

Test Procedures



1. Remove the lancing device cap. Insert a lancet into the lancet holder. Push down firmly until the lancet is seated in the lancet holder.

CAUTION: Only one person should use the lancing device.


CAUTION: Never share a lancet with another person. Lancets are for single use only. Always use a new, sterile lancet.


2. Carefully remove the lancet's protective cover. Replace the lancing device's cap.


3. If you are using an adjustable sampler, adjust the puncture depth setting as needed.


4. Cock the lancing device.

5. Use warm, soapy water to wash your hands. Dry them thoroughly. To increase blood flow to the fingertips, massage your hand two or three times from wrist to fingertips.
6. Press the blue power button to turn on the meter. The display check appears briefly. Carefully check that all symbols appear. Your last test result appears next. (If this is your first test, --- appears instead.)


(Example)

Make sure a test strip is not in the meter when you turn it on.


Next, the code number appears.


- Compare the code number on the meter display with the code number on the test strip vial. If they do not match, press the green **C** button until the codes match.

These two symbols will alternate on the display. The meter is now ready to perform a test. The meter will turn off automatically after 2 minutes if a test strip is not inserted.

Apply blood to
test strip

Insert test strip
into meter

8. Remove a test strip from the vial and replace the vial cap. Close it tightly.
- ◆ Check the expiration date and discard date. Do not use the test strips after these dates.


- ◆ Compare the colour of the confirmation dot on the test strip with the “Unused” colour dot shown on the Colour Chart on the test strip vial label. If the colour of

the confirmation dot on the test strip is darker than that shown on the Colour Chart, do not use the test strip. The result may be inaccurate.


CAUTION: Use a test strip as soon as possible after removing it from the vial. Prolonged exposure to moisture, light, or heat can cause inaccurate results. Store test strips tightly capped in the original vial in a cool, dry place away from direct heat and sunlight. Never transfer test strips to another container.

NOTE: If opening the test strip vial for the first time, write the discard date (4 months from day of opening) on the vial label and perform a control solution test following the instructions in **Checking the System.**


9. Hold the lancing device firmly against the side of your finger. Firmer pressure will give a deeper puncture. Use a different spot each time to help prevent soreness and callus formation. Press the release button. If needed, gently massage your finger from base to tip to help obtain the needed volume of blood.

10. Apply one drop of blood to the centre of the pink test square in one of the ways shown. **Never add a second drop as this may cause inaccurate results.**


- ◆ Turn the test strip so the pink test square faces your finger. Touch the centre of the pink test square on the test strip to the drop of blood on your finger. The pink test square will quickly absorb blood.


OR

- ◆ Hold the test strip under your finger and touch the drop of blood to the centre of the pink test square.


Be careful not to apply too much blood, especially if you tend to bleed easily or are taking blood-thinning medications.


If the entire white pad surrounding the pink test square is soaked, you have applied too much blood. Your test result may be inaccurate. You must start over with a new test strip.

11. Check the confirmation dot on the back of the test strip after a few seconds. If the confirmation dot is completely blue, go to step 12.


If white patches or streaks appear on the confirmation dot, you have not applied enough blood to obtain an accurate test. You must start over with a new test strip.

If there is no colour or if white patches or streaks appear after you have repeated the test, call the OneTouch® Customer Care Line at 1 800 663-5521.


12. Confirm that the meter is still on. Firmly insert the test strip, tapered-end first, confirmation dot side down, into the test strip holder until it stops.

CAUTION: *If you do not firmly insert the test strip into the meter until it stops, the test may start but you may get an inaccurately low result. Do not remove test strip until countdown is complete.*

IMPORTANT: *Immediately insert the test strip after applying the blood. If you do not insert the test strip within 2 minutes, discard the test strip and start over with a new test strip.*

A flashing clock symbol appears while the meter is reading your result. Your glucose result appears in about 30 seconds. Do not remove strip until countdown is complete. The wait may be shorter or longer depending upon the current temperature.


The test result is displayed alternating with the time, then with the date. The meter displays results between 0-27.8 mmol/L (0 and 500 mg/dL). If the result is above 27.8 mmol/L (500 mg/dL) the meter displays *HI* in virtually all cases.

(Examples)

IMPORTANT: *When proper procedures are followed, the meter will display HI in virtually all cases when your result is above 27.8 mmol/L (500 mg/dL). The confirmation dot colour will also be darker than the 19.4 mmol/L (350 mg/dL) sample colour dot on the Colour Chart. This indicates very high blood glucose (severe hyperglycemia). Follow the advice of your Healthcare Professional.*

If you get an error message, see “What to do” in **Troubleshooting**.

(Example)

If you get an error message instead of a number as a result, you can quickly check to see if your glucose is very low (2.8 mmol/L) or very high (19.4 mmol/L) prior to troubleshooting the error message. To quickly see if your glucose is very low (2.8 mmol/L) or very high (19.4 mmol/L), compare the colour of the test strip confirmation dot with the sample colour dots on the Colour Chart on the test strip vial label.

- ◆ If the test strip confirmation dot colour is approximately the same as or lighter than the “2.8 mmol/L” colour dot, your glucose level is very low.
- ◆ If the test strip confirmation dot colour is approximately the same as or darker

than the “19.4 mmol/L” sample colour dot, your glucose level is very high.


CAUTION: *The Sample colour dots are only approximate. They should be used only as a guide in the event you get an error message instead of a numerical result. They can give you an indication of whether your glucose is very low 2.8 mmol/L (50 mg/dL) or very high 19.4 mmol/L (350 mg/dL) while you troubleshoot the error message. Results obtained using the sample colour dots are not a replacement for a OneTouch® SureStep® meter test result. Do not use Colour Chart to verify control solution tests.*

- Record your test result in your logbook.

If a test result falls outside of the target range that you and your Healthcare Professional have established and you


cannot link it to diet, exercise, medication, or stress, repeat the test.

- Remove and discard the test strip.
- Turn off the meter by pressing the blue power button. The meter automatically turns off 3 minutes after it displays your test result.
- Remove the lancing device cap.

Removing Lancing Device Cap


17. Place the protective cover on a hard surface and push the lancet tip into it.

CAUTION: To reduce the risk of an accidental needle stick, carefully follow the lancet and lancing device instructions for safe removal and disposal of used lancets.

19. Return the lancing device to its original position, following the lancing device instructions.

18. Point the lancet down and away from you. Eject the used lancet into a container for sharp objects.

Understanding Your Results

The result is displayed on the meter as a plasma glucose value. The system's plasma calibration makes it easier for you to compare meter results with lab results. See **Comparing Meter and Laboratory Results** for more information.

You and your Healthcare Professional have established target ranges for your blood glucose levels and any actions you should take in response to your test results. Follow your Healthcare Professional's recommendations to receive the maximum benefit from testing.

If a test result falls outside of the target range that you and your Healthcare Professional have established and you cannot link it to diet, exercise, medication, or stress, repeat the test. If the result is similar and

still outside of your range, review **Testing Tips for Accurate Results**.

If you get an error message instead of a number as a test result, compare the test strip confirmation dot to the sample colour dots on the Colour Chart on the test strip vial label. If the comparison indicates your glucose is very low or very high and you have symptoms, follow the advice of your Healthcare Professional. If you do not have symptoms, follow the "What to Do" recommendations in **Troubleshooting** to resolve the problem. Repeat the test with a new test strip using the procedures outlined in **Troubleshooting**. If you continue to get an error or other message, call the OneTouch® Customer Care Line 1 800 663-5521.

Low or High Results

You have **low blood sugar** (hypoglycemia) if:

- ◆ Your meter result is 3.3 mmol/L (60 mg/dL) or below.
- ◆ The colour dot comparison shows approximately 2.8 mmol/L (50 mg/dL) (very low).

You have **high blood sugar** (hyperglycemia) if:

- ◆ Your meter result is 10 mmol/L (180 mg/dL) or greater.
- ◆ The colour dot comparison shows approximately 19.4 mmol/L (350 mg/dL) (very high).

If you have symptoms that match your test results, or if you continue to get results that fall below 3.3 mmol/L or above 10 mmol/L when you repeat the test, follow the advice of your Healthcare Professional.

CAUTION:

- ◆ *If your meter result is very low, 1.1 mmol/L (20 mg/dL) or less, and you **do not have symptoms** of low blood glucose, remove the test strip from the meter and look at the confirmation dot on the back of the test strip. You should do one of the following depending on the colour of the confirmation dot on the test strip:*
 - *If the confirmation dot is not completely blue or has white patches or white streaks (see “**How To Test Your Blood**”), repeat the test with a new test strip. If the problem continues, call the OneTouch® Customer Care Line at 1 800 663-5521.*
 - *If the confirmation dot is completely blue, compare the confirmation dot to the sample colour dots on the Colour Chart of the test strip label on the vial. If this comparison indicates that your*

blood glucose is very low (2.8 mmol/L), conduct a Control Solution Test (see **“Checking The System”**). If the Control Solution test indicates that your meter and test strips are working together properly, repeat your blood glucose test with a new test strip. If you continue to receive a very low blood glucose test result, the low glucose result is likely to be correct and you should follow your Healthcare Professional’s advice for treating low glucose levels. If the Control Solution test is out of range, call the OneTouch® Customer Care Line at 1 800 663-5521.

IMPORTANT: When proper procedures are followed, the meter will display **“HI”** in virtually all cases when your result is above 27.8 mmol/L (500 mg/dL). The confirmation dot colour will also be darker than the

“19.4 mmol/L” sample colour dot on the test strip vial label. This indicates very high blood glucose (severe hyperglycemia). Follow the advice of your Healthcare Professional. If you have never experienced blood glucose levels this high or have not been instructed in how to deal with them by your Healthcare Professional, you should contact your Healthcare Professional at once.

How to Check if Your Meter is Working Correctly

The easiest way to check if your meter and test strip are working together properly is to run a test with OneTouch® SureStep® Glucose Control Solution. See **Checking the System** for more information.

Comparing Meter and Laboratory Results

The OneTouch® SureStep® System is plasma calibrated. Because most laboratories measure plasma glucose, you can directly compare your meter test results with the lab results. However, your meter blood glucose test result may differ somewhat from your laboratory result due to normal sources of variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. Your meter glucose result should agree with a laboratory result within $\pm 20\%$ most of the time under normal conditions. A result within that range is considered accurate when testing with the OneTouch® SureStep® Meter. However, results can differ by more than $\pm 20\%$ in some situations. Such results are more likely to occur under the following conditions:

- ◆ If you have recently eaten, exercised, taken medication, or experienced stress.⁵ If you have eaten recently, the blood glucose level from a fingerstick can be up to 3.9 mmol/L (70 mg/dL) higher than that of a venous sample (blood drawn from a vein) used for a lab test.⁶
- ◆ For the most accurate comparison between meter and laboratory results, the meter test must be done within 10 to 15 minutes of the laboratory test. If samples are taken more than 15 minutes apart, the two tests could have very different results because blood glucose levels can change significantly within short periods of time.

- ◆ Factors such as a low or high hematocrit (the percentage volume of red blood cells in your blood). For example, if your hematocrit is too low (below 25%) or too high (above 60%).
- ◆ Severe dehydration (loss of body fluid) may also cause meter results to be different from a laboratory result. Repeated vomiting or diarrhea, both of which can occur when your glucose is very high, can cause this. Your blood may become thick, with not enough liquid to properly wet the test strip. This may cause your meter to give false low blood glucose results.
- ◆ If you test at a very low or high temperature or a low or high relative humidity or some other extreme environmental condition. If you test under these conditions, repeat the meter test in a climate-controlled place as soon as possible. The operating temperature is 10–35°C (50–95°F). The operating humidity is 10–90%. You should always test within these ranges.
- ◆ You may also get inaccurate results if you test when both the temperature and relative humidity are near the limits defined above. If you test under these conditions, repeat the meter test in a climate-controlled place as soon as possible.
- ◆ Conditions affecting circulation to the hands may also lead to inaccurate results. These include shock, some drugs called vasoactive agents, and Raynaud's disease.^{4, 7, 8}

See the **Important Technical and Healthcare Professional Information** section in the OneTouch® SureStep® Test Strip package insert for typical accuracy and precision data for the OneTouch® SureStep® system. Make

sure you read **Important things to know before testing** in the test strip package insert. You should also review “Testing Tips for Accurate Results” in the **Testing Your Blood** section.

Health Conditions Which Can Cause Out-of-Range Results

If you are getting results outside your target range and you believe the OneTouch® SureStep® Meter and OneTouch® SureStep® Test Strips are working together properly, ask yourself the following questions.

- ◆ “How do I feel?”
- ◆ “Do I have symptoms of low blood sugar (hypoglycemia): slurred speech, headache, tingling of lips, cool sweating, rapid heartbeat, disorientation, weakness,

dizziness, fainting, hunger, nervousness, irritability, or tremors?”

- ◆ “Do I have symptoms of high blood glucose (hyperglycemia): fatigue, dry mouth and skin, increased thirst or hunger, blurry vision, increased urination, rapid deep breathing, or a fruity smell to the breath?”
- ◆ “What have I eaten?”
- ◆ “How much have I exercised?”
- ◆ “Am I ill (common cold, flu, etc.)?”
- ◆ “Have I been under stress lately?”

If your Healthcare Professional has not already instructed you in dealing with these situations, you should contact your Healthcare Professional.

Caring for your Meter

General Care

- ◆ Keep your meter dry and avoid exposing it to extremes in temperature. Do not store your meter or test strips in your car or in the refrigerator or freezer.
- ◆ Dropping the meter can damage the electronics. If you drop your meter, do a control solution test (see **Checking the System**) and a display segment check (see **Display Symbols**) before using the meter to test your blood.
- ◆ Do not take the meter apart. Sensitive parts could be damaged causing inaccurate results. Taking the meter apart will also void your warranty.

Cleaning the Meter

CAUTION:


- ◆ Use water to clean the meter.
- ◆ Do not get water inside the meter case.
- ◆ Never immerse the meter or hold it under running water.
- ◆ **Do not** use alcohol, glass cleaners, or any other cleansers to clean the test strip holder or lens area.

NOTE: *For Healthcare Professionals or caregivers who assist multiple users, to **disinfect** the meter, cleanse the meter as described below using a 10% bleach solution (9 parts water, 1 part bleach) and then clean again with water to remove any remaining bleach. Thoroughly dry the meter. Failure to remove the excess bleach may result in an error message (Er5, Er6) or an inaccurately high result.*





Always clean the meter when this symbol appears.


1. Press down at the top end (near blue power button) of the test strip holder and slide it away from the meter.


2. Wipe the test strip holder gently with a cotton swab or soft cloth dampened with water. Thoroughly wipe the grey area on the inside cover. Clean both sides of the test strip holder base, especially around the hole. **Dry the test strip holder completely with a soft cloth or lint-free tissue.**

3. Dampen a cotton swab or soft cloth with water and wipe the lens area and contact points with it. Do not scratch the lens area. **Dry with a soft cloth or tissue that does not contain lint, lotion, or perfume. Do not use abrasive materials such as paper towels.**

4. Slide the test strip holder into the meter, press down, and push forward at the test strip insertion point until you hear a click. Make sure the test strip holder is firmly in place and does not extend beyond the meter.

5. Turn on the meter by pressing the blue power button. The display check appears briefly. Carefully check that all symbols appear.

6. If the *Er 5* symbol appears after cleaning, check to see if you inserted the test strip holder completely and securely (see **step 4** in this section).


Note: *Clean the outside of the meter by gently wiping it with a cloth dampened with water. Dry it with a soft cloth.*

Note: *If you disinfected the meter with a 10% bleach solution, you may not have completely removed the bleach residue from the grey area inside the test strip holder cover. Wipe this area with a cloth dampened with water and dry thoroughly. Failure to remove the excess bleach may result in an error message (Er-5, Er-6) or an inaccurately high result.*

Changing the Batteries

CAUTION: The OneTouch® SureStep® Meter uses two AAA alkaline batteries. **Do not use rechargeable batteries.**

1. When this symbol appears constant on the display, the batteries are low and need to be replaced as soon as possible with new AAA alkaline batteries. If this symbol starts flashing, battery power is insufficient to run the meter. You will not be able to test.

Note: Low batteries and battery removal will not affect test results stored in meter memory. However, the current time and date will be lost and must be reset each time you change the battery. See **Changing Meter Settings** for information on setting the current date and time.

2. Turn off the meter and remove the battery door by squeezing the tab and lifting it up.


3. Remove the old batteries with your fingers, or by tapping the meter gently against the palm of your hand. Never tap the meter against a hard surface.

4. Insert two new AAA alkaline batteries, matching the + end of each battery with the + sign in the battery compartment.

5. Insert the two battery cover hinges into the slots and press down gently until the tab catches.

6. Set the current time and date in the meter as described in **Changing Meter Settings**.

Storage Instructions

WARNING: *The system contains small parts such as lancets, vial caps, etc. Keep your system out of the reach of small children because these parts may be dangerous if swallowed. The test strip packaging contains a drying agent that may be harmful if inhaled or swallowed. It may also cause skin or eye irritation.*

- ◆ Store the OneTouch® SureStep® Meter at room temperature. If the meter is exposed to cold, direct heat, or sunlight (for example, when transported in a car), allow it to come to room temperature before testing your blood.
- ◆ Store your OneTouch® SureStep® Test Strips in a cool, dry place in the original moisture-resistant, light-protected vial. The storage temperature should not exceed 30°C (86°F). Do not refrigerate or

freeze test strips. Do not store test strips in any container other than the original vial. Do not expose the test strips to bleach fumes or hydrogen peroxide.

- ◆ Store OneTouch® SureStep® Glucose Control Solution in a cool, dry place at temperatures below 30°C (86°F). Do not refrigerate or freeze. When used for testing, the solution should be at normal room temperature above 15°C (59°F).
- ◆ Store the lancing device and lancets at room temperature out of reach of small children.

Troubleshooting

Error Messages

The meter displays specific error messages when it detects a problem with its internal functions, the test strip, or the testing procedure. However, the meter cannot detect all problems. For those problems the meter can detect, it may not produce an error message every time the problem occurs.

When an error message appears, note the error number, turn off the meter, and then resolve the error by following the instructions in the **“What to Do”** column. If you continue to get an error message after following the instructions, call the OneTouch® Customer Care Line at 1 800 663-5521.

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
	<p>The meter detected a problem with the colour of the test strip.</p>	<ol style="list-style-type: none"> Repeat the test with a new test strip and: <ul style="list-style-type: none"> ◆ Apply blood to the pink test square only. ◆ Make sure the confirmation dot turns completely blue. ◆ Insert the test strip within 2 minutes of applying blood.
		<p><i>Note: If you continue to get any error message, compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</i></p>
		<ul style="list-style-type: none"> ◆ Clean the meter and repeat the test with a new test strip. See Caring for Your Meter for further information. <ol style="list-style-type: none"> If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
<p>Er 1 <i>continued</i></p>	<p>The meter detected a problem with the colour of the test strip.</p>	<p>Professional Use:</p> <ul style="list-style-type: none"> ◆ Test the sample with an alternate method. While awaiting results: Compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial label. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), take appropriate precautions immediately. Before retesting verify that the sample hematocrit is within range. ◆ If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

The meter detected a problem with the test strip.

1. Repeat the test with a new test strip and:
 - ◆ Check the test strip expiration date. Do not use test strips beyond the expiration or discard date.


NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
<p>Er 2 <i>continued</i></p>	<p>The meter detected a problem with the test strip.</p>	<ul style="list-style-type: none"> ◆ Apply the sample to the pink test square only. ◆ Make sure the confirmation dot turns completely blue. ◆ Wait approximately 1 minute, but not more than 2 minutes, after applying sample before inserting the test strip into the meter. ◆ Do not alter the test strips. <p><i>Note: If you continue to get any error message, compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</i></p> <ul style="list-style-type: none"> ◆ Clean the meter and repeat the test with a new test strip. See Caring for Your Meter for more information. <p>2. If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663–5521.</p>

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
	<p>The meter detected a problem with the position of the test strip.</p>	<ol style="list-style-type: none"> Repeat the test with a new test strip and: <ul style="list-style-type: none"> Clean the meter, especially the contact points, before retesting. See Caring for your Meter for more information. <p><i>Note: Before cleaning the meter, compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</i></p> Do not move the test strip after you insert it into the meter. Remove the test strip only when the test is complete. If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
	The meter detected a temperature problem.	<ol style="list-style-type: none"><li data-bbox="608 186 1416 326">1. Repeat the test with a new test strip. Make sure the meter and test strip are at room temperature.<ul style="list-style-type: none"><li data-bbox="654 290 1262 326">◆ Do not apply heat or cold to the meter.
		<p data-bbox="685 352 1447 564">Note: If you continue to get any error message, compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</p> <ol style="list-style-type: none"><li data-bbox="608 590 1416 663">2. If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
	<p>The meter detected a problem with the optical system.</p>	<ol style="list-style-type: none"> 1. Repeat the test with a new test strip and: <ul style="list-style-type: none"> ◆ Clean the meter, especially the grey area inside the test strip holder, before retesting. See Caring for Your Meter.
		<p>Note: Before cleaning the meter compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</p> <ul style="list-style-type: none"> ◆ Make sure the test strip holder is fully inserted into the meter. See Caring for Your Meter for more information. ◆ Wait for the “Insert test strip” symbol to appear on the display before inserting the test strip. <ol style="list-style-type: none"> 2. If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.


NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

ERROR	POTENTIAL PROBLEM	WHAT TO DO
	<p>The meter detected a problem with the electronics.</p>	<ol style="list-style-type: none"> 1. Repeat the test with a new test strip and: <ul style="list-style-type: none"> ◆ If you use bleach to <i>disinfect</i> the meter, before retesting, wipe the meter and contact points with a cloth moistened with water to remove residual bleach and thoroughly dry the meter. Failure to remove the excess bleach may result in an error message (<i>Er5</i>, <i>Er6</i>) or an inaccurately high result. ◆ Recheck the code number entered in the meter and make sure it matches the code number on the test strip package. <p>Note: If you continue to get any error message, compare the colour of the test strip confirmation dot to the Colour Chart on the test strip vial. If the glucose level is very low (2.8 mmol/L) or very high (19.4 mmol/L), follow the advice of your Healthcare Professional.</p> 2. If these actions do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

Symbols

The following symbols appear on the meter display to help guide you through the testing process. If they occur during testing, they may indicate problems.

SYMBOL	WHAT IT MEANS	WHAT TO DO
--------	---------------	------------

Insert Test Strip

Firmly push the test strip into the test strip holder, with the tapered-end first and the confirmation dot facing down, until it stops. If you fail to completely insert the test strip, the test may start. However, you may receive an inaccurately low result.


Remove Test Strip

Remove the test strip from the meter. If white patches or streaks appear on the confirmation dot, you have not applied enough blood to obtain an accurate test. You must start over with a new test strip. When the “Insert Test Strip” symbol appears on the display, insert the test strip. Push it firmly and more slowly into the test strip holder until it stops.

Note: *It is possible to push the test strip into the test strip holder too fast. If this happens, the meter displays the “Remove Test Strip” symbol.*

NOTE: If the actions listed below do not resolve the error, call the OneTouch® Customer Care Line at 1 800 663-5521.

SYMBOL **WHAT IT MEANS** **WHAT TO DO**

Replace battery

Replace the batteries (see **Caring for your Meter**) as soon as possible. If this symbol appears constant, batteries are low. If this symbol is flashing, batteries are dead and the meter will no longer operate.

OTHER SYMBOLS

SYMBOL **WHAT IT MEANS**


Apply blood to test strip.

See p. 26-29

Meter and/or test strip holder need cleaning.

See p. 42-44

Meter is reading your result.

See p. 31

SYMBOL **WHAT IT MEANS**

Result from meter memory.

See p. 65-67

Result is over 27.8 mmol/L (500 mg/dL).*

See p. 36 **IMPORTANT**

Match meter code with the code on the vial. *See p. 16*

* *If you have never experienced blood glucose levels this high or have not been instructed in how to deal with them by your Healthcare Professional, you should contact your Healthcare Professional at once. If this reading does not reflect how you feel, repeat the test with a new test strip.*


Changing Meter Settings (Set-Up Mode)

The OneTouch® SureStep® Meter has been preset to meet the needs of most users in Canada. There are six features you can change in addition to the date and time.

FEATURE	PRESET	ALTERNATIVE	SEE PAGE
Set time and date	AM/PM month/day/year	Any	59
Automatic display of last result	On	Off	61
Unit of measurement	mmol/L	mg/dL	61
Beep signals during testing	On	Off	62
14- and 30-day averages	On	Off	63
Time format	AM/PM	24 hour	64
Date format	Month/Day	Day/Month	64

Entering Set-Up Mode

1. Turn off the meter.
2. Press and hold down the green **C** button.
3. Press and release the blue power button while holding down the green **C** button.
4. Release the green **C** button.

5. Each setting flashes on the display until you press the blue power button to move to the next feature.

6. You can change the setting while it flashes by pressing the green **C** button.

NOTE: *To move quickly through the features, press the blue button repeatedly or press and hold it down to scroll.*

Setting the Date and Time

The OneTouch® SureStep® Meter records the time and date of each glucose test you perform. You may have to change the time and date, for example, when the batteries are changed. You will also need to reset the time when daylight savings time starts or ends. The accuracy of your glucose results is not affected if the date and time settings are incorrect.

When setting the date and time, if you pass the desired number (for example, day, month year), continue pressing the green **C** button to cycle through the values and back to the desired number.

Setting the Time

To set the time, enter the set-up mode as described in “Entering Set-Up Mode” above. The time appears with the hour flashing.

Press and release the green **C** button to advance 1 hour. To move faster, hold down the green **C** button. With the correct hour on the display, press and release the blue power button and the minutes will start flashing.

Press and release the green **C** button to advance 1 minute. To move faster, hold the green **C** button down. With the correct minutes on the display, press the blue power button and the year setting will appear on the display.


Note: *Calculated 14-day and 30-day averages may change if you change the date and time. The averages are calculated from glucose results measured during the 14 or 30 days preceding the current date and time settings.*

Setting the Year

With the year on the display, press and release the green **C** button to advance 1 year. To move faster, hold down the green **C** button. With the correct year on the display, press and release the blue power button and the month and day will appear on the screen with the month segment flashing (month/day format is preset).

Setting the Month

Press and release the green **C** button until the correct month number appears. To move faster, hold down the green **C** button. With the correct month on the display, press the blue power button and the day segment will start flashing.

Setting the Day

Press and release the green button until the correct day appears. To move faster, hold the green button down. With the correct day on the display, press the blue power button and the display of last result setting will appear.

Display of Last Result

The meter is preset to display your last glucose result when you turn on the meter.

To turn this feature off, press the green **C** button. This symbol will flash, indicating that automatic display is turned off. Press the blue power button to move to the unit of measurement setting.

Unit of Measurement

The meter is preset to display test results in millimoles per liter (mmol/L), the unit of measure most often used in Canada. You may change the setting to milligrams per deciliter (mg/dL) the unit of measure used in the

United States. If you are not sure which unit of measurement to use, please check with your Health Care Professional. The mmol/L results include a decimal point and the mg/dL results do not.

To change the unit of measurement to mg/dL, press the green **C** button. This symbol will flash on the display.

Press the blue power button to move to the beep signals setting.

Beep Signals

The OneTouch® SureStep® Meter uses beeping tones to signal certain steps in the procedure. However, these “beep” signals do **not** mean that the step was performed correctly. Three types of beep signals sound.

1. A single, short beep identifies each step in the test procedure.
2. A single, long beep accompanies error messages, a “HI” result, and certain mistakes in doing the test procedure.
3. A three-tone beep sounds when the meter displays a test result and when the meter moves to the beep signal feature in set-up mode.

In set-up mode, a three-tone beep sounds and the following symbol flashes when you move to the beep signals feature.

To operate without beep signals, press the green **C** button. The following display will appear.

You will hear a three-tone beep when *OFF* is displayed. Press the blue power button to move to the 14- and 30-day test averages setting.

14- and 30-Day Test Averages

The meter is preset to display 14-day and 30-day averages of your blood glucose results in memory mode.


For accurate 14-day and 30-day blood glucose averages, you must mark control solution test results by pressing the green **C** button again while the glucose level is displayed. This will keep them from being included in your blood glucose test averages. If you wish to turn the averages off, press the green button. The following display will appear:

No averages will be displayed in memory mode. Press the blue power button to move to the time format setting.

Time Format

The meter is preset to display time in the AM/PM format. An alternative 24-hour time format is available.

To change to the 24-hour time format, press the green **C** button. Press the blue power button to move to the date format setting.

Date Format

The meter is preset to display the date in month/day format. An alternate day/month format is available.

To change the display of the date to the day/month format, press the green **C** button. Press the blue power button to turn off the meter.

Exiting the Set-Up Mode

You may exit the set-up mode at any time by pressing and holding the blue power button down until the display is blank.

Getting Results from Meter Memory

Background Information on Meter Memory

When you turn on the OneTouch® SureStep® Blood Glucose Meter you will see the result of your last test. The meter stores your last 150 results (blood glucose or control solution test) with the date and time in its memory. When the meter's memory is filled with 150 results, the oldest result is dropped and the newest one is added. The meter also gives 14-day and 30-day averages of your blood glucose test results. Results marked as control solution results will not be included in your averages.

Getting All Results from Meter Memory

To get the results stored in the meter's memory:






1. Press the blue power button and hold it down for about 5 seconds.
2. When the display flashes, let go of the blue button.

You are now in memory mode.

3. First, the 14-day average will appear. This is the average of your blood glucose results during the last 14 days. (Ask your Healthcare Professional how to use these averages.)
4. Then, when you press the green **C** button, your 30-day average appears. (Ask your Healthcare Professional how to use these averages.)


Note: A “HI” result is included in your 14-day and 30-day averages as 27.8 mmol/L.

- Each time you press and let go of the green **C** button, the meter displays the next most recent of your 150 stored test results. The date and time of each result is displayed alternately under the test result. If the test result was marked as a control solution result, **CEL** appears under the result.
- To scroll quickly through the results, press and hold the green **C** button.

(Examples)

- If a test resulted in an error message, the error message is stored in the memory.


(Example)

Note: A "HI" result is included in your 14-day and 30-day averages as 27.8 mmol/L.

First Time Turning on Meter

The first time you turn on your meter you will see the following symbol. This shows that there are no test results in memory.

You will also see this symbol after you turn on your meter when your last test ended in an error and the meter did not display a glucose result.

Uploading Memory Data

Test results stored in the memory may be uploaded to a personal computer using OneTouch™ Diabetes Management software.

See www.OneTouch.ca or call The OneTouch® Customer Care Line at 1 800 663-5521 for more information.

Specifications

The performance of the OneTouch® SureStep® Blood Glucose Monitoring System was tested under the conditions below. See **“Important Technical and Healthcare Professional Information”** in the OneTouch® SureStep® Test Strips package insert for typical accuracy and precision data.

Power Supply: Two AAA alkaline batteries

Battery Life: Approximately 18 months when doing two tests per day

Result Range: 0–27.8 mmol/L (0–500 mg/dL). Higher results are virtually always displayed as **“HI”**

Display: Liquid crystal

Blood Source: Whole blood

Calibration: Plasma equivalent

Hematocrit Range: 25–60%

Size: 8.9 cm x 6.1 cm x 2.0 cm (3.5" x 2.4" x 0.8")

Weight: 108 grams (with batteries) (3.8 ounces)

Operating Temperature Range*: 10–35°C (50–95°F)

Operating Humidity Range*: 10–90% relative humidity (noncondensing)

Memory: 150 results (blood glucose and control solution).

*See **“Comparing Meter and Laboratory Results”** for additional information.

Guarantee and Warranty

30-Day Money-Back Guarantee

If you are not fully satisfied with the OneTouch® SureStep® System, you may receive a full refund by calling a OneTouch® Customer Services Representative at 1 800 663-5521 (Canada) within 30 days of purchase. You will need to return the OneTouch® SureStep® Meter and a copy of your receipt to receive a refund. The Warranty Policy does not apply to the performance of the OneTouch® SureStep® Meter when used with any test strip other than OneTouch® SureStep® Test Strips from LifeScan, or when the OneTouch® SureStep® Meter or Test Strips are modified in any way.

This warranty is in lieu of all other warranties, express or implied, including any implied warranty of merchantability or fitness for a particular purpose.

Three-Year Warranty

The OneTouch® SureStep® Meter has a full 3-year warranty from the date of purchase. If, at any time during the first 3 years after purchase, the meter does not work for any reasons (except for obvious abuse), LifeScan will replace it with a new meter or equivalent product free of charge.

The Warranty Policy applies only to the original purchaser of this meter and does not include the batteries supplied with the meter.


Please complete the Warranty Registration Card and mail it to LifeScan.

Write your date of purchase here:


Before you return your meter or any product, first call your OneTouch® Customer Services Representative at 1 800 663-5521 (Canada) for instructions.

CODING THE METER


Before testing, make sure the meter is coded correctly.


1 Press blue power button to turn on meter. Make sure a test strip is not in the meter when you turn it on.

2 The meter automatically displays the last test result with the date. . .


3 . . . and then the code number is displayed afterward.


4 Compare the code number on the meter to the code number on the test strip vial. If the code numbers do not match, press the green **C** button until the code number on the meter matches the code number on the test strip vial.


PERFORMING A TEST


1 Press blue power button to turn on meter. Make sure a test strip is not in the meter when you turn it on.

2 Apply blood to the centre of the **pink test square**. Never add a second drop of blood.


3 Check confirmation dot on back of test strip and make sure it is completely blue.

4 Firmly push the test strip into the test strip holder, with the tapered-end first and the confirmation dot facing down, until it stops. Read result. Do not remove strip until countdown is complete.

BLOOD APPLICATION


How you can be sure that you have applied the right amount of blood:


Pink test square is where you apply a drop of blood.

White pad absorbs the excess blood that may extend beyond the pink test square. Do not apply blood directly onto the white pad.

Confirmation dot on back of test strip is where you check to confirm proper blood application.


Apply a drop of blood to the centre of the pink test square on the front of the test strip. The confirmation dot on the back should be completely blue for an accurate test.

If any white patches or streaks appear on the blue confirmation dot, you have not applied enough blood for an accurate test. You must start the test over with a new test strip to avoid inaccurately low results.

** Never add a second drop of blood.*

If there is no white visible on the white pad, you have applied too much blood. You must start the test over with a new test strip to avoid inaccurately high results.

** Be careful not to apply too much blood, especially if you tend to bleed easily or are taking blood-thinning medications.*

Numerics

14-and 30-day averages 19, 60, 63, 65

A

Accuracy 39

Applying sample

- Adequate amount 30
- Blood 28-30
- Control solution 17
- Not enough 30
- Pink test square 11, 17, 28-29
- Too much 17, 23, 29

B

Batteries

- Changing 44-46
- Symbol indicating replacement 12, 44

Battery compartment 7, 45-46

Beep signals 62-63

Blood

- Applying 28-30
- Obtaining 24-25, 28

Blood glucose

- Results 31-33
- Target ranges 35, 40
- Test 22-34

C

C Button 7

- Changing code number 14

- Changing meter settings 58-64

- Marking control solution result 19

Changing

- Batteries 44-46
- Code number 14
- Meter settings 58-64

Checking

- Code number 16, 26
- Confirmation dot 18, 30
- System 15-21

Cleaning

- Contact points 43
- Disinfecting 41
- Display symbol 12, 42
- Lens area 43
- Meter 41-44
- Test strip holder 42-43

Code number 11

- Changing 13-14
- Displayed at meter power up 13, 16, 26

Coding the meter 13-14

Colour Chart 9

- Checking confirmation dot against 32-33
- Unused test strip 7, 23, 27
- Very low or very high 9, 36

Comparing meter and lab results 3, 38-40

Confirmation dot 8

- Checking (control solution test) 18

- Checking (blood test) 30
 - Colour of unused test strip 17, 23, 27
 - Comparing to Colour Chart 9, 32-33
 - Facing down when inserting test strip . . . 11, 18, 31
 - Patches or streaks appear 30
- Control solution
- Applying 17
 - Discard date 4, 17
 - Expiration date 4, 17
 - Marking result 19
 - Out of range 17, 20-21
 - Results 20-21
 - Storage 47
 - Test 15-21
- Contact points 7
- Cleaning 43
- D**
- Data port 7
- Date and time
- Changing 59-61
 - Format 64
- Discard date
- Control solution 4, 17
 - Test strip 4, 16, 23, 27
- Display
- Of last result 61
 - Screen 7
- Symbols 10-12
- E**
- Error message
- Instead of blood glucose result 32-33, 35
 - Troubleshooting 48-57
- Expiration date
- Control solution 4, 17
 - Test strip 4, 16, 23, 27
- G**
- Glucose
- See also control solution
 - Levels 2
 - Results (blood) 31-33
 - Target ranges 35, 40
 - Test 22-34
- H**
- Hematocrit 39, 50, 68
- “HI” result 12, 32, 37, 66
- Humidity 39, 68
- L**
- Lancing device 24-25, 33-34
- Storage 47
- Lens area 7
- Cleaning 43
- Logbook
- Recording blood glucose results 33

- Recording control solution results 20

M

Matching code number 11, 16, 26

Memory 65-67

Meter

- Care 41
- Cleaning 41-44

Coding 13-14

- Description 7
- Display symbols 10
- Serial number 7
- Settings 58-64
- Storage 47
- Turning on 16, 25

Meter settings

- 14- and 30-day averages 63
- Beep signals 62-63
- Date format 64
- Date and time 59-61
- Display of last result 61
- Time format 64
- Unit of measurement 61-62

O

Optional Colour Chart 9

P

Pink test square 8, 11

- Applying blood 28-29

- Applying control solution 17

– Plasma glucose value 35, 38

– Power button 7, 16, 25

– Precautions 4-5

R

Result

– 14- and 30-day averages 19, 63

– Blood glucose 31-33

– Comparing to lab result 3, 38-40

– Control solution 20-21

– With date and time 12

– Getting results from meter memory 65-67

– High 12, 33, 36

– Low 36

– Out-of-range (control) 20-21

– Recording in logbook 20, 33

– Target ranges 35, 40

– Tips to ensure accuracy 23

– Turning on the meter 10

– Understanding 35-40

– Units 12

– Uploading to PC 67

– Using the Colour Chart 9, 32-33

S

Segment 10

Settings (meter) 58-64

Symbols

- Apply sample 11, 26

– Battery	12, 44, 57
– Clean meter	12, 42
– Clock	12, 19, 31
– Descriptions	56-57
– Empty memory	65
– Insert test strip	11, 26, 56
– Remove test strip	56

T

Temperature	39, 68
Test strip	8
– Code number	11, 14
– Colour of unused confirmation dot	17, 23, 27
– Discard date	4, 16, 23, 27
– Expiration date	4, 16, 23, 27
– Inserting firmly	18, 31
– Storage	46-47
– Storage precautions	27
Test strip holder	7
– Cleaning	42
– Removing	42
– Replacing	43
Troubleshooting	48-55

U

Unit of measurement	iv, 12, 61-62
---------------------------	---------------

W

Warnings	4-5
White pad	8

FOOTNOTES

- 1 American Diabetes Association, “Implications of the Diabetes Control and Complications Trial,” (Position Statement), *Diabetes Care* 1988; 21; (Supplement 1); S88-S90.
- 2 Wickham, NWR et al.: *Practical Diabetes*, Vol. 3, No. 2 (1986), 100.
- 3 Cohen, FE et al.: *Diabetes Care*, Vol 9, No. 3 (1986), 320–322.
- 4 Atkin, SH, et al. Fingerstick glucose determination in shock. *Ann of Intern Med.* 1991;114:1020–1024.
- 5 Surwit RS, Feinglos MN *Diabetes Forecast* 1988; April 49-57
- 6 Sacks, DB in “Carbohydrates”, Burtis, CA, Ashwood, ER (ed), *Tietz Textbook of Clinical Chemistry*. Philadelphia, WB Saunders Company (1994), p. 959
- 7 Sylvain, HF, et al. *American Journal of Critical Care.* 1995; 4:44–48.
- 8 Rushakoff RJ, Lewis SB. Case of Pseudohypoglycemia (letter), *Diabetes Care* 2001; 24:2157-2158.

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Owner's Booklet
Guide d'utilisation

Download from www.somanuals.com. All Manuals Search And Download.

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