**Heart Rate
Monitor**
Instruction Manual

- English
- Français
- Español
- Português

OMRON
A Good Sense of Health
HR-100C

Heart Rate Monitor

HR-100C

Omron Healthcare, Inc.
a world leader in the
development of wellness
products, offers a variety
of innovative blood pressure
monitors, respiratory
products, body fat analyzers,
massagers, heart rate
monitors and pedometers
to help improve personal
health and lifestyle.

Distributed by:
Omron Healthcare, Inc.
300 Lakeview Parkway
Vernon Hills, Illinois 60061

Toll Free Customer Service
1-800-634-4350

www.omronhealthcare.com

Made in China

HR-100CINST

Table of Contents

Warning	1
Introduction	2
Product Features	3
Components	3
Buttons and Functions	3
Attaching Your Heart Rate Monitor	5
General Precautions	6
Specifications	6
One Year Limited Warranty	7

Warning

Overexertion can cause serious injury, including heart attacks. Some individuals cannot safely elevate their heart rate to the levels of typically used heart rate training zones. No one should begin an exercise program without first obtaining medical clearance, especially if there is a personal or family history of heart disease, high blood pressure, if you are over age 40, have diabetes, high cholesterol, smoke cigarettes, are overweight or are taking certain medications. Stop exercising and seek medical attention if you notice signs of overexertion or heart problems, such as pain or pressure in the left or mid-chest area or left part of the neck, shoulder or arm, light-headedness, cold sweat, unusual paleness or fainting. Also note that the signals used by this monitor may interfere with a pacemaker or other implanted devices. Consult the manufacturer of the implant device and/or your physician prior to using this monitor.

Introduction

Congratulations!

Welcome to the future of wireless heart rate technology! Your purchase shows that you are serious about your exercise and are committed to achieving your best. Today, getting the most from your exercise doesn't necessarily mean you have to spend longer periods of time to achieve your goals. With Target Heart Rate Training, you will increase the effectiveness of your workout by monitoring and quantifying your results every step of the way. For most individuals as little as 20 minutes in your target heart rate zone each session is enough to achieve substantial health improvements. The Omron HR-100INT will act as your personal coach and guide you to your fitness goals more effectively.

Mission Statement

At Omron Healthcare we have listened to what people such as yourself have been asking for in a heart rate monitor and are committed to your satisfaction. We know you will enjoy the added features, along with those never found before in other heart rate monitors, and the ease of use incorporated into our products. If you ever have any questions, comments or suggestions, please write or call our staff. We greatly value your feedback and look forward to continue to provide you with products and support that is second to none. Thank you for your purchase and anticipated support.

Product Development Team


Product Features

- Shock resistant design
- Water resistance to 98 feet (30 meters)
- Wireless ECG accurate heart rate chest belt

Components

Your package of the HR-100INT includes:

- A. Pulse watch receiver (Heart Rate Monitor)
- B. Chest Transmitter (Transmitter Belt)

Buttons and Functions

A. Heart Rate Reading

Starting from the Time Mode:
Press MODE button to enter the
Heart Rate Mode.

Note: If no button is pressed or no heart rate reading is done for 5 minutes, the watch will automatically return to Time Mode.

B. Time and Alarm Setting

Start from the Time Mode:

1. Press & hold SET/SIGNAL button to enter the time setting.


2. Hour (blinking)
- Press MODE button to adjust.
- Press SET/SIGNAL button to proceed to the next item.

3. Minutes (blinking)
- Press MODE button to adjust.
- Press SET/SIGNAL button to proceed to the next item.

4. 12/24 HR Time Display (blinking)
- Press MODE button to adjust.
- Press SET/SIGNAL button to proceed to the next item.

5. Daily Alarm (blinking)
- Press MODE button to toggle ON/OFF.
- Press SET/SIGNAL button to proceed to the next item.

6. Alarm Hour (blinking)
- Press MODE button to adjust.
- Press SET/SIGNAL button to proceed to the next item.


7. Alarm Minutes (blinking)
- Press MODE button to adjust.
- Press SET/SIGNAL button to exit the setting mode when finished.

Note: By pressing the SET/SIGNAL button once in Time Mode can also turn the alarm ON or OFF.


Attaching Your Heart Rate Monitor

Step I: Attach the strap to the transmitter


Step II: Adjust the tension of the strap to fit comfortably around your chest. Be sure to wrap and center the belt, just below your muscles or breasts: (Fig.2)

Fig.2

Step III: Moisten the transmitter electrodes and wear the belt directly on your skin to get the best results (Fig.3)

Fig.3

Step IV: Wear your Heart Rate Monitor, within 3 feet of the transmitter.

Step V: If you prefer to wear the transmitter over a light shirt, be sure to moisten the shirt beneath the electrodes to insure proper signal pickup.

General Precautions

Should your heart rate monitor need battery replacement we recommend that you take it to a local jewelry store in order to guarantee that it will maintain its water resistance. Improper battery replacement or service will void the warranty. Repair should only be performed by Omron Healthcare, Inc., 300 Lakeview Parkway, Vernon Hills, IL 60061. (Customer Service Toll Free 1-800-634-4350).

Although your heart rate monitor is water resistant, be sure not to operate the device under water as pressure can cause it to leak.

The transmitter will automatically shut itself off once you have disconnected it from your chest. Dry it off after every use.

Be sure not to excessively twist or bend the electrode portions of the transmitter belt as it may cause mechanical damage to the internal circuit.

Specifications

Water Resistance 98 feet (30 meters)

Operating Temperature 32°F-122°F
(0°C-50°C)

Battery Type Watch and
Transmitter Belt CR2032

Battery Life Approximately 1.5 years,
average use of 1 hr./day, 7 days/wk.

One Year Limited Warranty

This product is warranted by Omron Healthcare, Inc. for a period of one year from date of original purchase (with the exception of battery and elastic strap) under normal use, to be free from defects in material and workmanship. This warranty does not cover damage caused by accidents, misuse, abuse, tampering or disassembling of the product. Omron Healthcare, Inc. reserves the right to determine the cause of malfunction and will repair or replace your Omron HRM at no charge except for shipping and handling. For warranty service, you must provide proof of purchase and date of purchase, along with a check or money order for \$5.00 to cover the costs of return shipping and handling. Warranty is voided if repairs have been made by individuals not authorized by Omron Healthcare, Inc.

No responsibility is assumed for any incidental or consequential damages, including, but without limitation, damages or injuries resulting from possible measurement inaccuracies or malfunctions.

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