# Reusable SpO<sub>2</sub> Sensors

*Instruction for Use*

## *A-Series Sensors*

**M1191A**  
**M1191AL**  
**M1192A**  
**M1193A**  
**M1194A**  
**M1195A**

## *Special Edition A-Series Sensors*

**M1191ANL**  
**M1192AN**  
**M1193AN**  
**M1194AN**  
**M1195AN**

## *B-Series Sensors*

**M1191B**  
**M1191BL**

# PHILIPS

# About This Edition

## Notice

### Copyright © 2006

Koninklijke Philips Electronics N.V.  
All rights are reserved.

## Printing History

### Current Edition/Revision:

Printed in Germany December 2006

### Previous Editions/Revisions:

Printed in Germany June 2006

Edition 1, Revision E.00.00

Printed in Germany August 2004

Edition 1, Revision D.00.00

Printed in Germany December 2003

## Contents

Instructions for Use.....	1	Návod na použitie .....	127
Manuel d'utilisation .....	10	Navodila za uporabo .....	136
Gebrauchsanweisung .....	19	Instruçiuni de utilizare .....	145
Gebruiksaanwijzing .....	28	Инструкции за употреба .....	154
Instrucciones de Uso .....	37	Kasutusjuhend .....	163
Instruções de Uso .....	46	Naudojimo instrukcija.....	172
Istruzioni d'uso .....	55	Lietošanas norādījumi .....	181
Bruksanvisning .....	64	Инструкция по эксплуатации	190
Bruksanvisning .....	73	Οδηγίες Χρήσης .....	199
Käyttöopas .....	82	Kullanım Talimatları.....	208
Brugerhåndbog.....	91	ユーザーズ・ガイド .....	217
Használati útmutató .....	100	使用手冊 .....	226
Instrukcja obsługi.....	109	使用说明书.....	236
Návod na použití .....	118		

United States law restricts this device to sale by or on the order of a physician.

OxiMax® is a registered trademark of Tyco Healthcare Group LP, Nellcor Puritan Bennett Division.

Download from [Www.Somanuals.com](http://www.Somanuals.com). All Manuals Search And Download.

## Product Overview

Three series of Philips reusable SpO<sub>2</sub> sensors offer monitoring capability for a broad range of patient types as shown in the table below. Reference *Specifications* later in this IFU for sensor accuracy.

- A-Series Sensors:** For use with any Philips/Agilent/HP instrument that lists them as an accessory. Refer to your instrument's Instructions for Use for A-Series sensor compatibility information.
- B-Series Sensors:** For use with the same Philips/Agilent/HP instruments that list A-Series M1191A and M1191AL sensors as compatible accessories; i.e. M1191B can replace M1191A and M1191BL can replace M1191AL.
- Special Edition A-Series Sensors (AN-Series):** For use with Philips brand OxiMax® compatible instruments M1020B or M3001A, Option A02. Also for use with the same Philips/Agilent/HP instruments that list A-Series sensors as compatible accessories. Each Special Edition A-Series sensor can replace its similarly numbered A-series sensor; i.e. M1191ANL can replace M1191AL; M1192AN can replace M1192A; M1193AN can replace M1193A; M1194AN can replace M1194A; and M1195AN can replace M1195A.

A-Series Sensors	B-Series Sensors	<sup>2</sup> Special Edition A-Series Sensors	Adult		Pediatric		Infant		Neonate	
			Finger	Ear	Finger	Ear	<sup>3</sup> Finger	Toe	<sup>3</sup> Hand	Foot
M1191A	M1191B	N/A	✓							
<sup>1</sup> M1191AL	<sup>1</sup> M1191BL	<sup>1</sup> M1191ANL	✓							
M1192A	N/A	M1192AN			✓					
M1193A	N/A	M1193AN							✓	✓
M1194A	N/A	M1194AN		✓		✓				
M1195A	N/A	M1195AN					✓	✓		

<sup>1</sup>These sensors have a Longer 3m (versus standard 2m) attached cable.

<sup>2</sup>These sensors extend compatibility to Philips brand OxiMax® compatible instruments.

<sup>3</sup>Preferred application sites for M1193A and M1195A sensors.

## Intended Use

Philips reusable SpO<sub>2</sub> sensors are for multi-patient use, when continuous non-invasive arterial oxygen saturation and pulse rate monitoring are required.


## Accessories Supplied with Sensors

- **M1627A Adult Wristband:** *Supplied only with adult finger sensors to secure sensor cable to back of adult patient's hand and wrist.*
- **Instructions for Use:** Describes sensor models, patient applications, Warnings against improper use, and specifications.

## Accessories Available for Separate Purchase

- **M1940A Adapter Cable:** 8-pin (round female) to 12-pin (round male) *adapter* cable. Adapts 8-pin sensor connector to 12-pin instrument jack. Extends sensor cable length by 2m. DO NOT USE with M1191AL, M1191ANL, M1191BL sensors.
- **M1941A Extension Cable:** 8-pin (round female) to 8-pin (round male) *extension* cable. Extends sensor cable length by 2m. DO NOT USE with M1191AL, M1191ANL, M1191BL sensors.

## Definition of Product Symbols

 Adult	 Pediatric	 Adult Finger	 Infant/Neo Hand	 Adult/Pediatric Ear
 Infant	 Neonate (Neo)	 Pediatric/Infant Finger	 Infant Toe	 Infant/Neo Foot
 Caution, See Documentation	<b>R<sub>x</sub> only</b> Prescription Use Only	 Storage Temp	 Contains no Latex	

---

## Warnings

- These sensors are for use only with Philips/Agilent/HP instruments. Before using, verify the combination of instrument/sensor is specified in the instrument's user documentation (e.g. Instructions for Use); otherwise patient injury may result.
- Pulse oximetry measurements are statistically distributed. Two-thirds of all pulse oximetry measurements can be expected to fall within the stated accuracy (refer to *Specifications* later in this IFU for stated sensor accuracy).
- Connect the sensor only to the SpO<sub>2</sub> connection or the SpO<sub>2</sub> adapter cable of the oximeter.
- Do not reuse the sensor on a different patient until it has been disinfected. The sensor can be reused on the same patient for the patient's entire stay.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41° C on the skin if the initial skin temperature does not exceed 35° C.
- Be sure to apply the sensor to the patient at the preferred or alternative application site, while following application instructions provided later in this document. Failure to do so can cause inaccurate measurements.
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, make sure you are using a sensor of the correct size and that the sensor is not too tight. If the sensor is too tight, because the application site is too large or becomes too large due to edema, the excessive pressure applied may result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.
- If the sensor is applied too loose it may fall off or compromise proper alignment of sensor optics and result in inaccurate readings.

## Warnings (cont.)

- Where possible, the application site for the sensor should be an extremity free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Avoid sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate measurements.
- Make sure the application site for the sensor is not deeply pigmented, or deeply colored. For example, nail polish, artificial nails, dye or pigmented cream may cause inaccurate measurements. In any of these cases reposition the sensor or choose an alternative sensor for use at a different site.
- Cover the sensor with opaque material under conditions of strong or excessive light (infrared lamps, OR lamps, photo therapy). Failure to do so can result in inaccurate measurements.
- Protect the connector from contact with any liquid.
- Inspect the sensor application site every 2 to 3 hours to ensure skin integrity, correct optical alignment, and circulation distal to the sensor site. Skin irritation or ulceration may occur if the sensor is attached to one location for too long. Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised. If the light source is not directly opposite the light detector, reapply the sensor, or choose an alternative sensor for use at a different site.
- Do not use a sensor during MRI scanning. This may cause burns or inaccurate measurements.
- **M1191AL/M1191ANL/M1191BL only:** Do not use these sensors with an extension or adapter cable.
- **M1195A/M1195AN only:** Infant patients are prone to movement; therefore, when applying a sensor to an infant, secure the sensor cable with tape to prevent it from coming off the infant's finger. If necessary, position the sensor on a different finger, or an alternative site (toe).

---

## Before Applying the Sensor

Be sure to read and understand all Warnings listed in your SpO<sub>2</sub> monitoring instrument IFU, as well as all Warnings described in this sensor IFU. Use sensors only with validated Instruments and at a recommended patient application area.

### Inspect Sensor for Damage

- Inspect the sensor outside and inside. To check the inside, gently open the sensor cavity and check for splits on or next to the transparent silicone that covers the optical elements. Be sure there are no blisters on the silicone, and no liquid leaking from sensor optics.
- Any sensor showing signs of damage or alteration must not be used for further patient monitoring; instead, dispose of it using proper disposal procedures (see below).

---

## Sensor Disposal

### Dispose Deteriorated Sensors

Any sensor showing signs of physical or electrical deterioration or failure should be disinfected, decontaminated, and disposed of according to local laws regarding the disposal of hospital waste.

---

## Change Application Site Periodically

Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised.

# Applying Finger (or Toe) Sensors

## Adult Finger Sensors

Any finger except thumb on patients more than 50kg.

- **A-Series:** M1191A or M1191AL
- **Special Edition (SE) A-Series:** M1191ANL
- **B-Series:** M1191B or M1191BL

*adult application  
shown*

## Pediatric Finger Sensors

Any finger except thumb for patients between 15kg and 50kg

- **A-Series:** M1192A
- **Special Edition (SE) A-Series:** M1192AN

## Infant Finger (or Toe) Sensors

Any finger or toe (except thumb) for patients between 4kg and 15kg.  
Finger or toe diameter should measure from 7-8mm (0.27-0.31”).

- **A-Series:** M1195A, or
- **Special Edition (SE) A-Series:** M1195AN

Step	Applying the Finger (or Toe) Sensor
1	Select appropriate sensor for patient size (defined above).
2	Position sensor over patient's finger. If applying an adult finger sensor, be sure the sensor cable is positioned <b>ON TOP</b> of the finger/hand as shown above.
3	Patient fingertip should touch but not protrude from end of sensor. If necessary, trim the fingernail to correctly position the sensor.
4	If applying an adult finger sensor, secure the cable to <b>BACK OF HAND</b> with M1627A wristband (supplied only with adult finger sensors.)
5	Plug sensor into instrument (or into adapter cable if needed).
6	Inspect and change the sensor application site periodically.

---

## Applying Hand/Foot Sensors

Neonatal hand or foot, for patients between 1kg and 4kg.

- **A-Series:** M1193A
- **Special Edition (SE) A-Series:** M1193AN

Step	Applying the Hand/Foot Sensor
1	Position the sensor over the hand or foot with optical components opposite each other.
2	Holding sensor, stretch strap slightly (not more than 2.5cm (1.0")).
3	Insert stretched strap into slot, hold it there while threading end through latch. If strap is too long, thread it through second latch.
4	Plug sensor into instrument (or into adapter cable if needed).
5	Inspect and change the sensor application site periodically.

---

## Applying Ear-Clip Sensors

Adult or pediatric ear, for patients more than 40kg.

- **A-Series:** M1194A, or
- **Special Edition (SE) A-Series:** M1194AN

Step	Applying the Ear-Clip Sensor
1	To improve perfusion, massage or warm the ear lobe.
2	Clip probe onto fleshy part of the lobe. Plastic fixing mechanism minimizes artifact generated by patient motion. Do not position probe on cartilage or pressed against the head.
3	Plug sensor into instrument (or into adapter cable if needed).
4	Inspect and change the sensor application site periodically.

---

## Cleaning and Low-Level Disinfection

The reusable sensors should be cleaned and disinfected, but never sterilized. Follow the procedure outlined below.

### Warnings

- Use only validated cleaning agents and disinfectants listed below; do not use any others. Failure to do so may damage the sensor or its connecting wires and shorten product lifetime or cause safety hazards.
- Select disinfectants carefully as some have very similar names but completely different compositions.
- Do not immerse the sensor connector in any of the cleaning solutions, disinfectant, or other liquid (only the sensor and cable housing may be immersed, not the connector).
- Do not soak sensors in disinfectants longer than specified by the disinfectant manufacturer.
- Do not sterilize the sensors.

### Validated Cleaning Agents

- Mild Detergent
- Salt Solution (1%)

### Validated Disinfectants

- |                                |                              |  |
|--------------------------------|------------------------------|--|
| Metricide <sup>®</sup> 28      | Cidex <sup>®</sup> Formula 7 | Kohrsolin <sup>®</sup> (2%)                    |
| Metricide <sup>®</sup> Plus 30 | Cidex <sup>®</sup> OPA       | Mucapur <sup>®</sup> -CD (1%)                  |
| Terralin <sup>®</sup> Liquid   | Cidex <sup>®</sup> Plus      | Isopropanol (70%) or<br>Isopropanol Wipe (70%) |
| Incidin <sup>®</sup> Liquid    | Omnicide <sup>®</sup> 28     |  |

Step	Cleaning and Low-Level Disinfection
1	Clean sensor following instructions supplied with cleaning agent.
2	Disinfect sensor following instructions supplied with disinfectant.
3	Rinse sensor in water, wiping it dry with a clean cloth, leaving to dry completely. If you notice any deterioration or damage to sensor or cable, dispose the sensor immediately.

---

## Specifications

### Sensor Accuracy

**A-Series Sensors:** For A-Series sensor SpO<sub>2</sub> and Pulse Rate accuracy specifications, reference your Philips/Agilent/HP instrument's Instructions for Use.

**B-Series Sensors:** SpO<sub>2</sub> and Pulse Rate accuracy specifications for M1191B and M1191BL sensors are the same as indicated for M1191A and M1191AL sensors in your Philips/Agilent/HP instrument's Instructions for Use.

**Special Edition A-Series Sensors (AN-Series):** SpO<sub>2</sub> accuracy (functional saturation) for Special Edition A-Series sensors used in combination with M1020B or M3001A, Option A02, specified as the root-mean-square (RMS) difference between the measured values and the reference values, in the range 70 to 100% is: 2% for M1191ANL and M1192AN; 3% for M1193AN, M1194AN, and M1195AN.

**Neonatal Applications:** The claimed accuracy numbers for neonatal application are already increased by an additional 1% to account for affect on oximeter measurements due to fetal hemoglobin in neonatal blood, as predicted by the literature.

### Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in these sensors are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

### Measurement Validation

The SpO<sub>2</sub> accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO<sub>2</sub> were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to black

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements can be expected to fall within the  $\pm$  Arms value measured by a CO-oximeter. Functional testers, such as an SpO<sub>2</sub> simulator cannot be used to assess the accuracy of pulse oximeter Sensors.

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For more information, please call your local Philips sales office listed in your telephone directory or a Philips regional office listed below for the location of your nearest sales office.

**United States:**

Philips Medical Systems  
North America Corporation  
3000 Minuteman Road  
Andover, MA 01810-1099  
(800) 225-0230

**Canada:**

Philips Medical Systems Canada  
281 Hillmount Road  
Markham, ON  
L6C 2S3  
(800) 291-6743

**Latin America Headquarters:**

Philips Medical Systems  
1550 Sawgrass Corporate Parkway #300  
Sunrise, FL 33323 USA  
Tel: (954) 835-2600  
Fax: (954) 835-2626

**Manufacturing Address:**

Philips Medizin Systeme  
Böblingen GmbH  
Hewlett-Packard Str. 2  
71034 Böblingen, Germany

**Europe, Middle East and Africa:**

Philips Medizin Systeme  
Böblingen GmbH  
Cardiac and Monitoring Systems  
Hewlett-Packard Str. 2  
71034 Böblingen, Germany  
Fax: (+49) 7031 463 1552

**Asia Pacific Headquarters:**

Philips Medical Systems  
30/F Hopewell Centre  
17 Kennedy Road  
Wanchai  
Hong Kong  
Tel: (852) 2821 5888  
Fax: (852) 2527 6727

<http://shop.medical.philips.com>

 0366 This product complies with the requirements of the Council Directive 93/42/EEC June 1993 (Medical Device Directive).

 Contains no Latex

Subject to modification

© 2006

Koninklijke Philips Electronics N.V.  
Printed in Germany December 2006  
Philips Part Number 453564052251

## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>