LER 2900 QUANTUM® BIG SCREEN MICROSCOPE


Instruction Guide

The Quantum® Big Screen Microscope is designed for use by children age 6 years and older. Its unique design invites young scientists to take a closer look at the wonders of the microscopic world. Using the Quantum Big Screen Microscope, children can focus in on objects they collect, or slides that have been prepared for them. The Quantum is an excellent tool to spark children's curiosity and foster an early appreciation for science. We hope your children enjoy their experiences.

Features and Accessories (Figure 1)

The Quantum® Big Screen Microscope includes these Features:



- A. **8.5 Inch/22 cm Screen** equipped with interior lamp for clear, crisp images
- B. Easy-turn **Focus Knobs**
- C. **Revolving Turret** for selecting **10x, 20x, or 40x** magnification
- D. **Slide Stage** for Learning Resources' Investigator® Slide Strips or prepared slides
- E. **Condensing Lever**
- F. **On/Off Switch**
- G. **Plug-In Outlet** for AC adapter (LER 2901, sold separately)
- H. **Snap-Tight™ Tray** for storing slide equipment

The Quantum Big Screen Microscope comes with these accessories: Four plastic slides, four slide covers, one stirring rod, one pair of tweezers and one spatula.


1) Inserting Batteries (Figure 2)

Requires 2 "C" type batteries (not included). Battery chamber is located at rear-left side of unit. Unscrew the battery chamber knob by turning it counterclockwise. To insert batteries, place each one in the chamber using printed symbols in back for reference. Replace the knob and turn it clockwise until it stops.

2) Changing the Bulb (Figure 3 & 4)

Allow the bulb to cool before replacing it. Unplug the Quantum and remove the faceplate and condensing lever; see Figure 3. Grasp the light module and turn to the left to unlock it. Lift up and the light module should come out with ease. Replace the light module with the one included in the set. Insert the new light module that is included in the set. Twist to the right to lock it in; see figure 4.

3) Making a Prepared Slide (Figure 5)

Prepared slides are used to view specimens you find. A specimen must be very thin so light will pass through it, such as hair, tissue paper or a leaf. If necessary, have an adult help slice a very thin piece from your specimen. If a specimen is too thick, it cannot be viewed because light will not pass through the specimen. Specimens should be no larger than a 3/8 inch or 1 cm square.

4) Mounting a Specimen onto a Slide

Wipe the plastic slide and slide cover clean. Pick up the specimen with the tweezers and place it in the center of the plastic slide. Place one or more drops of water on the sample, place the slide cover over it and press down. Be sure not to trap any air bubbles. Remove any excess water with blotting paper. The specimen is now ready for observation.

5) Viewing Slides

The Quantum projects the brightest and clearest image in a dimly lit room. To view specimens at **10x** power, slide the condensing lever to the **LEFT** until it clicks into position. To view specimens under **20x** and **40x** power, slide the lever to the **RIGHT**.

Turn on the lamp. Set the magnification to **10x**. Make sure the condensing lever is slid to the left position. Place the slide on the stage and secure it in place under the clips. Turn the focusing knobs until the image is clear. Choose a higher magnification by turning the revolving turret to the desired magnification. Don't forget to readjust the condensing lever to the right position when viewing specimens at **20x** or **40x**.

When viewing Investigator® Slide Strips, use the 10x magnification and position the condensing lever to the left position. Place a strip in upside down. Turn the focusing knobs until the image is clear. Please remember that some of these images have already been magnified tremendously. Further magnification may reduce image clarity.

6) Using the A/C adapter, LER 2901 (Figure 1)

Make sure Quantum is switched OFF. Insert adapter into plug on the right side of the Quantum. Plug adapter body into a household electrical socket. The Quantum is now ready to use.

Note: The light is brightest with the adapter; however, light bulb life is longer with batteries.
CAUTION!

1. Use only the recommended A/C adapter described above. The A/C adapter is NOT a toy. Adult supervision is recommended.
2. Examine A/C adapter cord, plug, body and other parts for damage before each use. If damaged, repair or replace adapter before using.

7) Handling Precautions!

1. The lenses are the most important parts of the Quantum. Please take care of them.
2. When the lenses get dusty or dirty, gently wipe with a clean lens tissue or a soft, clean lint free cloth. **DO NOT USE YOUR FINGERS.**
3. Place a plastic bag over your Quantum when it is not in use to protect it from dust.

4. Store your Quantum in a moisture-free place. Moisture on the lamp can reduce the light intensity.

8) CAUTION!

- Use only 2 "C" batteries.
- To remove batteries, see "**Inserting Batteries**". Remove old batteries and replace with new batteries. Be sure to insert batteries correctly according to (+) and (-) symbols on battery chamber.
- Be sure to insert batteries correctly (with adult supervision) and always follow the toy and battery manufacturer's instructions.
- Do not mix old and new batteries.
- Do not mix alkaline, standard (carbon-zinc) or rechargeable (nickel-cadmium) batteries.
- Non-rechargeable batteries are not to be recharged.
- Rechargeable batteries are to be removed from the toy before being charged.
- Rechargeable batteries are only to be recharged under adult supervision.
- Only batteries of the same or equivalent type as recommended are to be used.
- The supply terminals are not to be short-circuited.
- Always remove weak or dead batteries from the product.
- Remove batteries if product will be stored for an extended period of time.
- If using an AC adapter, check cord frequently and do not use if wire is damaged or exposed.

Please retain these instructions for future reference.

© Copyright 1999, Learning Resources, Inc.

Transformers or battery chargers used with the toy are to be regularly examined for damage to the cord, plug enclosure and other parts, and that, in the event of such damage, the toy must not be used with this transformer or battery charger until the damage has been repaired.

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