## ResMed Masks

### Disinfection/Sterilization Guide

English • Deutsch • Français • Italiano • Español • Português • Nederlands

#### **Mirage Vista™**

NASAL MASK

#### **Mirage Activa™**

NASAL MASK

#### **Mirage Kidsta™**

NASAL MASK

#### **Mirage Swift™**

NASAL PILLOWS SYSTEM

#### **Mirage Swift™ II**

NASAL PILLOWS SYSTEM

#### **Ultra Mirage™ II**

NASAL MASK

#### **Ultra Mirage™**

FULL FACE MASK

#### **Ultra Mirage™**

NON-VENTED NASAL MASK

#### **Mirage™**

NON-VENTED FULL FACE MASK SERIES 2

#### **Ultra Mirage™**

NON-VENTED FULL FACE MASK

**RESMED**

## Disinfecting/Sterilizing the Mask or Nasal Pillows System between Patients

This Disinfection/Sterilization Guide is intended for multipatient use of the mask/nasal pillows\* system in a sleep lab, clinic or hospital environment. If you use the mask as a single user in the home, refer to the User Guide for cleaning instructions.

The following are the recommended and validated procedures for cleaning, disinfection and sterilization of the mask. We advise that each healthcare facility consult its own procedures before carrying out the following instructions.

For details on the mask refer to the User Guide. Parts that cannot be disinfected/sterilized are listed in the table on page 5.

### WARNINGS AND CAUTIONS

- ResMed cannot give any assurance that deviations from the procedures listed below, and their effect on the performance of the product, will be acceptable.
- The mask components should not be subjected to autoclave or ethylene-oxide gas sterilization.
- When using detergents, disinfectants or sterilization agents, always follow the manufacturer's instructions.

#### Notes:

- *The headgear does not require disinfection or sterilization. However, it is necessary to thoroughly wash the headgear before using between patients. The cleaning procedure is detailed in the User Guide.*
- *The mask must be disassembled and cleaned prior to disinfection or sterilization. The mask must also be inspected after disinfection/sterilization. See the section "Inspection after Disinfection/Sterilization" for inspection details.*

---

\* References to mask/s in this guide refer to ResMed masks and nasal pillows system.

## Clinical Use Only

After disassembly, thoroughly disinfect or sterilize the mask components by following the procedure for *one* of these methods:

- Chemical Disinfection
- Thermal Disinfection
- STERRAD™ Sterilization.

### **Option 1: Chemical Disinfection**

#### **Pre-disinfection**

- 1 Disassemble the mask.
- 2 Clean the individual parts of the mask with a soft bristle brush for one minute while soaking them in detergent (eg, Aniosyme™ DDI or Alconox™). Pay particular attention to all crevices and cavities.
- 3 Rinse the components twice by agitating them vigorously in drinking quality water (5 litres per mask) and allow them to air dry out of direct sunlight.

*ResMed has tested the detergents Aniosyme DDI (diluted at 0.5% for 15 minutes) and Alconox (diluted at 1%).*

#### **Disinfection**

- 1 Soak the disinfectable mask components in a commercially available solution of peracetic acid 0.15% (eg, Anioxyde™1000), ortho-phthalaldehyde 0.55% (eg, CIDEX™ OPA) or glutaraldehyde 3.4% (eg, CIDEX Plus™).
- 2 Rinse the mask components in drinking quality water (5 litres per mask) and allow them to air dry out of direct sunlight.

*Testing by ResMed has shown that the mask components can withstand 15 cycles using solutions of peracetic acid (0.15% Anioxyde 1000 for 30 minutes), ortho-phthalaldehyde (0.55% CIDEX OPA for 20 minutes) or glutaraldehyde (3.4% CIDEX Plus for 20 minutes).*

**Note:** *If using the detergents and disinfectants listed above, the following combinations should be followed:*

- Detergent Aniosyme DDI with Anioxyde 1000 as the disinfectant  
OR
- Detergent Alconox with disinfectant CIDEX OPA or CIDEX Plus.

## Option 2: Thermal Disinfection

### Pre-disinfection

- 1 Disassemble the mask.
- 2 Clean the individual parts of the mask with a soft bristle brush for one minute while soaking them in detergent (eg, Neodisher™ Medizyme or Alconox). Pay particular attention to all crevices and cavities.
- 3 Rinse the components twice by agitating them vigorously in drinking quality water and allow them to air dry out of direct sunlight.

*ResMed has tested the detergents Neodisher Medizyme (concentration at 2% for 30 minutes) and Alconox (diluted at 1%).*

### Disinfection

Using a certified hot water disinfection system, soak the disinfectable mask components using **one** of the temperature-time combinations listed in the following sections.

#### European Standards

Based on prEN ISO 15883-1, the mask has been tested at the following temperature-time combinations:

- 70°C for 100 minutes
- 75°C for 30 minutes
- 80°C for 10 minutes
- 90°C for 1 minute.

#### Regional Requirements

Due to specific regional requirements, the mask has also been tested for disinfection at the following temperature-time combination:

- 93°C for 10 minutes.

*Testing by ResMed has shown that the mask components can withstand 15 cycles of thermal disinfection using the above temperature-time combinations. The device used for testing was Branson™ LTH 1828-24 hot water bath.*

## Clinical Use Only

### **Option 3: STERRAD Sterilization System**

**Note:** *The Ultra Mirage™ II Nasal Mask is not validated for the STERRAD Sterilization System.*

#### **Pre-sterilization**

- 1 Disassemble the mask.
- 2 Clean the individual parts of the mask with a soft bristle brush for one minute while soaking them in detergent (eg, Alconox). Pay particular attention to all crevices and cavities.
- 3 Rinse the components twice by agitating them vigorously in drinking quality water and allow them to air dry out of direct sunlight.

*ResMed has tested the detergent Alconox (diluted at 1%).*

#### **Sterilization**

Sterilize the sterilizable mask components as described in the manufacturer's instructions for the STERRAD Sterilization System (eg, Sterrad 100S or Sterrad NX).

*Testing by ResMed has shown that the mask components can withstand 15 cycles of the STERRAD 100S Sterilization System, when using the tray method.*

**Note:** *The use of pouches is not recommended.*

#### **Inspection after Disinfection/Sterilization**

Perform a visual inspection of each mask component. If any visible deterioration of a mask component is apparent (cracking, crazing, tears etc), the mask component should be discarded and replaced. Slight discoloration of the silicone components may occur and is acceptable.

**Note:** *If a healthcare facility requires an additional disinfection or sterilization process after reassembly, the number of validated cycles is halved.*

#### **Parts not Validated for Multipatient Use**

Some parts of the masks have not been validated for multipatient use. These parts need to be replaced with new parts when using between patients. Please check the Components Card on the website [www.resmed.com](http://www.resmed.com) for a list of available replacement parts for each mask system.

## Clinical Use Only

The parts that cannot be disinfected or sterilized are listed in the following table.

Masks/Pillows System	Parts that cannot be disinfected or sterilized
<ul style="list-style-type: none"> <li>• Mirage Vista™ Nasal Mask</li> <li>• Mirage Activa™ Nasal Mask</li> <li>• Mirage Kidsta™ Nasal Mask</li> </ul>	Swivel, short tubing
<ul style="list-style-type: none"> <li>• Mirage Swift™ Nasal Pillows System</li> <li>• Mirage Swift™ II Nasal Pillows System</li> </ul>	Short tube assembly
<ul style="list-style-type: none"> <li>• Mirage™ Full Face Mask Series 2</li> <li>• Ultra Mirage™ Full Face Mask</li> </ul>	Valve membrane
<ul style="list-style-type: none"> <li>• Non-vented Ultra Mirage™ Nasal Mask</li> <li>• Non-vented Ultra Mirage™ Full Face Mask</li> <li>• Ultra Mirage™ II Nasal Mask</li> <li>• Non-vented Mirage™ Full Face Mask Series 2</li> </ul>	None

**Note:** *Parts that cannot be disinfected or sterilized should be replaced when using between patients.*

---

**Note:** *Different disinfectants and detergents are available in each country. For this reason, ResMed has identified alternatives as examples and has validated the disinfectants and detergents listed in this guide.*

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**ResMed Ltd** (Manufacturer) 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia  
**ResMed Corp** (US Designated Agent) 14040 Danielson Street Poway CA 92064-6857 USA  
**ResMed (UK) Ltd** (EU Authorized Representative) 65 Milton Park Abingdon Oxfordshire OX14 4RX UK  
**ResMed Offices** Australia, Austria, Brazil, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, UK, USA (see [www.resmed.com](http://www.resmed.com) for contact details).

Mirage, Mirage Activa, Mirage Kidsta, Mirage Swift, Mirage Vista and Ultra Mirage are trademarks of ResMed Ltd and Mirage, Activa, Kidsta, Swift and Vista are Registered in U.S. Patent and Trademark Office.

CIDEX and STERRAD are trademarks of Johnson & Johnson. Alconox is a trademark of Alconox Inc. Neodisher Medizyme is a trademark of Chemische Fabrik Dr Weigert. Aniosyme and Anioxyde are trademarks of Laboratoires Anios. Branson is a trademark of Branson Ultrasonics.

## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>