Infrared Ear Thermometers iCcurate User Manual


Thank you for purchasing a quality Rosewill Product.

Please register your product at : www.rosewill.com for complete warranty information and future support for your product.

English

Operating

Product Description

Introduction

The eardrum is an ideal location for measuring the core temperature of the human body since it is close to, and shares blood vessels with the hypothalamus, the temperature control center of the brain. Medicare iCcurate is a compact and accurate infrared ear thermometer that detects the strength of thermal radiation emitted by the eardrum. The unique ergonomic design promotes ease of use.

Taking temperature with Medicare infrared ear thermometer is fast and safe.

Important Information : Please carefully read this instruction manual before using the product.

This product is designed to measure human body temperature in the ear. Do not use it for any other purpose.

This product is intended for home use only.

Use only genuine Medicare probe covers, models RPC-902(40 – pack).

Use of this thermometer is not intended as a substitute for consultation with your physician.

Do not use the thermometer if it malfunctions or has been damaged in any matter.

When the ambient temperature of the thermometer changes too much, such as moving the Thermometer from one place of lower temperature to another place of higher temperature, Allow the thermometer to remain in a room for 30 minutes where the temperature is between 10°C to 40°C (50°F- 104°F).

A warning to remove primary batteries if equipment is not likely to be used for some time.

How to Operate

1. Should make sure the probe cover is well attached, intact and undamaged before taking each measurement.
2. Turn on the thermometer by pressing the ON/M button. The thermometer will perform self-test with all segments displayed for 1 second. When the “READY” symbol appears, the thermometer is ready for use.
3. Stabilize the head of the person whose temperature is to be taken. To avoid injury and allow for easy use, the head must remain steady during the temperature taking process.
4. Perform an ear-tug to straighten the ear canal. For children under 1 year old, pull the ear straight back. For adults or children over 1 year old, pull the ear up and back.

5. Gently and carefully place the thermometer probe tip snugly into the ear canal while tugging the ear.
6. Press and hold the ACTIVATION button. A short beep will be heard to confirm the end of measurement.
7. Remove the temperature from the ear, and get the temperature reading on the LCD.
8. Press the BACKLIGHT button and hold it for easy nighttime reading, if desired.
9. The device is ready to take a new temperature only when the "READY" symbol appears steady.
10. Replace with a new probe cover after each measurement.
11. The thermometer will automatically turn off if no button is pressed within 60 seconds.

Recalling and Erasing Memory Data

The last temperature taken before the thermometer powers down is stored in memory. Up to 8

Temperature readings can be stored in memory. Up to 8 temperature readings can be stored into

The memory cells from M1 to M8. The latest reading is always at M1.

By pressing ON/M button, the cell number will be shown on the LCD.

After releasing the button, the stored temperature will be displayed with the memory mode Symbol M indicated.

Press and hold the ON/M button for at least 5 seconds to erase all the stored temperatures in Memory. "CLR" will be shown on the LCD.

An empty memory cell shows "---°C" or "---°F"

Changing the Temperature Scale

1. Wait until the thermometer powers down or make sure the thermometer is at off state.
2. Press and hold the ACTIVATION button first, and then press and release the "ON/M" button.
3. Wait until the display scale changes from °C to °F or vice versa. All the memory data will be changed accordingly.
4. Release the ACTIVATION button.

Important Notes about Taking Temperature

1. It is necessary to inspect the soft lens for damage each time the probe cover is changed. Do not use the thermometer if the soft lens is damaged.
2. It is important to gently position the probe into the ear canal snugly and aim the thermometer probe at the eardrum to obtain an accurate reading.
3. The same person may have different temperatures in the right and left ears. It is recommended that all temperatures be taken in the same ear when monitoring an individual's temperature.
4. Many people do not run a fever (a rise in body temperature above normal) even when they are ill. You must immediately consult a physician, regardless of temperature reading, if you feel an illness is present.
5. Under any of the following situations, you should take three consecutive measurements in the same ear and use the highest temperature of the three readings:
 - a. Newborn babies under ninety (90) days old.
 - b. Infants under three (3) years of age with immune disorders which are life-threatening at abnormal body temperatures.
 - c. You are unfamiliar with the temperature measuring process and notice a wide variance between consecutive readings.
6. When taking consecutive measurements, a variation of 0.2°C - 0.3°C (0.3°F - 0.5°F) can occur for any type of thermometer.
7. If an individual has been sleeping on one ear for a while, the temperature of that ear will be higher than normal. You should measure the other ear or wait for a few minutes for the ear temperature to normalize before measuring.
8. Excessive earwax will adversely affect the thermometer's accuracy. Making sure the ears are clean before taking a temperature measurement

should result in a more accurate reading.

9. The thermometer is designed to take accurate temperatures with a probe cover in place. Temperatures taken without a probe cover in place.

Temperatures taken without a probe cover attached will not be accurate, and can result in damage to the soft lens. Change the probe cover after each use to avoid infection and cross contamination.

Troubleshooting

MESSAGE	SITUATION	SOLUTION
E1	Operating temperature exceeds the range of 10°C - 40°C (50°F - 104°F)	Move to a room within the operating range wait 30 minutes before taking temperature
E2	Ambient temperature changes too fast	Wait until the ambient temperature is stable
E3	System Error	Contact the retailer or service center
HI LO	Temperature taken is not within Typical human temperature range. (34°C - 42.2°C or 93.2°F - 108°F)	Make sure a new and clean probe cover is attached to the thermometer Probe. Then, take a new temperature
	Battery capacity is too low. Taking Temperature is not allowed.	Install a new battery
	(1) Power is off. (2) Improper battery installation. (3) The battery is exhausted. (4) Display remains blank.	(1) Press ON/M button again. (2) Check the battery polarity. (3) Replace with a new battery. (4) Contact the retailer or service center.

Cleaning, Care and Storage

The probe tip is the most fragile portion of the thermometer and the soft lens is very delicate.

It is very important to always keep the probe cover attached to protect the soft lens from dirt and damage.

Whenever replacing the probe cover, inspect the soft lens and make sure it is clean and intact.

You may use cotton swab dampened with isopropyl (rubbing alcohol) or warm water to gently clean the soft lens, if dirty.

Use a clean, soft cloth to clean the surface of the device and LCD. Do not use solvents or immerse the device into water or other liquids.

Always keep the thermometer and probe covers within the storage temperature range (- 20°C to 50°C or - 4°F to 122°F) and humidity range (15% - 95% non-condensing).

It is recommended to store the thermometer and probe covers in a dry location free from dust.


When not in use, store the ear thermometer inside the carry case with a new probe cover attached to the probe.

Keep the thermometer, probe covers and battery out of reach of children.

When the ambient temperature of the thermometer changes too much, such as moving the thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 10°C to 40°C (50°F to 104°F).

Replacing the Battery

1. By using screw driver carefully to release battery cover. Before Changing the battery be sure the system is already power off.
2. Press the battery holder lock to release the battery.
3. Remove the battery and replace with a new one, type CR2032, Make sure the positive pole of the battery facing up.
4. Slide the battery cover back in until it snaps in place.
Do not dispose of used batteries in household waste. Take them to special local collection sites.
5. In case, if system is latched up after changing battery. You may not follow up the process of rule one. Just take off battery, waiting 30 sec, then load battery again.

Battery may explode if mistreated. Do not recharge, disassemble or dispose of in fire.

Product Specifications

Displayed temperature range: 34°C to 42.2°C (93.2°F to 108°F)

Operating ambient temperature range: 10°C to 40°C (50°F to 104°F)

Display Resolution: 0.1°C or 0.1°F

Storage & transportation temperature:

-20°C to 50°C (-4°F~122°F),85%RH,non-condensing for one month.

Not AP or APG equipment.

This product is IPX0 grade.

The manufacturer recommends to verify the accuracy by Rosewill Medicare's Customer Department

Every 2 years.

Accuracy:

Measurement Range	Error
<35.5 °C(89.6 °F) or > 42 °C(107.6 °F)	±0.3°C (±0.5°F)
35.5 to 42 °C(86.9 to 107.6 °F)	±0.2°C (±0.4°F)

Accuracy is determined under laboratory conditions in the calibration mode with a typical probe Cover installed.

Storage temperature ranges: -20°C to 50°C (-4°F to 122°F)

Storage humidity (max): 95%RH, non-condensing

Display Mode: Oral (The actual measured ear canal temperature plus mathematical adjustment to approximate the familiar oral range. However, this is not necessarily the same as an oral temperature measured at the same time).

EL back light (DX version only).

Data memory: Up to 8 temperature readings can be stored.


Power supply: One lithium battery, type CR2032

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