
ES - MANUAL DE INSTRUCCIONES **HU** - HASZNÁLATI ÚTMUTATÓ
PT - MANUAL DE INSTRUÇÕES **CZ** - NÁVOD K POUŽITÍ
EN - INSTRUCTIONS FOR USE **SK** - NÁVOD NA POUŽITIE
FR - MANUEL D'UTILISATION **PL** - INSTRUKCJA OBSŁUGI
DE - GEBRAUCHSANWEISUNG **BG** - ИНСТРУКЦИЯ ЗА УПОТРЕБА
IT - MANUALE DI ISTRUZIONI **RU** - ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ
EL - ΕΝΤΥΠΟ ΟΔΗΓΙΩΝ

**TENSIÓMETRO AUTOMÁTICO / AUTOMATIC TENSOMETER / TENSIOMÈTRE AUTOMATIQUE /
 AUTOMATISCHER BLUTDRUCKMESSER / TENSIOMETRO AUTOMATICO / ΑΥΤΟΜΑΤΟ
 ΠΙΕΣΟΜΕΤΡΟ / AUTOMATA VÉRNYOMÁSMÉRŐ / AUTOMATICKÝ MĚŘIČ KREVNÍHO TLAKU /
 AUTOMATICKÝ MERAČ KRVNÉHO TLAKU / CIŚNIENIOMIERZ AUTOMATYCZNY /
 АВТОМАТИЧЕН АПАРАТ ЗА ИЗМЕРВАНЕ НА КРЪВНОТО НАЛЯГАНЕ /
 АВТОМАТИЧЕСКИЙ ТОНОМЕТР**

MOD.: **TSB-10**

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
1. INTRODUCTION

This instrument is a fully automatic digital tensometer incorporating a method of oscillometric measurement for the swift, reliable measurement of systolic, diastolic arterial pressure and heartbeat frequency.

Read the instructions for use carefully before using the device and keep for future reference. Contact your doctor if you have any further queries about arterial pressure and measuring the same.

2. VERY IMPORTANT INFORMATION

- Do not forget: **self-measuring means control, not a diagnosis or treatment.** Unusual values should, at all times, be discussed with your doctor. Under no circumstances should the doses for each medication given by your doctor be changed.
- This device is designed exclusively for household use. Keep out of reach of children.
- The pulse beat indicator is not suitable for measuring the frequency of a pacemaker.
- In the event of heartbeat irregularity (arrhythmia), measurements taken with this instrument should only be evaluated with medical consultation.
- There are many different causes for the appearance of high arterial pressure values. For this reason, always consult your doctor since certain illnesses may cause deviations in the measurement of arterial pressure. **Never use the results of tests to change medically prescribed treatment on your own account.**

Electromagnetic interference:

This device contains sensitive electronic components (microcomputer). For this reason, keep this device away from electrical or strong electromagnetic fields, e.g., mobile telephones, microwaves may lead to the temporary loss of precise measurements.

3. STANDARD VALUES

Classification table of the arterial pressure values measured in a relaxed state, in accordance with the World Health Organization (WHO):

Values	Systolic arterial pressure (SYS)	Diastolic arterial pressure (DIA)
Hypotension	Below 100	Below 60
Normal values	Between 100 and 140	Between 60 and 90
Moderate hypertension	Between 140 and 160	Between 90 and 100
Moderately serious hypertension	Between 160 and 180	Between 100 and 110
Serious hypertension	Over 180	Over 110

4. PRODUCT DESCRIPTION

1. Connecting the armband
2. Liquid crystal display
3. O/I button and Memory
4. Batter compartment cover
5. Armband
6. Armband tube
7. AC/DC plug

5. INSERTING THE BATTERIES

Once the device has been unpacked, first insert the batteries. The battery compartment is located underneath the device.

- a) Remove the battery compartment cover, as shown in the figure (Fig. 1).
- b) Insert the 4 alkaline LR 06 (AAA) 2x1.5V type batteries, ensuring correct polarity (Fig. 2).
- c) The  icon on the screen indicates that the batteries are used and need replacing (Fig. 3).

Warning:

- When the  icon appears indicating that the batteries are used, **the device blocks until they are replaced**. All 4 batteries must be replaced.
- It is not advisable to use rechargeable batteries.
- If the device is not going to be used for a long period of time, remove the batteries from the compartment.
-  Used batteries should be deposited in special containers for the purpose.

Checking working order: Once the batteries are fitted in place, keep pressing the O/I button to check all the elements on the screen. If working correctly, all the segments should be displayed.

6. USE OF AC/DC POWER SUPPLY (not supplied)

The tensometer may be used with a power supply (6VDC/600 mA output with DIN connector). The supply used must comply with the appropriate legal requirements (EC marked on the label).

- a) Insert the DIN connector into the socket located on the left side of the device (Fig. 4).
- b) Plug in the transformer to an electrical 230 V socket.
- c) Check for power by activating the O/I key.

Warning:

- When the device is connected to the transformer, the batteries are automatically turned off.
- In the event of the power supply being cut off while measuring, e.g., if the power supply transformer is disconnected, the device should be "restarted". To do so, pull out the power connector and reinsert once the power supply is re-established.
- In the event of any queries about the transformer, consult a specialized seller.

7. CONNECTING THE TUBE

Insert the armband tube in the plug located on the left of the device, as shown in the illustration (Fig. 5).

8. BASIC RULES FOR MEASURING

- Avoid eating, smoking or any form of exercise prior to measuring. All these factors influence the result of the measurement.
- Take time to relax seated in an armchair, in a calm atmosphere, for ten minutes prior to measuring.
- While measuring, remain seated, calm and without moving or speaking.
- Remove any clothing that might put pressure on the arm.
- Always take the measurement on the same arm (normally the left) and, as much as possible, avoid moving it while measuring. If necessary, rest your arm on a cushion to keep it relaxed.
- Endeavour to take the measurements on a regular basis, at the same time of day and in the same conditions since arterial pressure changes in the course of the day.
- Do not take a measurement immediately after another. Wait for a few minutes, in a relaxed position, seated or lying down, before repeating measuring.
- If the armband is placed too far below or too far above the heart, the measurement will be altered showing a higher or lower pressure relative to the real values (for each 15 cm difference in height, the result of the measuring will alter by 10 mmHg).
- A poorly fitted armband or an air pocket along the sides leads to false measurements.

N.B.


To obtain comparable blood pressure measurement values, it is essential for them to be taken in identical conditions.

9. FITTING THE ARMBAND


- a) Slip the end of the armband (with the built-in rubber stop) through the metal arc to form a loop. The Velcro fastener is on the outside (if the armband is ready, as described above, overlook this step).

- b) Place the armband on the left arm so that the tube is facing the forearm.

- c) Place the armband on the arm, as shown in the figure, being careful to ensure that the lower edge of the armband is 2-3 cm over the elbow, with the end of the armband rubber on the inside of the arm.


- d) Hold out the free end of the armband and seal it with the Velcro fastening.


- e) Leave no gap between the arm and the armband, as that would change the reading. Avoid compressing the arm with any article of clothing (e.g., a jersey), in which case, it should be removed.

- f) Ensure that the armband with the Velcro fastener adheres comfortably to the arm but without being too tight. Hold out the arm over the table (the palm of the hand should be face up) so that the armband is at the same level as the heart. Be careful not to bend the tube.


- g) Keep calmly seated for two minutes before starting measuring.

N.B.

If it is not possible to place the armband on the left arm, it may be placed on the right arm. It is important for the measuring always to be taken on the same arm.

10. MEASURING


- a) Press the O/I button and the pump starts to inflate the armband. The viewer will indicate if the armband pressure is increasing.
- b) Once inflated, the pump stops and the pressure gradually decreases. The armband pressure and a control value are displayed.
- c) While checking the heart beat, the heart symbol starts to flash on the viewer and a "bip" is heard for each beat.
- d) As soon as measuring is completed, a long "bip" sounds out. The viewer will indicate the systolic and diastolic pressure as well as the patient's heart beat frequency.
- e) The result of the measuring will continue to be displayed until the device is turned off. If none of the buttons are pressed within 5 minutes, the device disconnects automatically to save battery power.

N.B.: If for any reason it is necessary to stop measuring blood pressure (e.g., if the patient is unwell), simply press the O/I button. The armband pressure will fall immediately.

11. MEMORY

The measurement data are memorized by the device until the following measuring operation or until the batteries are replaced. To recuperate the values for the last measuring, simply press the O/I button and keep pressing for a few seconds.


12. CLEANING AND CARE

- Clean the device with a slightly damp soft cloth. Do not use detergents, petroleum, solvents or similar fluids for cleaning.
- Stains on the armband or on the rubber can be removed with a damp cloth. Do not wash the armband.
- Do not expose the armband to extreme temperatures, to damp, dust or to direct sunlight.
- Avoid dropping, do not knock and treat with care. Avoid strong vibrations.
- Avoid folding over the armband and tube excessively.
- **Do not open the device**, otherwise the manufacture set calibration becomes invalid.

13. ERROR MESSAGES

If an error occurs while measuring, the operation is interrupted and the relevant error message is displayed.

Error n.	Possible cause(s)
ERR 1	The systolic and diastolic pressure have been determined and then the armband pressure has fallen below 20 mmHg. This situation can occur, e.g., when, after measuring the systolic, the armband tube has been disconnected. Other possible causes: it has not been possible to detect any beat.
ERR 2	Abnormal pressure impulses compromise the result of the measuring. Causes: the arm was moving during the measuring.
ERR 3	The armband has taken too long to inflate. The armband is incorrectly positioned or perhaps the tube connection is not fully airtight.
ERR 5	The readings measured indicate an unacceptable difference between the systolic and diastolic pressures. Take another reading, paying careful attention to the instructions. Contact your doctor if unusual readings continue.

14. OTHER POSSIBLE PROBLEMS AND THEIR SOLUTION

Faulty working order	Solution
The screen remains blank when connecting the device, despite having batteries.	<ol style="list-style-type: none"> 1. Ensure that the batteries are fitted according to the polarity and, if necessary, insert them correctly. 2. If the display data is unusual, reinsert the batteries or replace them.
No pressure despite the pump working.	<ul style="list-style-type: none"> • Check the armband tube connection and if it is necessary to connect it correctly.
The device frequently fails when measuring arterial pressure, or the values measured are too high/low.	<ol style="list-style-type: none"> 1. Check the position of the tensometer. 2. Measure the arterial pressure again in a state of rest and calm, observing the details indicated in "Basic Rules".
The results are different with each measurement despite the correct working order of the device and it indicating normal values.	<ul style="list-style-type: none"> • Pay attention to the following "Additional information" and to the points listed in common "Error Messages", then repeat measurement.
The values measured with this device differ from those taken by the doctor.	<ul style="list-style-type: none"> • Record the daily evolution of the values and consult your doctor.

15. ADDITIONAL INFORMATION

The level of arterial pressure is subject to fluctuations, even in healthy people. It is important to bear in mind that **comparable measurements always need the same conditions (of calm)**. If, despite observing all these factors, the fluctuations are higher than 15 mmHg and/or irregular beat tones are occasionally heard, consult your doctor.

Consult a specialized seller if you have problems with the device. **Do not attempt to fix the device yourself.**

Unauthorized opening of the instrument will leave any claim under the guarantee null and void.

16. RULES FOR REFERENCE

Device standard:	This device complies with the European standard requirements for non-invasive arterial pressure control instruments EN1060-1 / 12:95 EN1060-3 / 09:97 DIN 58130, NIBP – clinical research ANSI / AAMI SP10, NIBP – requirements
Electromagnetic compatibility:	This device complies with the European standard specifications EN 60601-1-2
Clinical test:	The clinical function test was conducted in Germany, in accordance with the DIN 58130 / 1997 N6 (sequential) procedure.

In compliance with the specifications in the EU 93/42/CEE directive for Class IIa medical products.

17. TECHNICAL SPECIFICATIONS

Storage temperature:	-5 °C to +50 °C
Humidity:	15 - 85% maximum relative humidity
Working temperature:	10 °C to 40 °C
Screen:	Liquid crystal display
Method of measurement:	Oscillometric
Measurement interval:	
SYS / DIA:	30 to 280 mmHg ± 3 mmHg
Beat:	40 to 200 beats per minute ± 5% of the reading
Tensometer pressure indication limits:	From 0 to 299 mmHg. As from 300 mmHg, "HI" is displayed.
Measurement resolution:	1 mmHg
Power supply:	4 x LR 06 type (AA) 1.5V batteries (supplied with the device).

18. WASTE ELIMINATION

The materials used in the packaging, product and accessories can be recycled. Correct separation of waste materials will allow recyclable materials to be re-used.

When the appliance is no longer useful, you should cut its cord to put it out of action and dispose of it using a suitable waste elimination procedure. For further information on this matter, consult your local authorities.

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