
Health o meter®

Your Measure for Good Health™

AUTOMATIC WRIST
**Blood Pressure
Monitor**

Model 7633

IMPORTANT NOTES

Download from www.Somanuals.com All Manuals Search And Download.
Please read this entire manual carefully before using this device.

Healthometer®

Automatic Wrist Blood Pressure Monitor

IMPORTANT SAFETY INSTRUCTIONS

DO NOT DESTROY

Features:

- **One-Touch Easy Operation**
- **Large, Easy-To-Read Digital Display**
- **90 Automatic Memory Feature**
- **MicroStep™ Technology Precisely Controls Deflation for Greater Accuracy**
- **Auto Shut-Off**
- **Convenient Storage Case**
- **Compact and Portable**

In accordance with International Standards:

CE ISO
0366 9001

Table of Contents

Introduction	4
Quick Reference Guide	5
Standards of Quality	6
About the Wrist Monitor	6
About Blood Pressure	6-7
Blood Pressure Standard	7
Blood Pressure Fluctuation	7
Carrying Case	8
Product Component Identification	8
Display Explanations	9
Installing Batteries	9-10
HOW TO USE	10-14
• Helpful Tips	10
• Applying the Cuff	11
• Proper Posture During Measurement	11-12
• Taking a Measurement	12-13
• Recalling Memory	14
Troubleshooting	14
Precautionary Notes	15
Specifications	16
Patents	16
Warranty Information	17
Blood Pressure Log	18-20

Introduction

Congratulations on your purchase of the new Health o meter® Wrist Blood Pressure Monitor Model 7633. This fully automatic, wrist mounted device is both easy to use and ideally suited for daily measurements. The large display includes systolic, diastolic and pulse measurements, all clearly shown at the completion of each reading.

In addition, you can store up to 90 measurements using the monitor's built-in memory feature, ideal for users who want to monitor and track their blood pressure on a regular basis. The Health o meter® blood pressure monitor is compact and portable, making it ideal for both home and travel usage.

Blood pressure measurements determined with the Health o meter® blood pressure monitor are equivalent to those using the cuff/stethoscope auscultation method and are within the limits prescribed by the American National Standards Institute (ANSI). This monitor is to be used by an adult and is intended for in-home use.


Note: Please read this manual carefully before use. For specific information about your own blood pressure, contact your physician. Please be sure to keep this manual for future reference.

Quick Reference

Please see the following pages for complete instructions.

A. Installing Batteries

D. Press ON/OFF/START button

Standards of Quality

This Blood Pressure Monitor conforms with the American National Standards Institute/American Association of Medical Instrumentation Standards for Non-Invasive Blood Pressure Monitors. In addition, the device complies with the European regulations and bears the CE mark "CE 0366". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC dated 14 June 1993 concerning medical devices, as well as the EMC directive 89/336/EEC:

ANSI/AAMI SB10	Electronic or automated Sphygmomanometers
EN 1060-1 (Dec. 1995)	Non-invasive blood pressure measuring equipment General requirements
EN 1060-3 (Feb. 1997)	Non-invasive blood pressure measuring equipment Supplementary requirements for electromechanical blood pressure measuring systems

About the Wrist Blood Pressure Monitor

Your Health o meter® blood pressure monitor uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure, which is equivalent to the air pressure. The measurement of your blood pressure is based on this baseline pressure. Once the cuff inflates to block the passage of blood in the artery, the deflation process begins. During the deflation of the wrist cuff, the unit is detecting the pressure oscillations generated between beats of the pulse. Any muscle movement during this period of time will cause measurement error. After detecting the amplitude and the slope of the pressure oscillations during the deflation process, your Health o meter® monitor will determine the systolic and diastolic pressures for you, your pulse rate is detected and displayed at the same time.

ABOUT BLOOD PRESSURE

What is Blood Pressure?


Blood pressure is the pressure exerted on the inside of the artery wall while blood flows through the artery. The pressure measured when the heart contracts and sends blood out of the heart is the "systolic" reading (higher number). The pressure measured when the heart dilates with blood flowing back into the heart is called the "diastolic" (lowest number) blood pressure.

Why Measure Your Blood Pressure?

Various health problems associated to high blood pressure are quite common with people in today's society. There is a significant correlation of high blood pressure with cardiovascular disease. This correlation has made measuring blood pressure a necessity for those individuals identified with high blood pressure, placing them into a high risk category pertaining to cardiovascular disease. Also, as a person ages, one's blood pressure can change over time, thus it makes good "health sense" to monitor your blood pressure on a regular basis.

Blood Pressure Standard

The World Health Organization (WHO) has developed a blood pressure standard, according to which areas of low and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups, etc.


It is important that you consult with your physician regularly. Your physician can tell you your optimize blood pressure range.

Blood Pressure Fluctuation


Blood pressure fluctuates all the time!

Your blood pressure will change over the course of the month, week and even throughout the day. It is also influenced by seasonal and temperature changes.


Carrying

Please store your
as shown below:

Display Explanations

Display:

Systolic Pressure

Diastolic Pressure

Pulse Rate

No. of Memory Sequence

Icons:

mem.

Appears when measurement value is stored in memory or is recalled from memory

Appears when batteries should be replaced


Shows the Pulse Rate per minute

Occurs when an error was made during measurement

Installing Batteries

1. Place your thumb on the side of the battery cover. Pull the battery cover in the direction of the cuff to open.
2. Install lower battery first in the bottom position with the battery ribbon underneath. **Battery Type: 2 Alkaline (AAA) size.**
3. Place in the other battery and pull up the battery ribbon over it.
4. Replace the cover by placing the end tabs of the battery cover in the slots on the top of the battery compartment. Click in the other end to secure the battery cover.

Replace the batteries if:

1. The weak battery mark appears in the display.
2. Nothing appears in the display when the power is switched on.

As the supplied batteries are for test only, they may discharge earlier than batteries you buy in stores. Replace the batteries in pairs.

Remove batteries when unit is not in use for extended periods of time.

Note:  Discard used batteries as directed by your local, environmental codes.

HOW TO USE:

Helpful Tips

Here are a few tips to help you obtain more accurate readings:


- Do not measure your blood pressure immediately after consuming a large meal. To obtain more accurate readings, please wait one hour before measuring.
- Do not smoke or drink alcohol before measuring your blood pressure.
- You should not be physically tired or exhausted while taking a measurement.
- It is important that you relax during measurement. Try to take a 15-minute rest before taking a reading.
- Do not take measurements if you are under stress or under tension.
- Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
- If the monitor is stored at very low temperature (near freezing), have it placed at a warm location for at least one hour before using it.
- **Wait about 5 minutes before taking a SECOND blood pressure measurement.**

Applying the Cuff

1. Remove all watches, jewelry, etc., prior to attaching the wrist monitor. Clothing sleeves should be loosely rolled up and the cuff should be wrapped on bare skin for correct measurements.
2. Apply cuff to **LEFT** wrist with palm facing up as shown in Figure A.
3. Make sure the edge of the cuff is about 3/8 inch from the palm as shown in Figure B.
4. In order to insure accurate measurements, fasten the fabric fastener strap securely around your wrist so there is no extra space between the cuff and the wrist as shown in Figure C.
Note: If the cuff is not wrapped tightly enough, the measurement values may be false. However, do not apply the cuff too tightly as it may constrict the smooth flow of blood to the hand.
5. If your physician has diagnosed you with poor circulation on your left arm, carefully place the cuff around your right wrist as shown in Figure D.

Proper Posture During Measurement

1. Place your elbow on a table so that the cuff is at the same level as your heart as shown in Figure A.
Note: Your heart is located slightly below your armpit and a bit to the left of the middle of your chest.
Relax your entire body, especially the area between your elbow and fingers.

2. If the cuff is not at the same level as your heart or if you can not keep your arm completely still throughout the reading, use a soft object such as a folded towel to support your arm as shown in Figure B. Do not allow hard objects to come in contact with the wrist cuff.
3. Turn your palm upwards.
4. Sit upright in a chair, and take 5-6 deep, relaxing breaths.
5. Avoid leaning back while the measurement is being taken as shown in Figure C.




FIGURE C

Taking a Measurement


1. After applying the cuff on the wrist and you are properly positioned and relaxed, press the **ON/OFF/START** button “⏻”. All digits will light up, checking the display functions. The checking procedure will be completed after about 2 seconds.


2. After all symbols appear, the display will show a blinking “0”. At this time, the monitor is “**Ready to Measure**” and will automatically inflate the cuff to start measurement.


- Once the monitor inflates the cuff to approximately 180 mmHg, measurement will begin.

This monitor will re-inflate automatically to approximately 220 mmHg if the system detects that your body needs a higher pressure to measure your blood pressure.

- When the measurement is completed, systolic, diastolic, and pulse will be shown simultaneously and be saved automatically into the memory. **The last 90 measurement readings will be saved into memory.**
- Measurement is now complete. Press the ON/OFF/START button “⏻” to turn the monitor off. The unit will automatically shut-off after 1 minute if no buttons are pressed.


Note:

- This monitor automatically switches off approximately **1 minute** after the last button operation.
- To interrupt the measurement, simply press any key. The cuff will deflate immediately after a key is pressed.
- During the measurement, do not talk or move your arm or hand muscles.

Recalling Memory

1. To recall stored blood pressure readings, simply press the “M” button. The memory values will be shown on the LCD display. Every new press of the “M” button recalls one set of previously stored values. The most recent memory readings will be recalled first.

Note:

1. The data in the memory will not be lost even when the Blood Pressure Monitor is switched off.
2. You can erase the memories stored by one of 2 ways:
 - remove the batteries
 - press the memory recall key for more than 5 seconds

Troubleshooting

If a problem should arise during use, please check the following points.

Symptom	Check Point	Correction
No display when the POWER button is pressed	Have the batteries run down?	Replace them with 2 new “AAA” alkaline batteries
	Have the batteries been properly positioned in the compartment?	Re-insert the batteries in the correct position
EE mark shown on the display or the blood pressure value is displayed excessively low (or high)	Is the cuff placed correctly?	Wrap the cuff properly so that it is positioned correctly
	Did you talk or move during the measurement?	Measure it again quietly
	Did you shake the arm with the cuff on it?	

Note: If the unit still does not work, please see the Warranty Information in this manual. Under no circumstance should you disassemble the unit and attempt repairs yourself.

Precautionary Notes

1. Your Health o meter® monitor contains high-precision components. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or sharp impacts to the main unit, and protect it from dust.
2. Clean the blood pressure monitor body and the cuff carefully with a slightly dampened soft cloth. Do not press on the unit too hard. Do not wash the cuff or use chemical cleaners on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
3. Leaky batteries can damage the unit. Remove the batteries when the unit is not used for an extended period of time.
4. Do not press the ON/OFF/START key when the cuff is not placed around the wrist.
5. To avoid hazardous situations, the unit should not be operated by children.
6. Do not disassemble the main unit or cuff.
7. If the unit is stored near a freezing temperature, allow it to acclimate at room temperature before use.
8. Your Health o meter® Blood Pressure Monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to make any adjustments to it. If the monitor is not working properly please contact Sunbeam Products, Inc. as outlined under the Warranty information in this manual.
9. For users diagnosed with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, poor circulation of blood, kidney problems, or for users who have suffered from stroke, or for unconscious users, this device may not be suitable for use. Please consult your physician prior to use.

Specifications:

Power Source	: 3V DC Two (AAA) Alkaline Batteries
Measurement Method	: Oscillometric
Measurement Range	: Pressure: 20~280mmHg; Pulse: 40~180 beats/ minute
Accuracy	: Pressure: ± 3 mmHg; : Pulse: $\pm 5\%$ of reading
Pressure Sensor	: Semiconductor
Inflation	: Pump Driven Inflation
Deflation	: Automatic Pressure Release Valve
Memory capacity	: 90-set memory
Auto shut-off	: 1 minute after last button operation
Operation Environment	: Temperature: 10°C - 40°C (50°F-104°F); : Humidity: 40 - 85%
Storage Environment	: Temperature: -10°C - 60°C (14°F-140°F); 10 - 95% : RH max.
Weight	: 125g (without batteries)
Wrist Circumference	: 13.5 - 22 cm (5.3 inches - 8.7 inches)
Limited patient	: Over 18 years old
Dimensions	: 79 (L) X 72 (W) X 71 (H) mm
K.B.	: Short time operation is 2 minutes : Classification of type BF

**To ensure accurate measurements, it is recommended that calibration be completed every two years.*

**Specifications are subject to change without notice.*

Patents:

The Health o meter[®] Model 7633 is patented in countries named below:

Germany :	# 295-11-738.9
	# 297-21-574.4
Japan :	# 3020497
Taiwan :	# 140134
Mainland China :	# ZL 95217421.9
	# ZL 97249758.7

The Health o meter[®] Model 7633 is patent pending in other countries.

Warranty Information

Your Health o meter® Blood Pressure Monitor is covered under the limited Warranty offered by Sunbeam Products, Inc. ("Sunbeam"). Sunbeam warrants that for a period of five years from the date of purchase, this product will be free from defects in material and workmanship. Sunbeam, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. Sunbeam dealers, service centers or retail stores selling Sunbeam Products do not have the right to alter, modify or in any way change the terms and conditions of this warranty.

This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, disassembly, repair or alteration by anyone other than Sunbeam or an authorized service center. Further, the warranty does not cover Acts of God or Nature, such as fire, flood, hurricane and/or tornado.

Sunbeam shall not be liable for any incidental or consequential damages caused by the breach of any express or implied warranty. Except to the extent prohibited by applicable law, any implied warranty of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty. Some states, provinces or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or province to province.

For more information, or to obtain warranty service, contact Sunbeam Products, Inc. at 800 672-5625 or visit the Health o meter® website at www.healthometer.com

©2002 Sunbeam Products, Inc. All rights reserved.

Health o meter® is a registered trademark of Sunbeam Products, Inc.

and MicroStep™ is a trademark of Sunbeam Products, Inc.

Distributed by Sunbeam Products, Inc., Boca Raton, FL 33431.

Manual printed in China.

BLOOD PRESSURE LOG

Name: _____

My Target Blood Pressure is: _____

I am to call my healthcare practitioner:

if my blood pressure goes above _____ or falls below _____.

if I have the following symptoms: _____

<u>DATE</u>	<u>TIME</u>	<u>BLOOD PRESSURE</u>	<u>COMMENTS</u>
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____
_____	_____	/	_____

BLOOD PRESSURE LOG

Name: _____

My Target Blood Pressure is: _____

I am to call my healthcare practitioner:

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
Download from Www.Somanuals.com. All Manuals Search And Download.

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©2002 Sunbeam Products, Inc. Todos los derechos reservados
Health o meter® es una marca registrada de Sunbeam Products, Inc.
y MicroStep™ es una marca comercial de Sunbeam Products, Inc.
Distribuida por Sunbeam Products, Inc., Boca Raton, FL 33431.