
WEIGHT & BODY FAT MONITORING SCALE


User's Manual
Model BFM081-63

www.healthometer.com

WELCOME!

Congratulations on your purchase of this Health o meter® wellness monitoring scale. It's so much more than just a scale. It helps monitor your weight, body fat, body hydration, body mass index (BMI) and suggested daily caloric intake (DCI) for up to four different users. It's easy to use. Once your profile is set, your scale is ready to use. Also, you can save up to 5 records for each user.

IMPORTANT!
Please read all instructions to familiarize yourself with features and operations before programming the scale.

CAUTION!

- The use of this scale is not recommended for persons with an electric implant (heart pacemaker).
- The scale will not calculate body fat unless you have bare feet.
- For an accurate body fat estimate, do not bend your knees and keep your legs/thighs apart and straight when standing on the scale.
- The body fat function is not recommended for use by pregnant women, persons with fever, swollen legs or other edemas, as well as a person who is overhydrated or dehydrated.
- Body fat percentages may be slightly higher or lower for children under 16, persons with diabetes and other health conditions.
- Women naturally store approximately 5% more body fat than men.

SCALE FEATURES

The Health o meter® Body Fat Monitoring Scale can help you in achieving your physical best by providing five critical measurement estimates: BODY WEIGHT, BODY FAT COMPOSITION, BODY HYDRATION, BODY MASS INDEX, and SUGGESTED DAILY CALORIC INTAKE; in addition, it helps you track your progress to a goal by allowing you to set your weight goal and calculating the

weight differences between your starting point, your previously saved weight and your set weight goal. Traditional bathroom scales indicate your overall weight—both lean mass and body fat together. The Health o meter® Body Fat Monitoring Scale lets you track your weight, body fat and body hydration separately.

SELECTING UNIT OF MEASUREMENT

NOTE: If you choose to see your weight in Kilograms, then the scale will also change the measurement of height from feet and inches to centimeters.

- While the scale is OFF, hold down the "SELECT" button until the unit of measurement icon starts blinking.
- Press the ⬆️ or ⬆️ button to choose the desired unit of measurement (LB for pounds or KG for Kilograms).
- Press the "SELECT" button to confirm your selection.
- Press the ⬆️ or ⬆️ button to select if you want to reset the scale.
- Press the "SELECT" button to confirm your selection.

NOTE: If you choose "YES", then all information for all users will be erased.

FIRST TIME PROGRAMMING & SETTING YOUR GOAL/PROFILE

- Press the desired USER button and the HEIGHT will start blinking.
- Press the ⬆️ or ⬆️ button to choose the desired AGE. Then press the "SELECT" button.
- Press the ⬆️ or ⬆️ button to choose the desired GENDER: FEMALE ♀, MALE ♂, ATHLETE FEMALE ♀ or ATHLETE MALE ♂. Then press the "SELECT" button.

NOTE: Athlete mode is intended for individuals that are very fit and who dedicate 12 hours or more per week to intense physical activity such as swimming, biking or running. It does not include individuals that recently started making a commitment to intense physical exercise.

The display will flash the height, age, gender & "SAVED". Then, it will display the words "YES" and "GOAL?", with the word "YES" blinking.

4. To set a GOAL, press the ⬆️ or ⬆️ button to select "YES" then press the "SELECT" button.

NOTE: You can set a goal (by selecting "YES") or just track your weight (by selecting "NO").


- Press the ⬆️ or ⬆️ button to select your desired GOAL weight, and press the "SELECT" button to store it. The display will flash your goal weight and the word "SAVED", followed by "----". Then the scale will calibrate to "00" with either "LB" or "KG".

Your profile and goal (if desired) has been saved and the scale is ready for use.

USING YOUR BODY FAT MONITORING SCALE


A. MONITORING YOUR WEIGHT, BODY FAT AND BODY HYDRATION OVER TIME

- Press one of the "USER" buttons. Display will show the user profile, then it will calibrate ("----") and will show "00" with either "LB" or "KG" and the user number.

NOTE: Make sure to stand on the scale with bare feet.

- Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. Stand still while reading the display.
- The display will first show your weight and body fat percentage.
- After you step off of the scale, the display will scroll through your water level, BMI (Body Mass Index) and suggested DCI (Daily Caloric Intake).

NOTE: If the user setup a GOAL weight, the display will show the estimated weeks to achieve your goal with the suggested Daily Caloric Intake. If the user did not setup a goal, the fourth screen will show the amount of calories needed to maintain the current weight.

NOTE: In order to preserve battery life, your scale will automatically turn off after a short period of inactivity.

B. WEIGHING WITHOUT SELECTING A USER BUTTON

- Place the scale on a flat surface. A hard surface is best for maximum accuracy.
- Step onto the scale and position your feet so that your weight is distributed evenly over the scale platform. Stand still while the scale calculates your weight. The LCD display will flash "00" several times for a few seconds before displaying your weight. Stand still while reading the display.

NOTE: This scale offers both "Instant On" (Step onto scale) or "Tap On" (Tap on lower right corner and wait for "00" and "Step On"). During the "Instant On" weighing, if you see "----", the scale has recalibrated. Please step on the scale to weigh yourself again. Whether you remain standing on the scale, or you step off the scale immediately after your weight is displayed, your weight will continue to display for a few seconds and then turn off automatically.

REVIEWING SAVED RECORDS (HISTORY)


- Press the desired "USER" button. The user profile will be displayed.
- Press the ⬆️ or ⬆️ button to review up to five (5) saved records per user.

NOTE: M1 is the oldest record. Once there are 5 records saved, each new record will replace older records in memory.

RESETTING AN INDIVIDUAL USER'S HISTORY DATA


WARNING: This process will erase the records and the profile for the selected user.

- While the scale is OFF, press the "USER" button.
- Press and hold the "SELECT" button immediately to enter the user resetting function.
- Press the ⬆️ or ⬆️ button to choose "YES".
- Press the "SELECT" button to confirm the operation.
- The horizontal bars "----" will flash twice. The scale will prompt for the user setup process.

CHANGING A USER'S HEIGHT AND GOAL WEIGHT


- Press and hold the desired "USER" button for 3 seconds.
- Press the ⬆️ or ⬆️ button to select "YES". Then press the "SELECT" button to enter the goal setup.
- To change your goal weight, press the ⬆️ or ⬆️ button to reach your desired weight and press the "SELECT" button to store it.

COMPLETELY RESETTING THE SCALE

WARNING: This process will erase ALL history for ALL users including user profiles and saved records.

- While the scale is OFF, hold down the "SELECT" button until the unit of measurement icon starts blinking.
- Press the "SELECT" button to acknowledge the current unit of measurement (LB for pounds or KG for Kilograms).
- To confirm your desire to reset the whole scale, press the ⬆️ or ⬆️ button to choose "YES" and press "SELECT".
- The horizontal bars "----" will flash twice. All history and profile are erased.


REPLACING THE BATTERIES

This unit includes four (4) AA batteries. When "Lo" appears on the display, replace all batteries.

CAUTION: Do not mix old and new batteries. Do not mix alkaline, standard (Carbon-Zinc) or rechargeable (Nickel-Cadmium) batteries. Removing the battery will not affect the recorded data and user profile.

Battery Replacement

- Locate the battery cover on the back of the scale. Push the locking tab and pull the battery cover up.
- Use the tip of a pen to loosen each battery. Hold the scale at an angle to allow the batteries to slide out.
- Insert four new batteries, pressing them firmly into their slot.
- Replace the battery cover.

FCC STATEMENT

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on or off, the user is encouraged to try to correct the interference by one or more of the following measures: a) Reorient or relocate the receiving antenna; b) Increase the separation between the product and the receiver; c) Connect the product into an outlet on a circuit different from that to which the receiver is connected; d) Consult the dealer or an experienced radio/TV technician for help. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

TROUBLESHOOTING


If "E" appears on the upper row of the display:

- The scale cannot detect a stable weight. Stand still while weighing.
- The weight has exceeded the scale's capacity (400 LB/181 KG).

If "E" appears on the lower row of the display:

- The scale cannot measure body fat.
- Ensure you stand on the scale with bare feet.

NOTE: Display will turn off automatically in 20 seconds if no key is pressed.

IT IS NORMAL FOR YOUR WEIGHT TO VARY DURING THE DAY AND FROM ONE DAY TO ANOTHER.

CARING FOR YOUR SCALE

Periodically wipe the scale surface down with a DRY cloth to remove dust.

CAUTION: Weight exceeding the capacity (400 LB / 181 KG) may damage your scale. Ensure there is no weight being applied to the scale while not in use to prevent battery drainage. This product is designed for personal use only. Not for commercial or industrial use. This unit is not waterproof; avoid contact with excessive moisture.

BODY HYDRATION

Body hydration percentage indicates the amount of fluids you have in your body. Water in your system is important because it regulates body temperature, converts food into energy and protects vital organs. Individuals with more body fat have proportionately less total body water and are more susceptible to fluid imbalances that cause dehydration.

- When tracking your body hydration level you should consider the following:
- Take your hydration measurements the same time of day.
 - Monitor your readings to establish a personal body hydration range.
 - Do not drink any caffeinated or alcoholic drinks before taking a measurement.
 - Drink water if dehydration is detected.

Extracted from Dialysis & Transplantation, © 2001; Sports Nutrition For Young Adults: Hydration; Crystal Springs: Hydration Information; The National Women's Health Information Center: 4woman.gov. Please consult with your physician to determine the appropriate hydration level for your body.

BODY MASS INDEX (BMI)

Body Mass Index (BMI) is one way to tell whether you are at a healthy weight. It measures your weight in relation to your height, and is closely associated with measures of body fat. The higher the BMI, the greater the risk of developing health problems.

BMI for Adults

For adults over 20 years old, BMI falls into one of the following:

BMI Weight	Status
Below 18.5	Low
18.5 - 24.9	Normal
25.0 - 29.9	High
30.0 and Above	Very High


To determine the BMI that is appropriate for your body, consult your physician.

NOTE: The materials in this manual are based on information from the National Institutes of Health and Center for Disease Control (CDC). They are not intended to replace advice from your doctor or fitness professional. Please consult with your physician before beginning any fitness program or fat or weight reduction program. Individual weight loss will vary. Jarden Consumer Solutions takes no responsibility for individual results or any claim made by a third party.

DAILY CALORIC INTAKE (DCI) FEATURE

Your scale includes a daily caloric intake (DCI) feature. The DCI number is derived from the information you programmed on your user profile and your activity level.

If you have not set a goal weight, then the DCI number will represent the number of calories required to maintain your current weight.


If you have set a goal weight, the DCI number will represent the suggested calories (500 calories per day which represent 3,500 calories per week, equivalent to 1 LB) needed to increase or decrease your weight. The scale will also show the approximate length of time to achieve your goal.

NOTE: As you continue to lose/gain weight, the scale will adjust the amount of time required to reach your goal. The suggested DCI number will remain constant.

BODY FAT COMPOSITION RANGES

Extracted from Principals & Labs for Fitness and Wellness, 1st Edition ©1999. To determine the percentage of body fat that is appropriate for your body, consult your physician.

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