

# ***QuickMedical***

The quickest way to measure your health

## **Model-QM9800**



## **USER'S MANUAL**

## **INTRODUCTION**

The electronic Baby Scale QM9800 has been designed to help you monitor the growth of your baby, accurately and reliably quantifying how much milk has been assimilated at each feed.

Taking up to 20kg, and with 10g divisions, the scale will also help you monitor your child's progress in terms of weight gain during all stages of growth.

As soon as the child can stand up on his own, weighing operations can be performed without the pan, transforming the unit into practical bathroom scale for your infant.

The TARE function allows you to weigh the baby using a soft cloth or towel.

The scale comes with a useful WEIGHT-BLOCK function that allows a steady weight to be displayed despite the child wriggling about.

**YOU ARE STRONGLY RECOMMENDED TO READ THE INSTRUCTIONS AND WARNING HEREIN BEFORE USING THE EQUIPMENT AND TO KEEP THE BOOKLET IN A SAFE PLACE.**

## **POWER SUPPLY**

The Baby Scale uses a 9V alkaline battery only.

Remove the battery cover on the underside of the scale. Insert the battery and refit the cover.

When "LO" (low battery indication) is displayed, please replace a new 9V alkaline.

## **SPECIFICATIONS**

- Maximum Capacity: 20 kg / 44.40 lb
- Minimum Capacity: 10 g / 1/2 oz
- Weight Units: Kg / Lb
- Power: 9V alkaline battery
- Display: 4 active digits
- Material: plastic
- Function Keys: ON/OFF, UNIT, HOLD, TARE

## **KEY FUNCTIONS**

### **1. ON/OFF**

- 1) Turn on the scale by pressing this key.
- 2) Turn off the scale by pressing this key.

### **2. UNIT**

- 1) Press this key to select the unit of measure you require: Kg or pounds.
- 2) After turning the scale on, if the display shows 0.00, the weight is expressed in kg and g. If the display shows 0:0, the weight is expressed in pounds and ounces.

### **3. HOLD (BLOCK WEIGHT FUNCTION)**

- 1) Due to the baby moving, the weight shown on the display changes all the time. When the baby is laid on the plate, press the "HOLD" key.
- 2) The display shows "hold" and the scale automatically calculate the average of the weights it has taken.
- 3) The display then shows the baby's weight and at the top the symbol "□" starts flashing.

### **4. TARE**

If the baby is laid on a towel or cover, use the TARE function:

- 1) Turn the scale on with the ON/OFF key.
- 2) The display shows "8888".
- 3) Wait until the display shows "0.00" then place the towel or cover on the plate
- 4) Press the TARE key and wait until the display shows "0.00"; the scale is now ready for use.

## **USING THE SCALE WITH THE PLATE**

1. Remove the 4 green rubber plugs.
2. Fix the plate to the scale using the supplied locking nut. (See the illustration)
3. Press the ON/OFF key.
4. The display shows "8888".
5. Wait until the display shows "0.00". The scale is now ready for use.
6. Lay the baby on the plate and wait until the display shows its weight.
7. Press the ON/OFF key to turn the scale off, otherwise the scale turn off automatically after 2 minutes without being used.

## **USING THE SCALE WITHOUT THE PLATE**

1. When the baby can stand up on his own and does not weigh more than 20kg.
2. Remove the plate and keep the locking nut for future use.
3. Replace the 4 green rubber plugs.
4. Turn the scale on and wait until the display shows "0.00".
5. Get the baby to stand on the scale to start weighting.

## **ERROR MESSAGES**

1. "LO" indicates that the batteries are flat. Replace them.
2. "ERR" shows the weight is too much.
3. The scale is ready for use when the display shows "0.00". If the instructions have not been followed correctly, the display shows "ERR".

## **MAINTANCE**

- Cleans scale with a damp cloth.
- Being a precision instrument, the scale must be treated with due care.
- We guarantee accuracy and repeatability with a tolerance within 3%.

## **WARRANTY**

This appliance is warranted for a period of 2 years from the date of purchase as shown by the **dealer's stamp and signature and by the purchaser's receipt for tax purposes which must be kept attached to this warranty.**

The warranty only covers defects in materials and workmanship and is not applicable if damage results from an accident, incorrect use, negligence or improper use of the product. Only use the accessories supplied: using any other accessories may cause the warranty to be void.

**Do not open the appliance for any reason: if it is opened or tampered with, the warranty shall be void for good.**

Use of the appliance in public shops is not permitted.

The manufacturer disclaims any liability for injury to persons or pets or damage to property, whether direct or indirect, resulting from non-compliance with all provisions indicated in the relevant instruction booklet, especially concerning warnings regarding installation, use and maintenance of the appliance.

**QuickMedical**

P.O. Box 1052

8300 Railroad Ave. SE

Snoqualmie, WA 98065

[www.quickmedical.com](http://www.quickmedical.com)

[sales@quickmedical.com](mailto:sales@quickmedical.com)

Phone: 425-831-5963

Toll Free: 888-345-4858

Fax: 425-831-6032

## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>