**TANITA®**  
**IRONMAN™**

*InnerScan™*  
*Body Composition Monitor*

**Model : BC-548**

- GB** A guide to your body composition monitor
- D** Ein Handbuch für Ihre Körperanalysewaage
- F** Guide de votre impédancemètre
- NL** Een gids bij uw lichaamssamenstellingsmonitor
- I** Una guida al vostro monitor della composizione corporea
- E** Guía para el uso de su monitor de composición corporal
- P** Um guia para o seu monitor de composição corporal

Ironman and M-dot are registered trademarks of World Triathlon Corporation.

**Tanita UK Ltd.**  
The Barn, Philpots Close, Yiewsley,  
Middlesex, UB7 7RY UK  
Tel: +44 (0)1895 438577  
Fax: +44 (0)1895 438511  
www.tanita.co.uk

**Tanita Europe GmbH**  
Dresdener Strasse 25  
D-71065 Sindelfingen  
Germany  
Tel: +49 (0)-7031-6189-6  
Fax: +49 (0)-7031-6189-71  
www.tanita.de

**Tanita France S.A.**  
Villa Labrouste  
68 Boulevard Bourdon  
92200 Neuilly-Sur-Seine France  
Tel: +33 (0)-1-55-24-99-99  
Fax: +33(0)-1-55-24-98-68  
www.tanita.fr

**Tanita Corporation**  
14-2, 1-Chome, Maeno-Cho  
Itabashi-Ku, Tokyo, Japan 174-8630  
Tel: +81(0)-3-3968-2123  
Fax: +81(0)-3-3967-3766  
www.tanita.co.jp  
ISO 9001 Certified

**Tanita Corporation of America Inc.**  
2625 South Clearbrook Drive  
Arlington Heights, IL 60005 USA  
Tel: +1-847-640-9241  
Fax: +1-847-640-9261  
www.tanita.com

**Tanita Health Equipment H.K. Ltd.**  
Unit 301-303 3/F Wing On Plaza  
62 Mody Road, Tsimshatsui East  
Kowloon, Hong Kong  
Tel: +852 2838 7111  
Fax: +852 2838 8667

©2006 Tanita Corporation

BC5487901  
Download from [www.Softmanuals.com](http://www.Softmanuals.com). All Manuals Search And Download.

## WHO CAN USE A BODY COMPOSITION MONITOR

This Body Composition Monitor is intended for adults aged 18 – 99 years. Children aged 7 – 17 years can use the monitor for weight and body fat percentage readings only; the other features are not applicable to children.

Your monitor is also equipped with an athlete mode for adults with athletic body types. Tanita defines an athlete as a person involved in intense physical activity of approximately 10 hours per week and has a resting heart rate of 60 beats per minute or less. Individuals who have been an athlete for a number of years but currently exercises less than 10 hours per week can also use athlete mode.

Pregnant women should use the weight function only. All other functions are not intended for pregnant women.

The Tanita Body Composition Monitor is intended for home use only. This unit is not intended for professional use including hospitals, medical or fitness facilities; it is not equipped with the adequate specifications required for heavy usage. Using the unit in a professional location will invalidate your warranty.

The Tanita InnerScan Body Composition Monitor product provides readings for informational purposes only. This product is not intended to diagnose or treat any disease or abnormalities. Please consult with your physician with any questions or concerns related to your health.

## WHY MONITOR BODY COMPOSITION?

Body composition monitors are designed for healthier living by giving you an insight into key health indicators that will enable you to monitor the impact of changes to your lifestyle:

- See the impact of a change in diet on the body to make sure you are dieting the healthy way
- Fine-tune your fitness programme by monitoring progress of muscle mass and Basal Metabolic Rate
- Monitor visceral fat, a key indicator of the risk of Type 2 diabetes and heart disease.
- Set a target for your physique and monitor your progress towards it.

## HOW DOES A BODY COMPOSITION MONITOR WORK?

Tanita Body Composition Monitors calculate your body composition using Bio -electrical Impedance Analysis(BIA). Safe, low-level electrical signals are passed through the body via the patented Tanita footpads on the monitor platform. It is easy for the signal to flow through fluids in the muscle and other body tissues but meets resistance as it passes through body fat , as it contains little fluid. This resistance is called impedance. The impedance readings are then entered into medically researched mathematical formulas to calculate your body composition.

## WHEN IS THE BEST TIME TO USE MY BODY COMPOSITION MONITOR?

Your body water levels naturally fluctuate throughout the day and night. Any significant changes in body water may affect your body composition readings; for example, the body tends to be dehydrated after a long night sleep so if you take a reading first thing in the morning your weight will be lower and your body fat percentage higher. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may also cause variations in your hydration levels.


To get the most reliable reading it is important to use your Body Composition Monitor at a consistent time of day under consistent conditions. We suggest taking a reading before your evening meal.

## WHAT IS BODY FAT PERCENTAGE?

**Body fat percentage is the amount of body fat as a proportion of your body weight.**

Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer. The chart below shows the healthy ranges for body fat.

<sup>1</sup> Based on NIH/WHO BMI Guidelines.

<sup>2</sup> As reported by Gallagher, et al, at NY Obesity Research Center.

To determine the percentage of body fat that is appropriate for your body, consult your physician.

### Healthy Range Indicator (Standard Adults only\*)

Your Body Composition Monitor automatically compares your body fat percentage reading to the Healthy Body Fat Range chart. After your body fat percentage has been calculated, a black bar will flash along the bottom of the display, identifying where you fall within the Body Fat Ranges for your age and gender.

- (-) : Underfat; below the healthy body fat range. Increased risk for health problems.
- (0) : Healthy; within the healthy body fat percentage range for your age/gender.
- (+) : Overfat; above the healthy range. Increased risk for health problems.
- (++) : Obese; high above the healthy body fat range. Greatly increased risk of obesity-related health problems.

#### \*Note :

If you select Athlete mode, the unit will not display the Healthy Range Indicator. Athletes may have a lower body fat range depending on their particular sport or activity. (Range of user age is between 18-99)

## WHAT IS TOTAL BODY WATER PERCENTAGE?

**Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight.** Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a long night and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may cause variations in your hydration levels.

Your body water percentage reading should act as a guide and should not be used to specifically determine your absolute recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track the relative change.

Every individual varies but as a guide the average total body water percentage ranges for a healthy adult are:

**Female :** 45 to 60%

**Male :** 50 to 65%

**Source :** Based on Tanita's Internal Research

**Note:** The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body fat the total body water percentage should gradually move towards the typical range given above.

## WHAT IS BASAL METABOLIC RATE (BMR)?

### WHAT IS BMR?

Your Basal Metabolic Rate (BMR) is the minimum level of energy your body needs when at rest to function effectively including your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You burn calories when sleeping.

About 70% of calories consumed every day are used for your basal metabolism. In addition, energy is used when doing any kind of activity however; the more vigorous the activity is the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the quantity of muscles you have, therefore increasing your muscle mass will help increase your basal metabolism.

By studying healthy individuals, scientists have found that as people age, their metabolic rate changes. Basal metabolism rises as a child matures. After a peak at the age of 16 or 17, it typically starts to decrease gradually.

Having a higher basal metabolism will increase the number of calories used and help to decrease the amount of body fat. A low basal metabolic rate will make it harder to lose body fat and overall weight.

### HOW DOES A TANITA AN IRONMAN BODY COMPOSITION MONITOR CALCULATE BMR?

The basic way of calculating Basal Metabolic Rate BMR is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).\*

\*Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week : A Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition

## WHAT IS MUSCLE MASS?

This feature indicates the weight of muscle in your body. The muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increases, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.

## WHAT IS PHYSIQUE RATING?

This feature assesses your physique according to the ratio of body fat and muscle mass in your body.

As you become more active and reduce the amount of body fat, your physique rating will also change accordingly. Even though your weight may not change, your muscle mass and body fat levels may be changing making you healthier and at lower risk of certain diseases.

Each person should set their own goal of which physique they would like and follow a diet and fitness programme to meet that goal.

Result	Physique Rating	Explanation
1	Hidden obese	<b>Small Frame Obese</b>
		This person seems to have a healthy body type based on physical appearance; however, they have a high body fat % with low muscle mass level.
2	Obese	<b>Medium Frame Obese</b>
		This person has a high body fat percentage, with a moderate muscle mass level.
3	Solidly-built	<b>Large Frame Obese</b>
		This person has both a high body fat % and a high muscle mass.
4	Under exercised	<b>Low Muscle &amp; Average Body Fat%</b>
		This person has an average body fat % and a less than average muscle mass level
5	Standard	<b>Ave. Muscle &amp; Ave. Body Fat %</b>
		This person has average levels of both body fat and muscle mass.
6	Standard Muscular	<b>High Muscle &amp; Ave. Body Fat % (Athlete)</b>
		This person has an average body fat % and higher than normal muscle mass level.
7	Thin	<b>Low Muscle &amp; Low Fat</b>
		This person has both a lower than normal body fat % and muscle mass level.
8	Thin and muscular	<b>Thin and muscular (Athlete)</b>
		This person has lower than normal body fat % while having adequate muscle mass.
9	Very Muscular	<b>Very Muscular (Athlete)</b>
		This person has lower than normal body fat % while having above average muscle mass.

Source : Data from Columbia University (New York) & Tanita Institute (Tokyo)

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