waterpik

Ultra DENTAL WATER JET





Prepare the Reservoir:



Remove the reservoir from base and fill with lukewarm water.



Replace the reservoir onto the base, **pressing down firmly** to ensure full seating.



Insert Tip:

Select the appropriate tip and insert into the handle. Press firmly until the tip clicks into place.

<u>Caution</u>: Do not press tip eject button while unit is running.





Adjust Water Pressure:

For first time use, adjust the pressure dial on the base to **lowest position**.

Gradually increase pressure as you prefer.





Pause Button Control:



You can temporarily stop water flow by pressing the pause button on the handle.

continued on back

Download from Www.Somanuals.com. All Manuals Search And Download



On/Off Procedure:



Lean over sink and close lips enough to prevent splashing, while still allowing water to flow from mouth into the sink, Turn unit ON.









Always **PAUSE** or turn unit **OFF** before removing tip from mouth.



Technique Tips for Best Results:

Aim the tip just above the gum line at a 90 degree angle.



Starting with the back teeth, follow the gum line and pause briefly between teeth.

Continue until all areas around and between teeth have been cleaned.



When finished:

For more details refer to the instruction booklet.



Free Manuals Download Website

http://myh66.com

http://usermanuals.us

http://www.somanuals.com

http://www.4manuals.cc

http://www.manual-lib.com

http://www.404manual.com

http://www.luxmanual.com

http://aubethermostatmanual.com

Golf course search by state

http://golfingnear.com

Email search by domain

http://emailbydomain.com

Auto manuals search

http://auto.somanuals.com

TV manuals search

http://tv.somanuals.com