BRAUN

Oral-B

Braun Oral-B Plak Control MD 9000 Oral Irrigator

Type 4714

Munddusche mit stufenloser Aufprallkraft- und Impulseinstellung
Oral irrigator with continuous water impact and impulse
Hydropulseur à pression d'eau et impulsions continues
Idropulsore con selezione continua del getto d'acqua
Monddouche met constante waterdruk en constante waterpulsen
Εκτόξευσης νερού με συνεχή πίεση και ψεκασμό νερού

Deutsch

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English

Our products are engineered to meet the highest standards of quality, functionality and design. We hope you thoroughly enjoy your new Braun Oral-B appliance.

The Braun Oral-B Plak Control oral irrigator removes remaining food particles from between your teeth, while it also effectively massages your gums. For this reason, the oral irrigator should always be used after brushing your teeth.

Prior to using any Braun Oral-B appliance we recommend that you carefully read the following use instructions so that you and your family know how to best operate the unit safely.

Caution

Check the cord occasionally for any damage. Should this be the case, unplug the motor part ② immediately and take the appliance to a Braun Service Centre. A damaged or non-functioning unit may no longer be used.

Description

- ① Water container
- ② Motor part
- ③ Removable compartment for attachable jets ④
- ④ Four attachable jets with different colour markings
- ⑤ Water impact control
- ⑥ Press button for maximum water impact in the on/off switch ⑦
- ⑦ On-/off switch
- ⑧ Irrigator handle
- ⑨ Cord storage
- ⑩ Connections for the charging unit of a plaque remover, Braun Oral-B Plak Control
- ⑪ Wall holder
- ⑫ Two screws with wall plugs

Wall holder

Your oral irrigator comes with a wall holder ⑪ for the motor part ②. Use the enclosed screws and wall plugs ⑫ to mount (a). The two grooves ⑬ on the back of the motor part ② are pushed down onto the hooks ⑭ of the wall holder ⑪. The cord storage compartment ⑨ offers practical storage of the cord. Only use as much cord as you need to reach the socket.

Specifications

Voltage: 220–240 Volt (~) ac, 50–60 Hz
Power intake: 20 Watt (therefore also suitable for shaver sockets)
Max. operating time: 20 minutes
Cooling off time: 2 hours
Water capacity: approx. 600 ml
Number of impulses: approx. 500–1000 impulses per minute, continuously adjustable

Switch

Turn on the oral irrigator by pushing the on/off switch ⑦ up. Turn it off by sliding the on/off switch back to its initial position (i.e. below the arrow).

The water impact control ⑤ regulates the impulse pressure: gentle = 2, strong = 8. The water impact is continuously variable and electronically controlled. You may press the button ⑥ in the on/off switch ⑦ to obtain maximum water impact, independent of the impact shown (2–8). By this, even the most resistant food particles can be removed from between the teeth.

Getting started

After connecting the unit to a socket, remove water container ① and fill with lukewarm water. Then place the water container firmly onto the unit (c). Before using for the first time, allow the water to run through completely. Take the handle

⑧ out of the magnetic holder (b), hold it over the washbasin and turn the unit on.

Use

Fill the water container ① with water again and add mouthwash if desired. Place it firmly onto the motor part (c). Attach one of the jets ④ to the handle ⑧ until it clicks in place (d). Turn the jet to the direction desired.

Select the water pressure with the water impact control ⑤; it is best to start with the gentlest pressure (setting 2), especially if your gums tend to bleed or if children are using the oral irrigator. To remove resistant food particles, you may use the press button ⑥ for maximum water pressure (corresponding to setting 8 on the water impact control ⑤).

Bend over the washbasin and guide the jet to your teeth. Then turn the unit on with switch ⑦. The water jet impulses will start in a few seconds. Keep your mouth slightly open to allow the water to run into the washbasin.

Single and multiple jet

The oral irrigator jet removes remaining food particles from between your teeth while also massaging your gums. There are two different jet functions:

Single jet: Point of the jet in position (e) to clean the teeth and the spaces between the teeth and to remove resistant food particles (materia alba).

Multiple jet: Point of the jet in position (f) to massage the gums and keep them healthy.

Begin with the single jet. Guide the water to the teeth and especially to the spaces between the teeth. Then massage the gums with the multiple jet.

The jet can easily be switched from one function to the other. Turn off the oral irrigator with switch ⑦ (to avoid uncontrolled water spraying) and press the point of the jet in or out. Then turn the unit on again with switch ⑦.

After use

Do not turn off the unit unless the water container is entirely empty. Bacteria could otherwise collect in the remaining water if the unit is not used for some time.

After use, return the handle ⑧ to the magnetic holder (b). The attachable jets ④ can be stored in the removable compartment ③.

Cleaning

Clean the motor part ② with a damp cloth. Wipe the handle ⑧ dry. The water container ① and the removable compartment for the jets can be washed in the dishwasher.

Dental centre

Your Braun Oral-B Plak Control oral irrigator can easily be turned into a complete dental centre: The connections ⑩ on the side of the motor part ② allow you to connect a Braun Oral-B Plak Control plaque remover. The plaque remover and the oral irrigator are then operated as a dental centre using the cord of the irrigator.

The charging unit for the plaque remover (type 4723) is available at your dealer.

Subject to change without notice.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC).

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English

Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either

by repairing or replacing the complete appliance as we may choose.

This guarantee extends to every country where this appliance is supplied by Braun or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Braun Customer Service Centre.

For UK only:

This guarantee in no way affects your rights under statutory law.

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Australia

Gillette Australia Pty. Ltd.,
Private Bag 10, Scoresby,
Melbourne, Victoria 3179,
☎ 1 800 641 820

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Hong Kong

Audio Supplies Company,
Room 506,
St. George's Building,
2 ICE House Street,
Hong Kong,
☎ 25 24-93 77

Hungary

Kisgep KFT,
Pannonia u. 23,
1136 Budapest,
☎ 13 494 955

Iceland

Verzlunin Pfaff h.f.,
Grensasvegur 13,
Box 714, 121 Reykjavik,
☎ 5 33 22 22

India

Braun Division,
c/o Gillette Div. Op. Pvt. Ltd.,
34, Okhla Industrial Estate,
New Delhi 110 020,
☎ 11 68 30 218

Iran

Tehran Bouran Company
Irtuc Building,
No 874 Enghelab Ave.,
P.O. Box 15815-1391,
Tehran 11318,
☎ 021 67 03 50

Ireland (Republic of)

Gavins Shaver Centre,
83/84, Lower Camden Street,
Dublin 2,
☎ 1800 509 448

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Philippines

Gillette Philippines Inc.,
Braun Service Centre
Villongco Road, Sucat,
Paranaque City,
☎ 2 8 42 31 54

Poland

Gillette Poland S.A.,
Budynek Orion,
ul. Domaniewska 41,
02-672 Warszawa,
☎ 22 548 89 74

Portugal

Gillette Portuguesa, Ltd.,
Braun Service,
Rua Tomás da Fonseca,
Torre G-9ºB,
1600-209 Lisboa,
☎ 808 2 000 33

Qatar

Al Baker Trading Co.,
P.O. Box 9, Doha,
☎ 04 15 0 57

Réunion

Dindar Confort,
58, rue Maréchal Leclerc,
P.O. Box 12,
97400 St. Denis,
☎ 40 12 02

Romania

Gillette Romania srl.
Calea Floreasca nr. 133-137
et 1, sect 1,
71401 Bucuresti
☎ 01-2319656

Russia

RTC Sovinservice,
Rusakovskaya 7,
107140 Moscow,
☎ (095) 264 96 02

Saudi Arabia

Salem M. Bakhshwain Sons,
P.O. Box 743,
Tabuk Street,
21421 Jeddah,
☎ 657 31 11

Schweiz/Suisse/Svizzera

Telion AG, Rütistrasse 26,
8952 Schlieren,
☎ 0844-88 40 10

Singapore (Republic of)

Beste (S) Pte. Ltd.,
6 Tagore Drive,
03-04 Tagore Industrial
Building,
Singapore 787623,
☎ 4 552 24 22

Slovakia

Techno Servis Bratislava spol.
s.r.o.,
Bajzova 11/A,
82108 Bratislava,
☎ (02) 555 68 161

Slovenia

Coming d.d.,
Slovenceva 024,
1000 Ljubljana,
☎ 386-01-560-40-00

South Africa (Republic of)

Fixnet After Sales Service,
17B Allandale Park,
P.O. Box 5716,
Cnr Le Roux and Morkels
Close,
Johannesburg 1685,
Midrand, Gauteng,
☎ 11 315 9260/1

St. Maarten

Ashoka,
P.O. Box 79,
Philipsburg,
Netherlands Antilles,
☎ 52 29 31

St. Thomas

Boolchand's Ltd.,
31 Main Street,
P.O. Box 5667,
00803 St. Thomas,
US Virgin Islands,
☎ 340 776 0302

Suomi

Gillette Group Finland Oy/
Braun,
P.O. Box 9,
Niittykatu 8, PL 9,
02200 Espoo,
☎ 09-45 28 71

Sverige

Gillette Group Sverige AB,
Dept. Sweden, Stockholm
Gillette
Räundavägen 12,
Box 702,
16927 Solna,
☎ 020-21 33 21

Syria

c/o Khalil Daoud
Amro Ibn Kalthoum Str.
Unisyria, P.O. Box 35002,
Damascus,
☎ 11 622 28 81

Taiwan, R.O.C.

Audio & Electr. Supplies Ltd.,
Brothers Bldg., 10th Floor,
85 Chung Shan N Rd., Sec. 1,
Taipei (104),
☎ 0 80 221 630

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Yemen (Republic of)

Saba Stores for Trading,
26th September Street,
P.O. Box 5278,
Taiz,
☎ 967 4-25 23 80

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>

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Driva S.A.,
Gillette de Venezuela S.A.,