Digital Thermometer

- **WARNING:** This product is not a toy. Store away from baby and young children's reach.
- **DO NOT** immerse any part of the thermometer in water.
- **DO NOT** drop the thermometer! Protect it from severe impact and shock.
- **DO NOT** place thermometer near water and moisture.
- Keep thermometer away from heat sources, such as fireplaces, radiators, stoves, and TVs or other appliances. Heat can damage the case or electrical parts and effect the temperature reading.
- **DO NOT** use the thermometer unit if:
 - The unit has been exposed to liquids.
 - Unit has been dropped or is damaged.
- Bending, biting, or chewing the thermometer will cause damage.
- Keep thermometer at room temperature for at least 15 minutes before use.
- **DO NOT OPEN THE THERMOMETER** as there are no user-serviceable parts inside.
- When the "▼" symbol (upside-down triangle) appears on the right side of the display, the battery is exhausted and the thermometer needs to be replaced. When discarding this product, be sure to dispose of it according to your state and local regulations.
- This thermometer is for oral use only and should not be used in any body location other than the mouth.

WARNING

Modifications not authorized by the manufacturer may void user's authority to operate this device.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:


- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Important Safety Instructions

- The instrument may be used only for measuring body temperature!
- Protect the instrument from impact and dropping!
- Avoid ambient temperatures above 60 °C/140°F. NEVER boil the thermometer!

Note:

- If water accidentally splashes on the device, please simply dry the device with any clean cloth before taking measurement. It won't affect the accuracy of the measurement.
- This product is for home use.

Turning on the Thermometer

To turn on the thermometer, briefly press the ON/OFF button.

1. A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C/90°F, an «L» and a flashing «°C» appear at the display field.
2. The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button. Otherwise the thermometer will automatically turn off after about 10 minutes.

Change the scale:

If the ON/OFF button is held down for more than 3 seconds when turning on the thermometer, then LCD will display the F/C, so you could press again to switch the scale. About 2 seconds after the button is released, the temperature value disappears and the thermometer is ready to take a measurement.

Measuring Method

Try to keep your mouth closed for 5 minutes before taking the temperature. Place the tip under the tongue at the back of the mouth. For the most accurate read, place the tip as far back under the tongue as possible. Press the tongue down over the tip, and close your mouth. Hold the thermometer steady keeping the tip in the same position for 1 minute.

- ✓ = correct areas
- X = incorrect areas

Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Do not submerge in any liquid.


Technical Specifications:

Type:	Maximum temperature thermometer
Temperature measurement range:	90 °F to 110 °F If measured temperature is below 90 °F, an "L" will display. If measured temperature is above 110 °F, an "H" will display.
Temperature measurement accuracy:	± 0.2 °F
Battery:	Non-replaceable. 1.5V. Type SR41.
Note:	This thermometer complies with ASTM E 1112, EN 12470-3, EC 60601-1, IEC 60601-1-2

Medicine spoon

1. Hold spoon and pour medicine to prescribed dosage according to the dosage directions on the bottle.
Note: Insert stopper into medicine spoon to prevent spills.
2. Hold baby upright when administering medicine. Unplug stopper then administer medicine slowly enough so your child may swallow comfortably.
3. Pause and repeat if necessary.



Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Wash thoroughly before first use and after each use. Do not clean in dishwasher.

Medicine syringe

1. Push plunger down and insert tip of syringe into medicine.
2. Slowly pull plunger up to fill with medicine to the prescribed dosage according to the dosage direction on the bottle. Line the level of the liquid medicine up to the prescribed dosage line on the syringe body.
3. Insert tip of syringe in child's mouth pointing towards cheek.
4. Push plunger in slowly to dispense medicine. Press slow enough to allow child to swallow comfortably.
5. Pause dispensing medicine to allow child swallow if necessary.

Care and Maintenance: Pull plunger out of syringe to clean. Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Repeat after each use. Do not clean in dishwasher.


Nasal aspirator

1. Secure baby in an upright position.
2. Squeeze and hold aspirator bulb and gently insert into baby's nostril pointed toward the side. **DO NOT** point tip toward center of baby's nostril. **DO NOT** insert tip too far into baby's nostril.
3. Release bulb slowly to remove mucus from baby's nostril.
4. Repeat on other nostril.
5. Repeat as necessary until congestion is relieved.

• Never push the aspirator too far into the nostril and always dry thoroughly after each use. Never use in baby's mouth or ears.


Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Squeeze and release the bulb several times to clean the inside of the bulb. Rinse with warm water. Squeeze out all excess water from the bulb. Do not clean in dishwasher.

Brush and Comb

Gently brush baby's hair with the soft bristles of the brush without worry for baby's sensitive scalp. Use either of the two teeth options to comfortably comb hair despite wetness or thickness.

Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Allow to drip dry. Do not clean in dishwasher.

Baby scissors

Our baby scissors make cutting your baby's delicate nails safe. The round tips ensure extra safety when handling your baby's nails.

NOTE: Cut fingernails following their natural shape and toenails straight across. Always leave a little white showing at the end of the nail.

Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Wash thoroughly before first use and after each use. Do not clean in dishwasher.

Nail clippers and Emery boards


WARNING: The nail clipper contains batteries and is not a toy. Keep product and batteries out of reach of children. Properly dispose of used batteries.

The curved clipping edge on our clippers makes clipping your baby's delicate nails safe and easy.

1. Lift the clip arm from hand grip and rotate 180°.
2. Rotate light source up until you hear a click and light comes on.
3. Gently push skin away from nail. Clip fingernails following their natural shape and toenails straight across. Always leave a little white showing at the end of the nails.
4. Use emery board to smooth baby's nails as needed.

Battery installation: Requires 2 -CR1616 3-volt lithium batteries (included). Replace when light begins to dim. Lift up clip arm and unfasten battery cover. Insert 2 new batteries, one on top of the other, positive (+) sides up, under the battery terminal. Resecure the battery cover.

Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Do not submerge in any liquid. Remove batteries during long periods of non-use. Do not clean in dishwasher.

Tweezers


Ideal for removing small, hard to grip items.

Care and Maintenance: Disinfect tweezers with rubbing alcohol.

Toddler toothbrush

The toothbrush's contoured handle is ideal for small hands to grasp. Clean toddler's teeth and gums comfortably with the toothbrush's soft bristles.

Care and Maintenance: Hand wash in warm water before and after use. Do not clean in dishwasher.

Fingerbrush

Can be used to massage sore teething gums.

1. Slide over parent's fingertip.
2. Lightly rub gums and teeth.

Note:

- Large bumps for massaging gums.
- Smaller bumps for cleaning out teeth.

Care and Maintenance: Hand wash in warm water with soap or non-citrus based detergents and rinse thoroughly. Wash thoroughly before first use and after each use. Do not clean in dishwasher.

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