

Invacare<sup>®</sup> Therapeutic Support Surfaces

# Invacare<sup>®</sup> Therapeutic Support Surfaces Group | Products

he Invacare<sup>®</sup> Group I Therapeutic Support Surfaces are designed to aid in the prevention and treatment of pressure ulcers. Three types of modalities allow you to choose the right product for the individual user. Product types include an alternating pressure pad system, gel overlay and a foam mattress replacement. All of these products are simple to use and easy to maintain.



## Invacare<sup>®</sup> CareGuard<sup>™</sup> 101 Home Care Foam Mattress

Model no. CG10180

- · Torso section of top layer anatomically sculptured for body conformity
- Head and heel section of top layer constructed of softer, pressure absorbing foam
- Multiple layers of foam combine to provide the maximum pressure reduction and durability - 1.85 lb., 29 ILD head and heel top layer, 1.8 lb., 36 ILD torso top layer and 46 ILD bottom layer
- · Reversible head-to-foot with contour cuts to reduce shear
- Articulation cuts on bottom layer facilitate gatching
- · Anti-microbial waterproof cover with a 70 denier top to reduce friction
- Non-slip bottom cover provides secure positioning on bed sleep surface
- · Foam meets California Technical Bulletin #117 for fire retardancy and mattress cover meets 16 CFR Part 1632, California Technical Bulletin #106
- Designed to meet EO184
- Does not contain latex

# Invacare<sup>®</sup> CareGuard<sup>™</sup> Gel Foam Mattress Overlay

Model no. CGGFMO

- Four large 2" thick gel bladders inside polyurethane foam
- · Reversible for increased service life
- Durable aggressive base material with strong elastic corner securing straps reduce overlay movement
- · Vapor permeable waterproof cover with low friction and shear (Do not machine wash)
- Meets California Technical Bulletin #117
- Designed to meet EO185
- Does not contain latex




### Invacare<sup>®</sup> CareGuard<sup>™</sup> APP Alternating Pressure Pad System

Model no CG9701

- Advanced pump technology alternately inflates and deflates air cells
- Variable pressure setting for added patient comfort
- · Pump can be mounted on the bed with built-in bracket
- Pad meets California Technical Bulletin #116
- Designed to meet EO180

# **Specifications**

Model no.	Weight	Length	Width	Thickness	Power	Code Approval	Weight Capacity	Limited Warranty*
CG10180	19 lb.	80"	35"	5"	n/a	n/a	250 lb.	2 years
CGGFMO	55 lb.	76"	34"	3.5"	n/a	n/a	250 lb.	6 months
	Pump: 3.5 lb.	10"	6"	4"	115V, 160 HZ	UL-1431		2 years
CG9701	Pad: 4.6 lb.	118"	34"	2.5"	n/a	n/a	250 lb.	30 days



- Does not contain latex



Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

### Canada

570 Matheson Blvd. E Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare is a registered trademark and Yes, you can., the blue medallion design and CareGuard are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2003 Invacare Corporation

Form No. 00-307 Rev. 12/03



"Invacare says 'Yes, you can™,' and I couldn't have said it better myself."

Download from Www.Somanuals.com. All Manuals Search And Download.

Free Manuals Download Website <u>http://myh66.com</u> <u>http://usermanuals.us</u> <u>http://www.somanuals.com</u> <u>http://www.4manuals.cc</u> <u>http://www.4manuals.cc</u> <u>http://www.4manuals.cc</u> <u>http://www.4manuals.com</u> <u>http://www.404manual.com</u> <u>http://www.luxmanual.com</u> <u>http://aubethermostatmanual.com</u> Golf course search by state

http://golfingnear.com Email search by domain

http://emailbydomain.com Auto manuals search

http://auto.somanuals.com TV manuals search

http://tv.somanuals.com