he Invacare® MNS400 is a full-depth (8") powered alternating pressure air mattress replacement system which reduces interface pressure for the treatment of Stage I to Stage IV pressure ulcers. The MNS400 reduces interface pressures below capillary closure by alternating pressure in the 8-inch air cell bladders every 5, 10 or 15 minutes. The power unit control combines simplified operations with effective results.



FEATURES

- Portable, lightweight mattress replacement system fits most standard medical beds
- Single pump 50 liter per minute air flow
- 16 Individual 8" nylon air cells reduce interface pressure below capillary closure while providing low interface pressure over the entire alternating cycle
- Mattress cells can be set to alternate every 5, 10 or 15 minutes to reduce interface pressures

- Waterproof, anti-microbial Ventex Recovery 5 therapy cover reduces friction and shear
- Durable ballistic nylon bottom protects the mattress cell from any abrasive surface
- Lower safety mattress supports the patient in the event of power failure
- · Durable lower transport mat for patient transport or power outages



User-Friendly Controls

- 10 pressure settings for patient comfort range
- Power failure alarm
- · Power unit is lightweight and compact
- · Soft touch, easy wipe, sealed keyboard
- · Low pressure alarm
- · Dual connector rapid attach and CPR deflate at the power unit



Model no. MNS400-S (Standard width) Model no. MNS400-N (Narrow width) Model no. MNS400-B (Bariatric width)

Electrical power	115vac/60Hz
Dimensions	
Power unit	11" W x 12.5" H x 5.25" D
Standard Mattress	80" L x 36" W x 8.5" H
Narrow Mattress	80" L x 32" W x 8.5" H
Bariatric Mattress	80" L x 42" W x 8.5" H
\ A /- ! - I - I	

Weight M

V

Power unit

Mattress	
Standard mattress	22 lb.
Narrow mattress	20 lb.
Bariatric mattress	24 lb.

Veight capacity	
MNS400-S	350 lb.
MNS400-N	350 lb.
MNS400-B	600 lb.

viateriai	
Mattress cells	70 denier nylon with a
	polyurethane coating
	CA# 117 Fire Retardency
Therapy pad	70 denier nylon top surface

with a polyester backing for comfort CA# 117 Fire Retardency

Safety code approval CE, UL2601, CSA, ETL Limited warranty

Power unit 1 year Mattress 6 months



Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada 570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare and the blue medallion design are registered trademarks and Yes, you can. is a trademark of Invacare Corporation. Specifications are subject to change without notification.

Manufactured for Invacare Corporation by Genadyne Biotechnologies, Inc.

© 2005 **Invacare Corporation**

Form No. 04-216



"Invacare says 'Yes, you can.™' and I couldn't have said it better myself.



Alternating Mode

In alternating mode, the 16 adjacent air cells of the Invacare MNS400 alternately inflate and deflate at 5, 10, or 15 minute adjustable intervals, which redistributes the pressure against the skin to promote capillary circulation.





Free Manuals Download Website

http://myh66.com

http://usermanuals.us

http://www.somanuals.com

http://www.4manuals.cc

http://www.manual-lib.com

http://www.404manual.com

http://www.luxmanual.com

http://aubethermostatmanual.com

Golf course search by state

http://golfingnear.com

Email search by domain

http://emailbydomain.com

Auto manuals search

http://auto.somanuals.com

TV manuals search

http://tv.somanuals.com