

Invacare®
Nutron™ R50LX
Power
Wheelchair

Invacare® Nutron™ R50LX Power Wheelchair



Yes, you can.™



Yes, you can.™

Is it possible for one power chair to be both stylish and powerful? The Nutron™ R50LX is both. Its compact design and dependable rear-wheel drive make it easy to maneuver around obstacles inside and outside the home. And when the consumer wants to travel, the R50LX disassembles and folds for easy transport.



The R50LX can be disassembled and folded for easy transport.

Specifications

Overall dimensions	
Length*	40.75"
Widths (excluding joystick)	20" - 26"
Heights	
17.75" seat-to-floor	35" - 38"
19.75" seat-to-floor	37" - 40"
Seat dimensions	
Widths	14" - 20"
Depths	16" - 18"
Back heights	16" - 19"
Adjustable-angle back	-5° to +20° (5° increments)
Seat angle	3° fixed
Weight capacity	250 lb.
Wheel and caster sizes	
Front casters (standard)	8" x 2"
Front casters (optional)	6" x 2"
Drive wheels	12.5" x 2.25"
Maximum speed	4 mph
Incline capability	9°
Turning radius**	
Front (footrest in place)	33"
Rear (anti-tippers in place)	12"
Weight	
With 22NF batteries	155 lb.
Without batteries and front riggings	85 lb.
Ground clearance	3.5"
Warranties	
Base frame	60 months
Seat frame	36 months
Electronics	12 months
Motor/gearbox	12 months

* Measured with PW93 front riggings and anti-tippers with 16" seat depth.

** Measured with 16" W x 16" D seat.

Features

- Rear-wheel drive and compact size provide optimum maneuverability and accessibility
- Oversize, high-strength tubular frame with bold styling
- Seat-to-floor heights as low as 17.75" permit easier transfers
- Invacare MKIV™-R11 LP electronics for precise driver control
- Reinforced crossbrace design folds to 13" for transport
- Back angle and height adjust for maximum comfort and support
- Telescoping front riggings allow optimal foot placement
- Flat-free wheel inserts mean worry-free operation
- Powerful 350-watt, high-torque motors ensure dependable performance
- Regenerative braking and automatic motor locks for smooth stopping
- Nylon upholstery adds style and comfort
- 8-amp charger provides maximum power
- Waterfall arm pads offer user support
- Optional fold-down back permits maximum portability
- 250-lb. weight capacity
- Chair is available in 5 frame colors

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare is a registered trademark and Yes, you can., Nutron and MKIV are trademarks of Invacare Corporation. Specifications are subject to change without notification.

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>