



eed a lightweight folding wheelchair with simple adjustability and reliable performance at an economical price? The Invacare® Patriot™ is the model of choice. The classic manual rehab wheelchair is a clean design offering valuable features including: a lightweight, aluminum frame; adjustable seat to floor heights; durable center mount front riggings; and adjustable angle head tubes all in a compact folding design.



Simple, triangular axle bracket offers 8 different axle positions – 4 axle positions in font of back canes and 4 axle positions behind back canes.

## **Highlights**

- · Lightweight, aluminum frame
- Adjustable seat-to-floor height and rear axle positions
- · Adjustable angle head tubes
- Durable center mount front riggings
- Extra smooth precision head tube bearings
- Compact for transport (11.75" wide folded with wheels)
- Adjustable-height back available in 1" increments, 15"-19"
- Multiple arm, footplate and tire options
- 14 color options
- Quick-release axle option
- · Seat pouch standard



Precision head tube bearings allow casters to turn smoothly.

# **Specifications**

- peemeasiem				
Overall				
Width open	7.50" plus seat width*			
Width closed	11.75"*			
Depth (w/riggings)	43" or 45"			
Seat widths	14", 16", 17", 18", 20"			
Seat depths	16", 18"			
Front Seat-to-Floor		Front Caster	Rear Wheel	

it seat-to-i looi			Tront Caster	iteai vviicei		
			Mounting	Mounting Po	osition	
Rear Wheel Size	Seat-to-Floor	Caster	Position	Axle Mount	Axle Mounting Bracket	
20"	15.50"	6"	top	middle	upper	
	16.50"	6"	t	refer to ser	refer to service manual	
	17.50"	6"	bottom	bottom	lower	
	17.50"	8"	top	bottom	lower	
	18.50"	8"	t	refer to service manual		
22"	15.50"	6"	†	refer to ser	vice manual	
	16.50"	6"	†	refer to service manual		
	17.50"	6"	bottom	bottom	lower	
	17.50"	8"	top	bottom	lower	
	18.50"	8"	†	refer to ser	vice manual	
	19.50"	8"	†	refer to ser	refer to service manual	
24"	16.50"	6"	†	refer to ser	refer to service manual	
	17.50"	6"	bottom	middle	upper	
	17.50"	8"	top	middle	upper	
	18.50"	8"	†	refer to ser	vice manual	
	19.50"	8"	bottom	bottom	lower	
nee-to-heel**	Min.	Max. (range adjusts in 1" increments)				
70 degree tapered	14.75"	18.75"				
70 degree	12.50"	17.50"				
60 degree	13.75"	17.50"				
Elevating	14.50"	20.25"				

#### Arm Types

Conventional dual-point

Desk or Full Length

Fixed or Adjustable Height

Space-Saver dual-point

Desk or Full Length

Fixed or Adjustabl	e Height	
Back heights	15"-19"	
Weight capacity	250 lb.	
Weight	29 lb.	
Limited warranty	Lifetime on frame and crossbraces	

 $<sup>^*</sup>$ Measurements taken from outside of handrim to outside of handrim with 0 degrees camber



Invacare Corporation www.invacare.com

#### **USA**

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

### Canada

570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare is a registered trademark and Yes, you can., the blue medallion design and Patriot are trademarks of Invacare Corporation.

Specifications are subject to change without notification.

© 2003 Invacare Corporation

Form No. 03-051



"Invacare says 'Yes, you can",' and I couldn't have said it betten myself."

Carola Valmer

<sup>\*\*</sup>Measurement taken from top of seat rail to back of footrest

 $<sup>^{\</sup>dagger Not \ available \ as \ a}$  factory setting Download from Www.Somanuals.com. All Manuals Search And Download.

Free Manuals Download Website

http://myh66.com

http://usermanuals.us

http://www.somanuals.com

http://www.4manuals.cc

http://www.manual-lib.com

http://www.404manual.com

http://www.luxmanual.com

http://aubethermostatmanual.com

Golf course search by state

http://golfingnear.com

Email search by domain

http://emailbydomain.com

Auto manuals search

http://auto.somanuals.com

TV manuals search

http://tv.somanuals.com