

**Five-Year Commercial Warranty**



## Rehabilitation Treadmill


# PACKED WITH REHABILITATION

**Treadmill always starts at 0.0 mph and resets to 0.0 mph.**

**Reversing belt option is ideal for practitioners who prescribe downhill and retro-walking therapies.**

**AccuTrack contact heart rate monitoring option complements cardiopulmonary and respiratory therapies and conditioning without the need for wearing a chest strap.**

# TION SPECIFIC FEATURES



**15% grade electric elevation with variable elevation over-speed protection**

**Ultra high-torque drive moves users weighing up to 400 pounds at a true 0.1 to 12.0 mph.**

**Patient-attachable safety lanyard (treadmill stops when lanyard is pulled)**

**Remote start-stop switch option is ideal for modalities in which a therapist must manipulate the patient's legs while remaining in control of treadmill operation.**

## L7 Rehabilitation Treadmill Specifications

- Zero starting speed
- 3 HP continuous-duty motor
- Actual 0.1 to 12.0 mph speed range
- Ultra high-torque drive
- 400-pound user weight capacity
- 15% grade electric elevation
- Variable elevation over-speed protection
- 20" x 58" four-ply treadbelt
- Low 5¾" step-up height
- VFX<sup>Plus</sup> shock absorption system
- Extended parallel medical handrails
- Patient-attachable safety lanyard
- 2½" diameter steel rollers
- 1-inch thick reversible deck
- Rust-free aluminum frame
- Size: 32" x 77"
- Weight: 300 lbs.
- Choice of Titanium or Matte Textured Black finish

**Commercial Warranty:** 5 years parts / 1 year labor\*



**Electrical Requirements:** 120VAC, 60HZ, 15A (International voltages available)

**Certifications:** ETL certified to UL-1647 and CAN CSA 22.2 no. 68

\*Designed for use in rehabilitation settings. Heavy use pay-for-membership facilities such as health clubs, gyms, and YMCAs are not eligible for this warranty program.

## Rehabilitation Treadmill Options

- Reversing belt direction for downhill and retro walking
- Hand-held remote stop-start switch
- AccuTrack contact heart rate monitoring
- L8 frame with 22" x 63" four-ply treadbelt and 500-pound user weight capacity
- 100 micro-amp hospital isolation transformer
- **ElliptiMill ORTHOPEDIC SHOCK ABSORPTION SYSTEM**  
Provides the ultimate reduction of impact on joints and knees.



Orthopedic 4-ply, cushioned treadbelt  
5X magnification



## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>