

RESMED

ResScan™ Patient Management System

Greater insights, better outcomes




ResScan™ Patient Management System

Welcome to the world of ResScan

Compatible with all ResMed treatment devices, ResScan™ is much more than a software program. It's an **easy** and **flexible** patient management system providing excellent insights into therapy to enable better outcomes.

This unique patient management system offers:

- a clear view of your patients
- the right quality of objective information at the right time
- the tools to easily access, manage and use that information.



ResScan software with unique start screen



can USB adapter (also available as serial port)



Easy: user-friendly experience

The ResScan system brings intuitive design to patient management.

- Use one single software program for all ResMed treatment devices.
- Transport data easily, either via the ResScan Data Card or via the ResLink™ module.
- Let the unique start screen guide you through the software.
- Download all data in one step.
- Customize and save your preferred onscreen view of data.
- Create custom report templates only the first time, to see what you want every time.
- Email reports from within the software.
- Flag patient follow-ups and review at a glance on the start screen.

Flexible: choose the solution to suit your needs

Access information directly or remotely.

- ResScan Data Card*—with convenient mailback envelope
- Cable—with either the ResScan serial or USB adapter*
- ResLink module*—clips onto premium ResMed devices

Choose the level of monitoring your patients need.

- Standard monitoring: view trends up to 365 nights—summary data in ResScan software or on the flow generator LCD.
- Nightly monitoring: view therapy changes up to five nights—detailed flow generator data in ResScan software.
- Advanced nightly monitoring: view therapy changes with additional parameters including oximetry up to 30 nights—summary and detailed ResLink data in ResScan software.



ResScan Data Card* and ResScan Data Card reader*



ResScan Data Card mailback envelope*



ResLink*



ResScan reports



Oximetry option with ResLink

Complementing ResMed's wide range of treatment devices, the ResScan system includes patient management software as well as comprehensive data access options.

* Certain features and functions apply to specific ResMed flow generators and are based on regional availability.

Download from [Www.Somanuals.com](http://www.Somanuals.com). All Manuals Search And Download.

ResScan Patient Management System

Your window to valuable therapy benefits with ResMed devices



CPAP systems

- View usage data.
- Verify therapy success with apnea detection, leak and AHI data (premium systems).
- Fine-tune treatment with easy setting changes via the ResScan Data Card (for ResMed S8™ devices).

AutoSet™ systems

- Obtain objective evidence of how your patients respond to AutoSet technology.
- View how the device responds to your patients' changing needs.
- Capture graphic representations of AutoSet therapy for your records, your referral sources and patient education.

VPAP Adapt SV™ system

- Monitor the effectiveness of adaptive servo-ventilation therapy by comparing target ventilation against minute ventilation.
- Use detailed data to identify whether end expiratory pressure (EEP) requires adjustment to manage upper airway obstruction.

VPAP™ Auto system

- Monitor your respiratory insufficiency patients. Manage their changing needs using detailed ResLink data, which includes oximetry data.
- Review compliance and therapy outcomes for obstructive sleep apnea (OSA) patients.
- Remotely access VPAP Auto settings and download data via the ResScan Data Card.

RECOMMENDED COMPUTER REQUIREMENTS FOR RESSCAN SOFTWARE

- Windows® 98 SE, Windows NT 4.0 SP6, Windows 2000 SP2 or Windows XP
- Microsoft Internet Explorer 5.5
- 700 MHz Pentium II Processor (minimum 350 MHz) computer
- 192 MB RAM (min 128 MB RAM)
- 100 MB HDD
- Super VGA Monitor/Graphic Card (1024x768, 256 colours, 8MB video RAM)
- RS 232 Serial and USB port
- Printer
- CD ROM drive

PRODUCT CODES

USA

ResScan Software	31301
ResScan USB Adapter	22203
ResScan Serial Adapter	22201
Serial Cable	17952
ResScan Data Card Reader	22207
Spare S8 Blank Module	22904
Spare ResScan Data Card Module	22915
ResScan Data Card 10 Pack	22903
ResScan Data Card Envelopes	22906
ResLink	30924
S8 ResLink	22206
Oximeter (with Flex sensor cable)	1430050
Oximeter finger clip sensor	70275
USB SmartMedia Card Reader	22212
Card Reader USB hub	22211

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>