

6RUNNER® 10 (888WNLM)

SHOPRIDER®



The 6Runner® series powerchairs, all with 6 points of contact with the ground, provide the desired stability of traditional rear driven powerchairs, along with the maneuverability of modern mid wheel driven powerchairs. The 6Runner® series powerchairs deliver the best features available in today's powerchair technology all in one finished package. Tomorrow's products today, only from Shoprider.

Available Colors



Red



Blue



6RUNNER® 10 (888WNLM)

- Luxurious leatherette 20" high back recline executive seat
- Stable 4 post seating system
- Adjustable length, height and angle foot plate
- Adjustable width, height and flip-up arm rest
- Midwheel drive maneuverability
- 10" wheel
- 300 lb. capacity

6RUNNER® 10 SPECIFICATIONS

(888WNLM)

Overall Dimensions (L*1 x W x H*2)	in	30 x 23 x 44
Number, Size of Tires (front, middle, rear)	ea, in	2, 7"; 2, 10"; 2, 7"
Total Weight (w/36 Ah batteries, seat)	lb	191
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger	Amp	4 amp off board
Output Power of Motor (Reference Only)	hp	0.8 x2
Max. Single User Weight (On the Level Road)	lb	300
Max. Speed*3 (On the Level Road)	mph	5
Suggested Climbing Angle, Max. Single User Weight	%, lbs	Shown as Illustration
Range*4 (Per Charge with 36Ah battery/90Kg User Weight) (After the battery ad mechanical moving parts fully break in)	mile	16.2
Turning Radius	in	15.75



Adjustable height




Many seating options available



Flip-up arm rest



Reclining seat



All specification subject to change without prior notice.

WARNING! Please make sure that you are physically fit to operate this unit. Consult your dealer /doctor for advice.

- *1. Include the footplate.
- *2. Depending on the height of the seat back.
- *3. When proceeding up any incline, please move your seat to the most forward position or if you have the deluxe seat, make sure that it is in the 90 degree (upright) position.
- *4. The actual driving range varies with the factors shown as below.

(A) The weight of occupant	(E) Ambient Temperature
(B) Level Ground	(F) The Way of Driving
(C) Battery Conditions	(G) ETC.
(D) Type of Charger	
- *5. The Centers for Medicare and Medicaid Service (CMS).
The Healthcare Common Procedure Coding System (HCPCS)

All specifications are subject to change without prior notice.

Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

- 1) Include the anti-tip wheel or the rear castor.
- 2) Driver weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown as below:
 - a) The weight of occupant
 - b) Ground surface
 - c) Battery capacity and conditions
 - d) Type of charger
 - e) Ambient temperature
 - f) The way of driving
 - g) If the battery and mechanical moving parts fully break in
 - h) Etc.

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>