

Patient Alarm For Bed and Chair

- Alarm includes ON/OFF switch and dual volumes.
- Audio alarm sounds when patient gets out of chair or bed.
- Uses one 9v battery (included).
- Six month limited warranty.

GF13605	Patient Alarm for Wheelchair	Pad Size 10"x15"	1 ea
GF13606	Patient Alarm for Bed	Pad size 11"x30"	1 ea



Patient Alarm With Pull Cord

- Attaches to bed/chair with self contained clip.
- Audio alarm sounds when pin connector is pulled from the unit.
- Includes ON/OFF switch and dual volumes.
- Uses one 9V battery (included).
- One year limited warranty.

GF13602	Patient Alarm with Pull Cord, Small	1 ea
----------------	-------------------------------------	------



Patient Alarm with Magnetic Pull Cord

- Attached to bed/chair with self contained clip and adjustable cord.
- Audio alarm sounds with magnetic pull switch is removed.
- Includes ON/OFF switch and dual volumes.
- Uses two AAA batteries (included).
- One year limited warranty.



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2006 © Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>